REF. 33160 · REF. 33163  
REF. 33104 · REF. 33116

## Knee support

genusKILL 63

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## Knee support

## DECLARATION OF CONFORMITY

As the manufacturer, ORTHOSERVICE AG states under its sole responsibility that this product is a Class I medical device that has been manufactured in accordance with the requirements of Directives 93/42/EEC and 2007/47/EC. These instructions have been written in conformity with these directives, and their aim is to ensure the safe and proper use of this medical device.

## TRADEMARKS OF MATERIALS

Velcro® is a registered trademark of Velcro Industries B.V.

## PRECAUTIONS

It is advised that the pressure exerted by the product does not act on parts of the body suffering from wounds, swelling or lumps. The product should not be tightened too firmly so as to avoid the development of pressure points or compression of the nerves or blood vessels underneath. If in doubt about how to apply, contact an orthopaedic technician. Please read the composition of the product on the inside label carefully. The device should not be worn in the vicinity of open flames or strong electromagnetic fields (only 33104, 33116).

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## WARRANTY

This warranty is governed by the statutory provisions of the country in which the product was purchased. Where relevant to that country's regulations concerning warranties and the relationship between the seller and buyer shall apply. If you believe that the terms of the warranty may be applicable, first contact the retailer where the product was purchased. Orthoservice will repair or replace all or part of the device and its accessories if defects in materials and workmanship arise during the period of six months from the date of purchase.

## SELECTION/SIZE

Code	REF. 33160 · REF. 33163 · REF. 33104 · REF. 33116					
Size	XS	S	M	L	XL	XXL
Thigh circ. taken 15 cm above mid-patella cm	33/39	39/45	45/52	52/58	58/66	66/75
Length cm Genuskill 04/16	34	34	34	34	34	34
Length cm Genuskill 60/63	27	27	27	27	27	27
Colour	grey					

Genuskill 04/60/63: fits right and left · Genuskill 16: specify right or left

## MAINTENANCE

- Do not bleach
  - No chemical cleaning
  - Do not iron
  - Do not tumble-dry
  - Waschrichtung: Handwasche bis 30° C mit neutraler Seife. Dry away from heat sources.
- Dispose of this product and its parts responsibly

## INDICATIONS GENUSKILL 04/60/63

- Gonalgies due to contusions
- Early stage arthrosis
- Solo Genuskill 04/63: Femoro-patellar chondropathy and chondromalacia in the knee
- Solo Genuskill 04/63: Prevention of displacement of the patella during sporting activity
- Solo Genuskill 04: Légère instabilité rotulienne
- Solo Genuskill 60: 1<sup>re</sup> grade sprains to collateral ligaments
- Solo Genuskill 60: Tendinitis

## INDICATIONS GENUSKILL 16

- Lateral-medial instability of the patella
- Femoro-patellar arthrosis
- Gonarthrosis

## CONTRAINDICATIONS

Currently no known

## CHARACTERISTICS AND MATERIALS

- BreathSkill fabric support:
  - soft and breathable
  - 3D weave
  - a special woven mesh covering for a stronger, more durable design
- CompreSkill on the back of the support:
  - Stretchy multi-fibre yarn for compression
  - Ultra-low-profile, hassle-free design (knee flexion)
  - Fully adaptable
  - Perfect snug fit around the joint
  - No open kneecap Genuskill 60 only
  - Shorter design Genuskill 60 / 63 only
  - Open kneecap Genuskill 63 only
  - FlexSkill on the front of the support (bonded with a technical fabric) with special low-profile bonded seams. Genuskill 04 only
  - Patella stabiliser included Genuskill 04 only
  - Lateral and medial stays keep the knee support correctly positioned during wear. Genuskill 04 / 16 only
  - Patella alignment system Genuskill 16 only
  - lateral J-shaped foam stabiliser
  - support and pressure for the kneecap with hook-and-loop straps (customisable sizing)

## PUTTING ON THE APPLIANCE

## Genuskill 60

- 1 Slide the tubular support onto the knee with the orange label towards the top (fig. A).
- 2 Position the knee support so that the CompreSkill fabric is behind the knee.

## Genuskill 04 - Genuskill 63

- 1 Slide the tubular support onto the knee with the orange label towards the top (fig. B).
- 2 Position the knee support so that the kneecap is perfectly centred in the opening provided (fig. C).

## Genuskill 16

- 1 Undo the patella straps and fasten them back on themselves for the time being.
- 2 Slide the tubular support onto the knee with the orange label towards the top (fig. D) taking care to align the opening with the kneecap. Lateral instability: the stabiliser must be lateral to the patella. Mediale instability: the stabiliser must be medial to the patella (fig. E).
- 3 Pass the upper strap through the buckle provided, tighten, and fasten with the Velcro® (fig. F).
- 4 Repeat this procedure with the lower strap (fig. G).

N.B.: if the straps are too long, they can be cut to size and the corresponding Y-shaped Velcro® can be reattached (fig. H.).

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