
Invacare® Kirilos / Atlante

Kirilos, Atlante Electric, Atlante Hydraulic

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This manual **MUST** be given to the user of the product.
BEFORE using this product, read this manual and save for future reference.

Yes, you can.®

©2017 Invacare Corporation

All rights reserved. Republication, duplication or modification in whole or in part is prohibited without prior written permission from Invacare. Trademarks are identified by ™ and ®. All trademarks are owned by or licensed to Invacare Corporation or its subsidiaries unless otherwise noted.

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Contents

This manual MUST be given to the user of the product.
BEFORE using this product, read this manual and save for future reference.

1 General	4	5.4 Closing/Opening Legs	16
1.1 Introduction	4	5.4.1 Closing/Opening legs manually	16
1.2 Symbols in this manual	4	5.5 Performing an Emergency Stop	16
1.3 Intended use	4	5.6 Activating an Emergency Release	16
1.4 Scope of delivery	5	5.7 Charging the battery	17
1.5 Service Life	6	6 Lifting the Patient	19
1.6 Warranty Information	6	6.1 Safe Lifting	19
1.7 Limitation of liability	6	6.2 Preparing to Lift	21
2 Safety	7	6.3 Attaching the Sling to the Lift	23
2.1 Safety information	7	6.4 Lifting and Transferring the Patient From a Bed	24
2.2 Operating Information	7	6.4.1 Floor Transfers (Lifting from the Floor)	28
2.2.1 General	7	6.4.2 Commode Transfer Guidelines	30
2.2.2 Pinch Points and Positioning	8	6.4.3 Wheelchair Transfer	32
2.3 Radio Frequency Interference	8	7 Maintenance	34
2.4 Product labeling	8	7.1 Maintenance and Safety Inspection	34
7.1.1 Safety inspection checklist		7.2 Lubricating the Lift	37
3 Components	10	7.3 Cleaning the Sling and Lift	37
3.1 Main parts of the lift	10	7.4 Checking and tightening bolts	38
3.2 Accessories	11	7.5 Maintaining the Hydraulic Pump	38
4 Setup	12	8 After Use	39
4.1 Safe Assembly	12	8.1 Transportation and Storage	39
4.2 Assembling the Lift	12	8.2 Disposal	39
4.3 Disassembling the Lift	13	8.3 Reuse	39
5 Usage	14	9 Troubleshooting	40
5.1 Introduction	14	9.1 Identifying and repairing faults	40
5.2 Locking/Unlocking the Rear Casters	14	10 Technical data	42
5.3 Raising/Lowering the Lift	14	10.1 Dimensions and weight	42
5.3.1 Raising/Lowering an Electric Lift	14	10.2 Electrical system	43
5.3.2 Raising/Lowering a Hydraulic Lift	15	10.3 Hydraulic system	44
		10.4 Environmental conditions	44
		10.5 Materials	44
		10.6 Electromagnetic compliance (EMC) information	44
		10.7 Electromagnetic compliance (EMC)	46

1 General

1.1 Introduction

Thank you for choosing an Invacare product.

This user manual contains important information about the handling of the product. In order to ensure safety when using the product, read the user manual carefully and follow the safety instructions.

Please note that there may be sections in this user manual, which are not relevant to your product, since this manual applies to all available models (on the date of printing). If not otherwise stated, each section in this manual refers to all models of the product.

1.2 Symbols in this manual

Symbols and signal words are used in this manual and apply to hazards or unsafe practices which could result in personal injury or property damage. See the information below for definitions of the signal words.

WARNING

Indicates a hazardous situation that could result in serious injury or death if it is not avoided.

CAUTION

Indicates a hazardous situation that could result in minor or slight injury if it is not avoided.

IMPORTANT

Indicates a hazardous situation that could result in damage to property if it is not avoided.

Tips and Recommendations

Gives useful tips, recommendations and information for efficient, trouble-free use.

This product complies with Directive 93/42/EEC concerning medical devices. The launch date of this product is stated in the CE declaration of conformity.

Manufacturer of the product.

1.3 Intended use

Mobile patient lifts are indoor transfer devices, designed to be used in most of the common lifting situations in hospitals, nursing facilities and domestic areas, for example:

- Between the bed and a wheelchair
- To and from the toilet
- Lowering and raising patients to/from the floor


The mobile patient lift can be used to transfer and position completely or partially immobile patients, who cannot be transferred with other types of lifts or transfer aids. All position changes are possible without assistance of the patient. The mobile patient lift is only intended to lift patients up to the maximum weight limit stated in technical data. There are no known contraindications for this product.

Selecting the appropriate slings and accessories for each individual is important to assure safety when using a patient lift. See Invacare's sling and accessory user manuals for further information on those devices.

Invacare recommends that the patient be transferred to a shower chair or other means for bathing.

The mobile patient lift can be turned (rotated) in place for transfers with limited floor space.

1.4 Scope of delivery

Ⓐ	Mast assembly including boom, hanger bar and actuator or hydraulic pump (1 piece)
Ⓑ	Battery (1 piece, Kirilos and Atlante Electric only)
Ⓒ	Hand control (1 piece, Kirilos and Atlante Electric only)
Ⓓ	Leg spreader handle (1 piece)
Ⓔ	Base with legs (1 piece)
Ⓕ - not shown	Cable for charging (1 piece, Kirilos and Atlante Electric only)
Ⓖ - not shown	Sling (1 piece)
Ⓗ - not shown	User Manual (1 piece)

The items listed in the table are included in your package.

1.5 Service Life

WARNING!

Risk of Injury or Damage

The product has been tested for the service life stated in this manual. Use of the product beyond this time period may cause product damage and injury.

- Only use the product for the service life stated in this manual. DO NOT exceed the service life of the product.
- Remove the product from service when the service life has been met. DO NOT continue to use the product.
- Perform all maintenance according to the recommended schedule in this manual.

The product has an expected lifetime of 8 years when used in accordance with safety instructions, maintenance intervals and correct use, stated in this manual. The effective product lifetime can vary according to frequency and intensity of use.

Actuator Service Life

Number of lifts per day	Service life of actuator (in years)
1–2	10
3	9
4	6
5	5
6	4

Number of lifts per day	Service life of actuator (in years)
7–9	3
10–13	2
14–27	1

1.6 Warranty Information

Terms and conditions of the warranty are part of the general terms and conditions particular to the individual countries in which this product is sold.

Contact information for your local Invacare office is located inside the back cover of this manual.

1.7 Limitation of liability

Invacare accepts no liability for damage arising from:

- Non-compliance with the user manual
- Incorrect use
- Natural wear and tear
- Incorrect assembly or set-up by the purchaser or a third party
- Technical modifications
- Unauthorized modifications and/or use of unsuitable spare parts

2 Safety

2.1 Safety information

WARNING!

– Do not use this product or any available optional equipment without first completely reading and understanding these instructions and any additional instructional material such as user manuals, service manuals or instruction sheets supplied with this product or optional equipment. If you are unable to understand the warnings, cautions or instructions, contact a healthcare professional, dealer or technical personnel before attempting to use this equipment - otherwise, injury or damage may occur.

The information contained in this document is subject to change without notice.

Check all parts for shipping damage before using. In case of damage, Do not use the equipment. Contact the dealer or Invacare representative for further instructions.

2.2 Operating Information

This section of the manual contains general safety information about your product. For specific safety information, refer to the appropriate section of the manual and procedures within that section. For instance, for safety information related to assembling the lift, refer to the Setup section of the manual.

2.2.1 General

WARNING!

Risk of Falling

The Invacare mobile patient lift is NOT a transport device. It is intended to transfer an individual from one resting surface to another (such as a bed to a wheelchair).

Invacare slings and patient lift accessories are specifically designed to be used in conjunction with Invacare patient lifts.

WARNING!

Risk of Falling

Do not attempt any transfer without approval of the patient's physician, nurse or medical assistant. Thoroughly read the instructions in this User Manual, observe a trained team of experts perform the lifting procedures and then perform the entire lift procedure several times with proper supervision and a capable individual acting as a patient.


- Use common sense in all lifts. Special care must be taken with people with disabilities who cannot cooperate while being lifted.
- Use steering handle on the mast at all times to push or pull the patient lift.
- Be sure to check the sling attachments each time the sling is removed and replaced, to ensure that it is properly attached before the patient is removed from a stationary object (bed, chair or commode).

WARNING!

- The patient lift can be used indoors on a level surface. If the patient lift is used in the area of a shower or bath, ensure that the patient lift is wiped clean of any moisture after use.
- Periodically inspect all components of the patient lift for signs of corrosion. Replace all parts that are corroded or damaged.

2.2.2 Pinch Points and Positioning

WARNING!

Risk of Injury

- Pinch points are present in several locations on the lift and fingers could be pinched. The hanger bar can move suddenly and cause injury.
- Always keep hands and fingers clear of moving parts.
 - When positioning lift, be aware of the position of the hanger bar and the patient.

2.3 Radio Frequency Interference

WARNING!

Risk of Injury or Damage

- Most electronic equipment is influenced by Radio Frequency Interference (RFI). CAUTION should be exercised with regard to the use of portable communication equipment in the area around such equipment, otherwise injury or damage may occur. If RFI causes erratic behavior:
- PUSH the Red Power Switch OFF IMMEDIATELY.
 - DO NOT turn the Power Switch ON while transmission is in progress.

2.4 Product labeling

Label location

The product label is placed on the mast of the patient lift and contains the main product information, including technical data.

Symbols on the product label

Manufacturer address

Date of manufacture

Type Reference Number

Serial Number

Safe Working Load

Class II equipment
(Kirilos and Atlante Electric only)

Type B applied part
(Kirilos and Atlante Electric only)

This product complies with Directive 93/42/EEC concerning medical devices. The launch date of this product is stated in the EC declaration of conformity.

WEEE conform.
(Kirilos and Atlante Electric only)

Read Manual
(Atlante Hydraulic only)


Other labels and symbols

Read Manual
(Kirilos and Atlante Electric only)

3 Components

3.1 Main parts of the lift

(A)	Base
(B)	Leg
(C)	Caster
(D)	Rear caster with brake
(E)	Boom
(F)	Mast
(G)	Hanger bar
(H)	Hook for sling
(I)	Actuator (Kirilos and Atlante Electric only)
(J)	Control unit with battery (Kirilos and Atlante Electric only)
(K)	Emergency Stop (Kirilos and Atlante Electric only)
(L)	Steering handle
(M)	Hand control (Kirilos and Atlante Electric only)
(N)	Leg spreader handle
(O)	Hydraulic pump – replaces actuator (I) for Atlante Hydraulic
(P)	Manual emergency lowering device (Kirilos and Atlante Electric only)

3.2 Accessories

CAUTION!

Compatibility of slings and hanger bars / strap hooks

Invacare® uses a "Loop and Coat Hanger Bar / Strap Hook System" as do many other manufacturers. Therefore other suitable patient transfer systems (slings), manufactured by other companies, can be used on the Invacare patient lift range as well.

However we do recommend:

- A risk assessment is always to be carried out by a professional prior to issuing lifting equipment. It is important that the Task, Individual, Load, Environment and Equipment are considered in the risk assessment.
- Always choose the sling design and size according to the patient's weight, size and physical ability whilst considering the type of transfer to be carried out.
- Do only use slings that are suitable for a "Loop and Coat Hanger Bar / Strap Hook System".
- Do not use slings for "Keyhole Hanger Bar" or for "Tilting Frame Hanger Bar" designs.

Available Accessories

Sling models for "Loop and Coat Hanger Bar System":

- Full body support slings – without head support
- Full body support slings –with head support
- Slings for dress/toileting – with or without head support
- Slings for amputee

4 Setup

4.1 Safe Assembly

WARNING! **Risk of Injury**

- Improper assembly may cause injury or damage.
- Assembly **MUST** be performed only by qualified personnel.
 - Use only Invacare parts in the assembly of this patient lift. The base legs, the mast, boom, pump assembly and the hanger bar are manufactured to specifications that assure correct alignment of all parts for safe functional operation.
 - **DO NOT** overtighten the mounting hardware. This will damage the mounting bracket.

If there are any issues or questions during assembly, contact a local Invacare representative. Refer to the contact information on the last page in this manual.


4.2 Assembling the Lift

WARNING!

- Check all parts for visible defects or damage before assembly. In case of any damage, do not use the product and contact Invacare.
- Make sure the emergency stop is activated before assembly or disassembly.
- Take care when lifting components during assembly. Some parts are heavy. Always remember to adopt the correct lifting position.

Assembling the mast to the base


1. Put the base **A** on a level surface and lock the rear castors **B**.
2. Loosen the mast knob **C**, but leave it screwed into the base **A**.
3. Insert the tube of the mast **D** into the hole of the base **A**.
4. While supporting the mast assembly **D**, tighten the mast knob **C**.
5. Place the leg spreader handle **E** onto its support **F** and push the spring pin, to lock the handle **E** on its support **F**.

Installing the actuator or hydraulic pump to the boom


2 x 17 mm wrench

1. Loosen the nut **A** and remove bolt **B** from the boom mounting bracket **C**.
2. Place the actuator **D** or hydraulic pump **E** in the boom mounting bracket **C** and align the holes.
3. Reinstall the bolt **B** and tighten the nut **A**.

Wiring

(Kirilos and Atlante Electric only)

1. Plug the hand control cable on socket **A** and motor cable on socket **B** control unit.
2. Place the battery **C** on the control unit. Push the top of the battery against the mounting bracket until there is an audible click.

4.3 Disassembling the Lift

1. Put the base on a level surface so that all castors make contact with the floor.
2. Lock the rear castors.
3. Reverse the relevant steps of the assembling procedure for your model.


5 Usage

5.1 Introduction

The operation of the patient lift is an easy and safe procedure.

-  Before using the lift with a patient, refer to the following procedures for safety information and instruction:
- Operating Information in the Safety section
 - Lifting and Transferring the Patient in the Lifting the Patient section

5.2 Locking/Unlocking the Rear Casters

1. To Lock the Caster (A): Press tab (B) down.
2. To Unlock the Caster (A): Press tab (B) up.

5.3 Raising/Lowering the Lift


WARNING! **Risk of Injury**

The lift could tip and endanger the patient and assistants.

– Invacare does recommend that the rear casters be left unlocked during lifting procedures to allow the patient lift to stabilize itself when the patient is initially lifted from a chair, bed or any stationary object.

5.3.1 Raising/Lowering an Electric Lift

The hand control is used to raise or lower the lift.


1. To raise the lift — Press and hold the UP (A) button to raise the boom and the patient.
2. To lower the lift — Press and hold the DOWN (B) button to lower the boom and the patient.

Release the button to stop raising or lowering the lift.

5.3.2 Raising/Lowering a Hydraulic Lift

The pump handle and control valve are used to raise or lower the hydraulic lift.

Refer to the safety information in the Introduction of the Usage section before performing this procedure.

There are two controls on the pump assembly:

- The control valve (A)
- The pump handle (B)

To raise the lift:

1. Turn the control valve (A) clockwise to the CLOSED position.
2. Move the pump handle (B) up and down to elevate the boom (C).

To lower the lift:

1. Turn the control valve (A) counterclockwise to the OPEN position.
2. The weight of the patient in the sling (D) will cause the boom to descend.
3. If there is no patient in the sling, pull down on the boom (C).

The rate of descent of the boom is controlled by the amount that the control valve is open.

A safety gate is part of the hydraulic system that controls the maximum descent of the boom regardless of how far the control valve is opened.

5.4 Closing/Opening Legs


WARNING! **Risk of Injury**

The lift could tip and endanger the patient and assistants.

- The legs of the lift must be in the maximum open position for optimum stability and safety. If it is necessary to close the legs of the lift to maneuver the lift under a bed, close the legs of the lift only as long as it takes to position the lift over the patient and lift the patient off the surface of the bed. When the legs of the lift are no longer under the bed, return the legs of the lift to the maximum open position.

5.4.1 Closing/Opening legs manually

1. To open the legs, pull the leg spreader handle Ⓐ off and then to the left until it locks in the notch of the bracket.
2. To close the legs, push the leg spreader handle Ⓐ off and then to the right until it locks in the notch of the bracket.

The shifter handle **MUST** lock into its mounting slot to lock the legs in the full closed position.


WARNING! **Risk of Injury**

The lift could tip and endanger the patient and assistants.

- If the leg spreader handle is **NOT** positioned completely into its mounting slot, **DO NOT** use the patient lift until the leg spreader handle is properly seated and the legs of the patient lift are locked in place.

5.5 Performing an Emergency Stop

(Kirilos and Atlante Electric only)


1. Press the red emergency button Ⓐ on the control unit to stop the boom and patient from raising or lowering.
2. To reset, rotate the emergency button clockwise.

5.6 Activating an Emergency Release

Electrical emergency release with control unit


(Kirilos and Atlante Electric only)

1. Insert a pen into the hole labeled Emergency Down (A) on the control unit (B)

Secondary (manual) emergency release with actuator

(Kirilos and Atlante Electric only)

In case of partial or total power failure, or if the battery runs down while using the lift, it is equipped with a manual emergency release system located at the bottom of the actuator.

It is recommended that the primary emergency release be used. The secondary (manual) emergency release is only a back-up to the primary emergency release.


1. Rotate the emergency release knob (A) counterclockwise.

The boom will lower slowly without additional downward force on the boom.

Emergency release with hydraulic pump

(Atlante Hydraulic only)

1. Turn the control valve (A) counterclockwise to the OPEN position.

5.7 Charging the battery

(Kirilos and Atlante Electric only)


! **IMPORTANT!**

- Make sure the emergency stop is not activated while charging the battery.
- Do not use the lift while charging the battery.
- Make sure that charging takes place in a room with good air ventilation.
- Do not use or move the lift without unplugging from the wall outlet after charging.
- Do not attempt to use the lift if the battery housing is damaged.
- Replace a damaged battery housing before further use.

It is recommended to charge the battery daily to ensure optimal use of the lift and prolong the life of the battery. Furthermore, it is recommended to charge the battery before first use.

Control unit

The control unit is equipped with a sound signal. A beep indicates that the battery has low capacity, but lowering the patient is still possible. It is recommended to charge the batteries as soon as the sound signal is heard.

1. Plug the cable for charging © into a power outlet.

The battery will charge in approximately 4 hours. The charger stops automatically when the batteries are fully charged. The right yellow charge diode ① will continuous light during charging and switch off when fully charged. The left green diode ② will light continuously while the control unit is connected to the mains, and lights up when any button on the hand control is pressed.

2. Disconnect the power cord from the power outlet after the battery has been fully charged.

6 Lifting the Patient

6.1 Safe Lifting

**WARNING!****Risk of Death, Injury or Damage**

Improper use of this product may cause death, injury or damage.

The lift could tip and endanger the patient and assistants.

The Invacare mobile patient lift is NOT a transport device. It is intended to transfer an individual from one resting surface to another (such as a bed to a wheelchair).

- Wheelchair and bed wheel locks MUST be in a locked position before lowering the patient onto or lifting the patient off of the wheelchair or bed to prevent the wheelchair or bed from moving during transfer.
- Before transferring, check that the wheelchair (bed, commode or other surface) weight capacity can withstand the patient's weight.

**WARNING!****Risk of Death, Injury or Damage**

Improper use of this product may cause death, injury or damage.

The lift could tip and endanger the patient and assistants.

The Invacare mobile patient lift is NOT a transport device. It is intended to transfer an individual from one resting surface to another (such as a bed to a wheelchair).

- The legs of the lift must be in the maximum open position for optimum stability and safety. If it is necessary to close the legs of the lift to maneuver the lift under a bed, close the legs of the lift only as long as it takes to position the lift over the patient and lift the patient off the surface of the bed. When the legs of the lift are no longer under the bed, return the legs of the lift to the maximum open position.
- Invacare recommends locking the rear casters of the lift ONLY when positioning or removing the sling from around the patient.
- Invacare recommends that the rear casters be left unlocked during lifting procedures to allow the patient lift to stabilize itself when the patient is initially lifted from a chair, bed or any stationary object.

WARNING!

Risk of Death, Injury or Damage

Improper use of this product may cause death, injury or damage.

- Use the steering handle on the mast at all times to push or pull the product.
- Avoid using the product on an incline. Invacare recommends that the product only be used on a flat surface.
- During transfer, with the patient suspended in a sling attached to the lift, DO NOT roll lift over uneven surfaces that could cause the lift to tip over.

WARNING!

Risk of Injury or Damage

Damage to parts of the lift (pendant, casters, etc.) caused by impact with the floor, walls or other stationary objects may cause damage to the product and lead to injury.

- DO NOT allow parts of the lift to impact the floor, walls or other stationary objects.
- ALWAYS store the pendant properly when not in use.

WARNING!

Risk of Death

The pendant cord can cause injury if improperly positioned and secured.

- ALWAYS be aware of the location of the pendant cord in relation to the patient and caregivers.
- DO NOT allow the pendant cord to become entangled around the patient and caregivers.
- The pendant must be secured properly. ALWAYS store the pendant properly when not in use.

WARNING!

Risk of Entrapment or Strangulation

Items in the patient's surroundings can cause entrapment strangulation during lifting. To avoid entrapment or strangulation:

- Before lifting, check that the patient is completely free of his/her surroundings.

WARNING!

Risk of Entrapment

There is a risk of entrapment between the hanger bar hooks and the sling.

- Use caution when lifting.
- NEVER put hands or fingers on or near the hooks when lifting.
- Ensure the patient's hands and fingers are away from the hooks before lifting.

6.2 Preparing to Lift


WARNING! **Risk of Injury**

During transfers and lift operation, the boom can impact the patient or caregivers and cause injury.


- ALWAYS be aware of the position of the boom during transfers.
- Ensure the boom is positioned in a way that it cannot impact the patient or bystanders.
- ALWAYS be aware of your body position in relation to the boom during transfers.

1. Refer to the Safety section of the manual and review the information in Safe Lifting in the Lifting the Patient section of the manual before proceeding further. Observe all warnings indicated.
2. Position the patient onto the sling (A). Refer to your sling user manual.

3. Unlock the rear casters. Refer to Locking/Unlocking the Rear Casters in the Usage section of the manual.

4. Open the legs of the lift. Refer to Closing/Opening Legs of the Lift in the Usage section of the manual.

5. Use the steering handle **B** to push the patient lift into position.


WARNING!


– When using the lift in conjunction with beds or wheelchairs, be aware of the position of the lift in relationship to those other devices so that the lift does not become entangled.

Before positioning the legs of the patient lift under a bed, make sure that the area is clear of any obstructions.

6. Lower the patient lift for easy attachment of the sling.

7. Lock the rear casters. Refer to Locking/Unlocking the Rear Casters in the Usage section of the manual.
8. Attach the sling to the lift. Refer to Attaching the Sling to the Lift in the Lifting the Patient section of the manual.

6.3 Attaching the Sling to the Lift

**WARNING!****Risk of Injury or Death**

Improperly attached, improperly adjusted, or damaged slings can cause the patient to fall or cause injury to assistants.


- Use an Invacare approved sling that is recommended by the individual's doctor, nurse or medical assistant for the comfort and safety of the individual being lifted.
- Invacare slings and patient lift accessories are specifically designed to be used in conjunction with Invacare patient lifts.
- After each laundering (in accordance with instructions on the sling), inspect sling(s) for wear, tear, and loose stitching.
- Bleached, torn, cut, frayed, or broken slings are unsafe and could result in injury. Discard immediately.
- DO NOT alter slings.
- Be sure to check the sling attachments each time the sling is removed and replaced, to ensure that it is properly attached before the patient is removed from a stationary object (bed, chair or commode).
- Position the patient in the sling as directed by the instructions provided with the sling.
- Adjustments for patient safety and comfort should be made before moving the patient.

**WARNING!****Risk of Injury or Death**

Improperly attached, improperly adjusted, or damaged slings can cause the patient to fall or cause injury to assistants.

- DO NOT use any kind of plastic back incontinence pad or seating cushion between patient and sling material that may cause the patient to slide out of the sling during transfer.
- When connecting slings equipped with color coded straps to the patient lift, use shorter straps at the back of patient for upright support. Using longer straps will leave less support for the patient's back and the patient will be in a more reclined position. The loops of the sling are color coded and can be used to place patient in various positions. The colors make it easy to connect both sides of the sling equally. Make sure that there is sufficient head support when lifting a patient.
- The hanger bar **MUST** be attached to the lift **BEFORE** attaching the sling.

1. Place the straps (A) of the sling (B) over hooks (C) of the hanger bar (D). The slings may be equipped with color coded straps to assist with proper attachment. Match the corresponding straps on each side of the sling for an even lift of the patient.

2. Lift and transfer the patient. Refer to Lifting and Transferring the Patient in the Lifting the Patient section of the manual.

6.4 Lifting and Transferring the Patient From a Bed

WARNING!

Risk of Injury

The lift could tip and endanger the patient and assistants.

- Refer to the safety information and instructions in the following procedures BEFORE performing this procedure:

Safe Lifting in the Lifting the Patient section of the manual

Preparing to Lift in the Lifting the Patient section of the manual

Attaching the Sling to the Lift in the Lifting the Patient section of the manual


Raising/Lowering the Lift in the Usage section of the manual


Invacare recommends at least two (2) individuals assist in transferring the patient into and out of a bed using this product.

1. Move the lift to the patient area, open legs and prepare to lift. Refer to Preparing to Lift in the Lifting the Patient section of the manual.
2. Attach the sling to the lift. Refer to Attaching the Slings to the Lift in the Lifting the Patient section of the manual.
3. Unlock the rear casters.

4. Lower the bed to the lowest position.

5. Lift the patient high enough to clear the stationary object with their weight fully supported by the lift. Refer to Raising/Lowering the Lift in the Usage section of the manual.


The boom will stay in position until the DOWN (


) button is pressed.


6. Before moving the patient, check again to make sure that the sling is properly connected to the hooks of the hanger bar. Refer to Attaching the Sling to the Lift in the Lifting the Patient section of the manual. If any attachments are not properly in place, lower the patient back onto the stationary object and correct the problem.


7. Using the steering handle, move the lift away from the stationary object.

8. Using the handles on the sling, turn the patient so that he/she faces the assistant operating the patient lift.


9. Lower the patient so that his/her feet rest on the base of the lift, straddling the mast.

The lower center of gravity provides stability making the patient feel more secure and the lift easier to move.

10. Move the patient lift with both hands firmly on the steering handle.


11. Read and understand the information in the Lifting the Patient section of the manual that pertains to transfer to or from specific types of surfaces BEFORE proceeding:

- Bed Transfer
- Floor Transfers (Lifting from the Floor)
- Commode Transfer Guidelines
- Wheelchair Transfer

12. Raise or lower the lift to position the patient over the stationary surface. Be sure to raise or lower the patient enough to clear the sides of the stationary object.


13. Lower the patient onto the stationary surface.


14. Ensure the patient is fully supported by the surface to which you are transferring.
 15. Lock the rear casters.
 16. Raise the bed to a good working height (usually hip height of caregivers).

17. Detach the sling from the hanger bar.

18. Unlock the rear casters.
 19. Move the lift away from the area.

6.4.1 Floor Transfers (Lifting from the Floor)


WARNING!
Risk of Injury


The lift can injure a patient who is on the floor or an assistant working on the floor.

– Always use two assistants when transferring a patient from the floor.


1. Perform these steps in addition to those in Lifting and Transferring the Patient From a Bed when transferring from the floor.
2. Determine if the patient has suffered any injuries from a fall. If no medical attention is needed, proceed with the transfer.
3. Position the sling **A** under the patient. Refer to the sling user manual for more information about positioning slings.

4. One assistant should have the patient bend his/her knees and raise his/her head off of the floor. This assistant should support the patient's head with a pillow **B**.


5. The other assistant should open the legs of the lift. Refer to Opening/Closing Legs of the Lift in the Usage section of the manual.


6. Unlock the rear casters.


7. Position the lift with one leg under the patient's head and the other leg under the patient's bent knees. Keep the sling straps inside of the legs of the lift.

8. Lower the boom so the hanger bar is directly over the patient's chest.

9. Attach the sling. Refer to Attaching the Sling to the Lift in the Lifting the Patient section of the manual.

10. Proceed with the transfer. Refer to Lifting and Transferring the Patient From a Bed in the Lifting the Patient section of the manual.

6.4.2 Commode Transfer Guidelines


 The Invacare patient lift is NOT intended as a transport device. If the bathroom facilities are NOT near the bed or if the patient lift cannot be easily maneuvered towards the commode, then the patient MUST be transferred to a wheelchair and transported to the bathroom facilities before using the patient lift again to position the patient on a standard commode.

 The slings with commode openings are designed to be used with either a commode chair or standard commode.


 Invacare recommends at least two (2) individuals assist in transferring the patient onto and off of a commode using this product.

1. Perform these steps when transferring to a commode in addition to those in Lifting and Transferring the Patient From a Bed.
2. Before transferring the patient, the patient lift should be guided to the bathroom facilities to check that it can be easily maneuvered towards the commode.


3. Attach the slings to the lift. Refer to Attaching the Slings to the Lift in the Lifting the Patient section in the manual.


4. Unlock the rear casters.
5. Elevate the patient high enough to clear the commode chair arms and have their weight supported by the patient lift. Refer to Raising/Lowering the Lift in the Usage section in the manual.

6. Position the lift so the legs are outside of the commode legs. The push handles of the lift should be opposite of the commode, as shown in the figure. Both assistants should help guide the patient over the commode.

7. Lower the patient onto the commode, leaving the sling attached to the hanger bar hooks. Invacare recommends that the sling remain connected to the hanger bar hooks during the patient's use of either the commode chair or standard commode.


8. When complete, recheck for correct sling attachment.

9. Raise the patient off of the commode.

10. When the patient is clear of the commode surface, use the steering handles to move the lift away from the commode.

11. To return the patient to the bed, wheelchair or other surface. Reverse the following procedures:
 - Lifting and Transferring the Patient from a Bed in the Lifting the Patient section of the manual
 - Raising/Lowering the Lift in the Usage section of the manual
 - Wheelchair Transfer in the Lifting the Patient section of the manual

6.4.3 Wheelchair Transfer


WARNING! **Risk of Injury**

- Before transferring, check that the wheelchair weight capacity can withstand the patient's weight.
- The wheelchair wheel locks **MUST** be in a locked position before lowering the patient into the wheelchair for transport.


Invacare recommends at least two (2) individuals assist in transferring the patient into and out of a wheelchair using this product.

1. Perform these steps in addition to those in Lifting and Transferring the Patient From a Bed when transferring to or from a wheelchair.
2. Engage the wheel locks of the wheelchair to prevent movement of the wheelchair.

3. Position the patient over the seat with their back against the back of the chair. Use the straps or handles on the side and the back of the sling to guide the patient's hips as far back as possible into the seat for proper positioning.

4. Begin to lower the patient.
5. With one assistant behind the chair and the other operating the patient lift, the assistant behind the chair will pull back on the grab handle (on selected models) or sides of the sling to seat the patient well into the back of the chair. This will maintain a good center of balance and prevent the chair from tipping forward.

7 Maintenance

7.1 Maintenance and Safety Inspection

WARNING!

Risk of Falling

Maintenance **MUST** be performed only by qualified personnel.

Improper assembly may cause injury or damage.

- Regular maintenance of patient lifts and accessories is necessary to assure proper operation.
- **DO NOT** overtighten the mounting hardware. This will damage the mounting bracket.

Service Interval

At normal daily operation, a service check-up should take place every year, according to the Safety Inspection Checklist. When performing annual or regular maintenance, all parts designed to carry load must be, as a minimum, tested with maximum load. All safety features must be checked according to EN ISO 10535:2006 Annex B.

General Maintenance

Regular cleaning will reveal loose or worn parts, enhance smooth operation and extend the life expectancy of the lift.

Follow the maintenance procedures described in this manual to keep your patient lift in continuous service.

The Invacare Patient Lift is designed to provide a maximum of safe, efficient and satisfactory service with minimum care and maintenance.

It is important to inspect all stressed parts, such as slings, hanger bar and any pivot points for signs of wear, cracking, fraying, deformation or deterioration. All parts of the Invacare patient lift are made of the best grades of steel, but metal to metal contact will wear after considerable use. Replace any defective parts immediately and ensure that the lift is not used until repairs are made. Refer to the Safety Inspection Checklist for specific information regarding wear items.

There is no adjustment or maintenance of either the casters or brakes, other than cleaning, lubrication and checking axle and swivel bolts for tightness. Remove all debris, etc. from the wheel and swivel bearings. If any parts are worn, replace these parts immediately.

If you question the safety of any part of the lift, contact your Dealer or Invacare representative immediately and advise him/her of your problem.

Daily Inspections

The patient lift should be checked each time it is used. Perform the following checks in addition to those listed in the Safety Inspection Checklist. If you question the safety of any part of the lift, do not use. Contact your Dealer or Invacare representative immediately.

- Visually inspect the patient lift. Check all parts for external damage or wear. If damage is found, do not use. Contact your Dealer or Invacare representative immediately.
- Check the emergency lowering function (both electrical and/or mechanical). Check all parts for external damage or wear. If damage is found, do not use. Contact your Dealer or Invacare representative immediately.
- Check that all hardware and attachment points for damage or wear. Check all parts for external damage or wear. If damage is found, do not use. Contact your Dealer or Invacare representative immediately.
- Verify that the hand control is functional (lifting and leg movements).
- Charge the battery every day the lift is used.
- Check the emergency stop function.

7.1.1 Safety inspection checklist

A person who is suitably and properly qualified and well acquainted with the design, use and care of the lift should perform periodic inspections.

Date of Inspection:	Initials:
THE CASTER BASE <ul style="list-style-type: none"> <input type="checkbox"/> Inspect for missing hardware. <input type="checkbox"/> Base opens/closes with ease. <input type="checkbox"/> Inspect casters and axle bolts for tightness. <input type="checkbox"/> Inspect casters for smooth swivel and roll. <input type="checkbox"/> Inspect and clear wheels of debris. <input type="checkbox"/> Inspect pivot joints for wear. 	
SLINGS AND HARDWARE <ul style="list-style-type: none"> <input type="checkbox"/> Check all sling attachments each time it is used to ensure proper connection and patient safety. <input type="checkbox"/> Inspect sling material for wear. <input type="checkbox"/> Inspect straps for wear. <input type="checkbox"/> Inspect stitching. 	
THE HYDRAULIC PUMP/ELECTRIC ACTUATOR ASSEMBLY <ul style="list-style-type: none"> <input type="checkbox"/> Check for leakage. <input type="checkbox"/> Inspect hardware on mast, boom and base. <input type="checkbox"/> Check for wear or deterioration. If damaged, return to factory. <input type="checkbox"/> Cycle to ensure smooth quiet operation of the electric actuator. <input type="checkbox"/> Check for smooth operation of the pump handle for hydraulic pumps. 	

THE BOOM

- Check all hardware and hanger bar supports.
- Inspect for bends or deflections.
- Inspect bolted joints of boom for wear.
- Inspect to ensure that the boom is centered between the base legs.
- Check the mast pivot bolt. Ensure that the bolt is tightly secured.
- Inspect pivot joints for wear.

THE MAST

- Mast must be securely assembled to boom.
- Inspect for bends or deflections.
- Inspect pivot joints for wear.

THE HANGER BAR

- Check the bolt/hooks for wear or damage.
- Check sling hooks for wear or deflection.
- Inspect pivot joints for wear.

LEG SPREADER HANDLE

- Ensure leg spreader handle operates smoothly.

CLEANING

- Whenever necessary.

7.2 Lubricating the Lift


WARNING!

Risk of Falling

Hydraulic oil or lubricant on the floor can cause a fall and injury.

- Wipe excess lubricant from the lift after lubrication.
- Use a facial tissue to wipe excess oil from around the hydraulic piston.
- If excess oil is leaking from the hydraulic pump, contact a dealer or qualified technician for service.

The Invacare lift is designed for minimum maintenance. However, a six month check and lubrication should ensure continued safety and reliability.

Keep lift and slings clean and in good working order. Any defect should be noted and reported to your Dealer or Invacare representative as soon as possible.

Refer to the figure for lubrication points. Lubricate all pivot points with a light grease (waterproof auto lubricant). Wipe all excess lubricant from lift surface.

1. Hanger Bar
2. Boom Mounting Bracket
3. Boom/Mast Mount
4. Mast Mounting Bracket

7.3 Cleaning the Sling and Lift

Cleaning the Sling

Refer to the washing instructions on the sling and to the sling manual for cleaning details.

Cleaning and Disinfecting the Lift

CAUTION!

Risk of Damage

- Motors, control unit and mounting parts can be damaged if the lift is cleaned improperly.
- Never use acids, alkaline or solvents for cleaning the lift.
 - Dry the lift carefully after cleaning.

To prevent cross-infection, the hoist must be cleaned and disinfected after each use.


A soft cloth, dampened with water and a small amount of mild detergent, is all that is needed to clean the patient lift. The lift can be cleaned with non-abrasive cleaners.

Never use acids, alkaline or solvents for cleaning the lift. Dry the lift carefully after cleaning.

Motors, control unit and mounting parts can be damaged if the lift is cleaned any other way than stated above.

The lift must be wiped with a moistened, firmly wrung cloth with ordinary household disinfectants. Only use disinfection detergents approved by the facility and follow the facility policy. For more information about the residence time and concentration of disinfectants, please contact your disinfectant dealer or the manufacturer of the disinfectant.

7.4 Checking and tightening bolts

1. Check that the bolts are through the brackets and the locknuts are tight and secure.
2. If needed, do one or more of the following:
 - Tighten locknut and back-off the locknut 1/8 of a turn.
 - Replace the locknut.

7.5 Maintaining the Hydraulic Pump

(Atlante Hydraulic only)

WARNING!

Risk of Damage

The pump is sealed at the factory. To avoid voiding the warranty or causing product damage:

- If service is required, the pump unit **MUST** be returned to Invacare for repair.
- **DO NOT** attempt to open the pump or obtain local service. This will void the warranty and might result in damage and costly repair. Consult your dealer or Invacare for further information.

Risk of Falling

Hydraulic oil or lubricant on the floor can cause a fall and injury.

- Wipe excess lubricant from the lift after lubrication.
- Use a tissue to wipe excess oil from around the hydraulic piston.
- If excess oil is leaking from the hydraulic pump, contact a dealer or qualified technician for service.

All parts of the hydraulic pump are precision machined, then carefully assembled and tested to ensure reliable service. The pump assembly is completely enclosed and sealed with neoprene rings to prevent leakage of hydraulic oil. A small amount of oil (about a drop) will accumulate around the piston from time to time and should be removed with a tissue.

8 After Use

8.1 Transportation and Storage

During transportation, or when the patient lift is not to be used for some time, the emergency stop button should be pushed in.

The patient lift must be stored at normal room temperature. If it is stored in a damp, cold or wet environment then the motor and other mounting parts may be prone to corrosion. Refer also to Environmental Conditions.

8.2 Disposal

**WARNING!****Environmental Hazard**

This product has been supplied from an environmentally aware manufacturer that complies with the Waste Electrical and Electronic Equipment (WEEE) Directive 2012/19/EU.

Device contains lead acid batteries.

This product may contain substances that could be harmful to the environment if disposed of in places (landfills) that are not appropriate according to legislation.

- DO NOT dispose of batteries in normal household waste. They MUST be taken to a proper disposal site. Contact your local waste management company for information.
- Please be environmentally responsible and recycle this product through your recycling facility at its end of life.

8.3 Reuse

This product is suitable for reuse. The maximum number of times it can be reused is dependent upon product condition. To prevent the transmission of infection, the patient lift and slings must be cleaned after each use. Before reuse or refurbishment of the lift, refer to Cleaning the Sling and Lift in the Maintenance section of the manual. Always provide the user manual with the reused or refurbished lift.

9 Troubleshooting

9.1 Identifying and repairing faults

WARNING!

– Only personnel having received the necessary instruction or training by Invacare must perform service and maintenance on the product.




Symptoms	Faults	Solution
Patient lift feels loose.	Mast/Base joint loose.	Refer to Assembling the Mast to the Base in the Setup section.
	Tie - Rods are loose.	Contact your Dealer or Invacare representative.
Casters/Brakes noisy or stiff.	Fluff or debris in bearings.	Clean casters from fluff and debris.
Noisy or dry sound from pivots.	Needs lubrication.	Refer to Lubricating the Lift in the Maintenance section.
Electric actuator fails to lift.	Hand-control or actuator connector loose.	Connect hand control or actuator connector. Ensure connectors are seated properly and fully connected.
	Battery low.	Charge batteries. Refer to Charging the battery in the Usage section.
	RED emergency stop button pressed IN.	Rotate RED emergency stop button CLOCKWISE until it pops out.
	Battery not connected properly to control box.	Reconnect the battery to the control box. Refer to Charging the battery in the Usage section.

Symptoms	Faults	Solution
	The connecting terminals are damaged.	Replace the battery pack. Refer Charging the battery in the Usage section.
	Actuator needs service or load exceeds max. user weight.	Contact your Dealer or Invacare representative.
Unusual noise from actuator.	Actuator is worn or damaged or spindle is bent.	Contact your Dealer or Invacare representative.
Boom will not lower in uppermost position.	Boom requires a minimum weight load to lower from the uppermost position.	Pull down slightly on the boom.
Boom will not lower during a power retraction.	Shoulder bolt at the junction of the boom and mast may not be properly installed.	Refer to Checking and Tightening Bolts in the Maintenance section.
The control unit emits a beeping sound during lifting, and the motor stops.	Max. load is exceeded	Reduce the load (and the lift will function normally).
Manual/Hydraulic pump fails to lift when pumped.	Control valve not fully closed.	Close the control valve.
	Manual/Hydraulic pump in need of replacement	Contact your Dealer or Invacare representative.
Oil leaking from hydraulics.	Manual/Hydraulic pump in need of replacement.	Contact your Dealer or Invacare representative.

 Contact your dealer or Invacare representative if the above does not solve your problems.

10 Technical data

10.1 Dimensions and weight

Dimensions [mm]	Kirilos	Atlante Electric	Atlante Hydraulic
Front / rear castor diameter	75	100	
Max. reach at 600 mm (a)	480		
Max. reach from base (b)	520		
Reach from base with legs spread to 700 mm (c)	260		
CSP* max. height / highest position (k)	1780	1820	
CSP* min. height / lowest position (l)	850	890	
Lifting range (m)	930		
Total width (open) centre to centre of castor (p)	960		
Internal width at maximum reach (q)	810		
Min. internal width (r)	580		
Min. free height (s)	30	60	
Height to upper edge of legs (t)	110	150	
Base length (u)	1170	1200	
Max. internal base length (v)	870	890	

Dimensions [mm]	Kirilos	Atlante Electric	Atlante Hydraulic
Total width (closed), external measure (w)	680	700	
Total width (closed), internal measure (x)	530	510	
Total width (open) external measure (y)	1030	1060	
Total width (open) internal measure (z)	900	870	
Turning radius	1230		
Min. height at Sling Hook-up	575	610	
Max. height at Sling Hook-up	1490	1530	

* CSP = Central Suspension Point

Weights [kg]	Kirilos	Atlante Electric	Atlante Hydraulic
Maximum lifting capacity	140		
Total weight incl. hanger bar	33		37

Weights [kg]	Kirilos	Atlante Electric	Atlante Hydraulic
Weight, mast assembly, incl. battery, boom, hanger bar and actuator or pump	17		21
Base section with legs	16		

10.2 Electrical system

	Kirilos	Atlante Electric
Voltage output	24 V DC, max. 250 VA	
Voltage supply	100 – 240 V AC, 50/60 Hz	
Maximum current input	400 mA	
Protection class (entire device)	IPX4	
Insulation class	Class II equipment	
	Type B applied part Applied part complying with the specified requirements for protection against electrical shock according to IEC60601-1.	
Sound level	45 – 50 dB (A)	
Working ability	40 full lifts without battery charge with batteries at 50% of full capacity	

	Kirilos	Atlante Electric
Intermittent (periodic motor operation)	10%, max. 2 min. / 18 min.	
Battery capacity	2.9 Ah	
Manual emergency lowering	Yes	
Electric emergency lowering / lifting	Yes	

10.3 Hydraulic system

Technical specifications of hydraulic pump

	Atlante Hydraulic
Push rod max. stroke	300 mm
Push force (dynamic)	310 daN. at 95 Bar
Maximum static load	390 daN. at 120 Bar
Rod load / lever load ratio	20 / 1
Weight	9 Kg. (9.2 daN)

10.4 Environmental conditions

	Storage and transportation	Operation
Temperature	-10°C to +50°C	+5°C to +40°C

Relative humidity	20% to 75%	20% to 90% at 30°C, not condensing
Atmospheric pressure	795 hPa to 1060 hPa	

10.5 Materials

Component	Material
Base, legs, mast and boom	Steel, powder coated
Hanger bar	Steel, powder coated and foam
Actuator housing, hand control, mast protector, casters and other plastic parts	material according to marking (PA, PP, PE)
Carabiner, bolts and nuts	Steel, rust protected, zinc-plated

10.6 Electromagnetic compliance (EMC) information

Medical Electrical Equipment needs to be installed and used according to the EMC information in this manual.

This equipment has been tested and found to comply with EMC limits specified by IEC/EN 60601-1-2 for Class B equipment.

Portable and mobile RF communications equipment can affect the operation of this equipment.

Other devices may experience interference from even the low levels of electromagnetic emissions permitted by the

above standard. To determine if the emission from the lift is causing the interference, run and stop running the lift. If the interference with the other device operation stops, then the lift is causing the interference. In such rare cases, interference may be reduced or corrected by the following:

- Reposition, relocate, or increase the separation between the devices.

10.7 Electromagnetic compliance (EMC)

Guidance and manufacturer's declaration – electromagnetic emission

The patient lift is intended for use in the electromagnetic environment specified below. The customer or the user of the patient lift should assure that it is used in such an environment.

Emissions test	Compliance	Electromagnetic environment - guidance
RF emissions CISPR 11(partly)	Group I	The patient lift uses RF energy only for its internal function. Therefore, its RF emissions are very low and are not likely to cause any interference in nearby electronic equipment.
RF emissions CISPR 11(partly)	Class B	The patient lift is suitable for use in all establishments including domestic establishments and those directly connected to the public low-voltage power supply network that supplies buildings used for domestic purposes.
Harmonic emissions IEC 61000-3-2	Class A	
Voltage fluctuations / flicker emissions IEC 61000-3-3	Complies	

Guidance and manufacturer's declaration – electromagnetic immunity

The patient lift is intended for use in the electromagnetic environment specified below. The customer or the user of the patient lift should assure that it is used in such an environment.

Immunity test	IEC 60601 test level	Compliance level	Electromagnetic environment – guidance
Electrostatic discharge (ESD) IEC 61000-4-2	± 6 kV contact ± 8 kV air	± 6 kV contact ± 8 kV air	Floors should be wood, concrete or ceramic tile. If floors are covered with synthetic material, the relative humidity should be at least 30 %.

<p>Electrostatic transient / burst IEC 61000-4-4</p>	<p>± 2 kV for power supply lines ± 1 kV for input/output lines</p>	<p>± 2 kV for power supply lines ± 1 kV for input/output lines</p>	<p>Mains power quality should be that of a typical commercial or hospital environment.</p>
<p>Surge IEC 61000-4-5</p>	<p>± 1 kV line(s) to line(s)</p>	<p>± 1 kV line(s) to line(s)</p>	<p>Mains power quality should be that of a typical commercial or hospital environment. Product is double-insulated. There are no other possible connections to earth</p>
<p>Voltage dips, short interruptions and voltage variations on power supply input lines IEC 61000-4-11</p>	<p>< 5% U_T (>95% dip in U_T) for 0,5 cycle 40% U_T (60% dip in U_T) for 5 cycles 70% U_T (30% dip in U_T) for 25 cycles < 5% U_T (>95% dip in U_T) for 5 sec</p>	<p>< 5% U_T (>95% dip in U_T) for 0,5 cycle 40% U_T (60% dip in U_T) for 5 cycles 70% U_T (30% dip in U_T) for 25 cycles < 5% U_T (>95% dip in U_T) for 5 sec</p>	<p>Mains power quality should be that of a typical commercial or hospital environment. If the user of the patient lift requires continued operation during power mains interruptions, it is recommended that the patient lift be powered from an un-interruptible power supply or a battery. U_T is the a. c. mains voltage prior to application of the test level.</p>
<p>Power frequency (50/60 Hz) magnetic field IEC 61000-4-8</p>	<p>3 A/m</p>	<p>30 A/m</p>	<p>Power frequency magnetic fields should be at levels characteristic of a typical location in a typical commercial or hospital environment.</p>

<p>Conducted RF IEC 61000-4-6</p> <p>Radiated RF IEC 61000-4-3</p>	<p>3 V</p> <p>3 V/m</p>	<p>3 V</p> <p>3 V/m</p>	<p>Portable and mobile RF communications equipment should be used no closer to any part of the patient lift including cables, than the recommended separation distance calculated from the equation applicable to the frequency of the transmitter.</p> <p>Recommended separation distance:</p> $d = \left[\frac{3,5}{V_1} \right] \sqrt{P}$ $d = \left[\frac{3,5}{E_1} \right] \sqrt{P} \quad 80 \text{ MHz to } 800 \text{ MHz}$ $d = \left[\frac{7}{E_1} \right] \sqrt{P} \quad 800 \text{ MHz to } 2,5 \text{ GHz}$ <p>where P is the maximum output power rating of the transmitter in watts (W) according to the transmitter manufacturer and d is the recommended separation distance in metres (m).^b</p> <p>Field strengths from fixed RF transmitters, as determined by an electromagnetic site survey,^a should be less than the compliance level in each frequency range.^b</p> <p>Interference may occur in the vicinity of equipment marked with the following symbol:</p> 
----------------------------------------------------------------------------	-------------------------	-------------------------	-----------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------

^a Field strengths from fixed transmitters, such as base stations for radio (cellular/cordless) telephones and land mobile radios, amateur radio, AM and FM radio broadcast and TV broadcast cannot be predicted theoretically with accuracy. To assess the electromagnetic environment due to fixed RF transmitters, an electromagnetic site survey should be considered. If the measured field strength in the location in which the patient lift is used exceeds the applicable RF compliance level above, the patient lift should be observed to verify normal operation. If abnormal performance is observed, additional measures may be necessary, such as reorienting or relocating the patient lift.

^b Over the frequency range 150 kHz to 80 MHz, field strengths should be less than [V1] V/m.

At 80 MHz and 800 MHz, the higher frequency range applies.

Recommended separation distances between portable and mobile RF communications equipment and the patient lift

The patient lift is intended for use in an electromagnetic environment in which radiated RF disturbances are controlled. The customer or the user of the patient lift can help prevent electromagnetic interference by maintaining a minimum distance between portable and mobile RF communications equipment (transmitters) and the patient lift as recommended below, according to the maximum output power of the communications equipment

Rated maximum output of transmitter [W]	Separation distance according to frequency of transmitter [m]		
	150 kHz to 80 MHz $d = \left[\frac{3,5}{V_1}\right]\sqrt{P}$	80 MHz to 800 MHz $d = \left[\frac{3,5}{E_1}\right]\sqrt{P}$	800 MHz to 2,5 GHz $d = \left[\frac{7}{E_1}\right]\sqrt{P}$
0.01	0.12	0.12	0.23
0.1	0.37	0.37	0.74
1	1.17	1.17	2.33
10	3.69	3.69	7.38
100	11.67	11.67	23.33

For transmitters rated at a maximum output power not listed above the recommended separation, distance d in metres (m) can be estimated using the equation applicable to the frequency of the transmitter, where P is the maximum output power rating of the transmitter in watts (W) according to the transmitter manufacturer.

At 80 MHz and 800 MHz, the separation distance for the higher frequency range applies.

These guidelines may not apply in all situations. Electromagnetic propagation is affected by absorption and reflection from structures, objects and people.

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Invacare Poirier SAS