

# USER MANUAL

MS5750/ MS5751

Stand-on Floor Scale



## **Explanation of Graphic Symbols on Label/Packaging**

Caution, consult accompanying documents before use		X	Separate collection for waste of electrical and electronic equipment, in accordance with Directive 2002/96/EC
Manufacturer of medical device			Manufacturing year of medical device
Carefully read user manual before installation and usage, and follow instructions for use.		↟	Medical electrical equipment with Type B applied part
REF	Device catalogue number	EC REP	Authorized representative in the European Community
LOT	Manufacturer's batch or lot number	MD	Device is a medical device
SN Serial number		UDI	Unique Device Identifier
<b>( €</b> 2460		Device conforms to 93/42/EEC as amended by 2007/47/EC Medical Device Directive. Four digit number refers to Notified Body.	
		Device complies with International Organization of Legal Metrology (Class III) requirements (verified models only)	
<b>CE</b> M18 0122		Device complies with EC directives (verified models only)	
		M: Conformity label Directive 2014/31/E weighing instrument	U for non-automatic
			onformity verification was CE label was applied. (ex:
		0122: Refers to Not	cified Body for metrology

#### Copyright Notice Charder Electronic Co., Ltd.

No.103, Guozhong Rd., Dali Dist., Taichung City 41262 Taiwan

Tel: +886-4-2406 3766 Fax: +886-4-2406 5612

Website: www.chardermedical.com E-mail: info\_cec@charder.com.tw

Copyright® Charder Electronic Co., Ltd. All rights reserved. This user manual is protected by international copyright law. All content is licensed, and usage is subject to written authorization from Charder Electronic Co., Ltd. (hereinafter Charder) Charder is not liable for any damage caused by a failure to adhere to requirements stated in this manual. Charder reserves the right to correct misprints in the manual without prior notice, and modify the exterior of the device for quality purposes without customer consent.



Charder Electronic Co., Ltd. No. 103, Guozhong Rd., Dali Dist., Taichung City, 41262 Taiwan

# CONTENTS

I. Safety Notes	
A. General Information	
B. EMC Guidance and Manufacturer's Declara	
II. Installation	12
A. Assembly	
B. Inserting Batteries	
C. Using Adapter	
III. Indicator	
A. Indicator and Key Functions	
B. Display layout	
IV. Using Device	
A. Basic Operation	
B. Hold	
C. Tare	
D. Body Mass Index (BMI)	
E. Body Surface Area (BSA)	
F. Print	
V. Device Setup	19
VI. Setup USB Connection to PC	21
VII. Wireless Connection	
VIII. Troubleshooting	
IX. Product Specifications	27
X. Declaration of Conformity	28

# **⚠**I. Safety Notes

#### A. General Information

Thank you for choosing this Charder Medical device. It is designed to be easy and straightforward to operate, but if you encounter any problems not addressed in this manual, please contact your local Charder service partner.

Before beginning operation of the device, please read this user manual carefully, and keep it in a safe place for reference. It contains important instructions regarding installation, proper usage, and maintenance.

#### **Intended Use**

This device is intended to measure the weight of subjects who can stand unassisted, for diagnosis of weight-related issues by professionals.

#### **General Handling**

- Device should be placed on stable, flat, solid, non-slippery surface.
- Usage on soft surfaces (ex: carpet) may result in inaccurate results.
- Ensure all parts are properly locked and tightened before operating the device.
- Device is intended to measure one subject at a time.

#### **Safety Instructions**

- Batteries should be kept away from children. If swallowed, promptly seek medical assistance.
- Expected service life: 5 years.
- Always comply with appropriate regulations when using electrical components under increased safety requirements.
- Ensure voltage marked on power supply matches mains power supply.
- The device is intended for indoor use only.
- Observe permissible ambient temperatures for use

#### **Environmental**

All batteries contain toxic compounds; batteries should be disposed of via designated competent organizations. Batteries should not be incinerated.

#### Cleaning

Device surface should be cleaned using alcohol-based wipes. Corrosive cleansing liquids should not be used. Pressure-washers should not be used.

- Do not use large amounts of water when cleaning the device, as it may cause damage to the internal electronics.
- Always disconnect device from mains power before cleaning.

#### Maintenance

■ Device does not require routine maintenance. However, regular checking of accuracy is recommended; frequency to be determined by level of use and state of device. If results are inaccurate, please contact local distributor.

#### Warranty/Liability

- The period of warranty shall be eighteen (18) months, beginning on the date of purchase. Please retain your receipt as proof of purchase.
- No responsibility shall be accepted for damage caused through any of the following reasons: unsuitable or improper storage or use, incorrect installation or commissioning by the owner or third parties, natural wear and tear, changes or modifications, incorrect or negligent handling, chemical, electrochemical, or electrical interference.
- All maintenance, technical inspections, and repairs should be conducted by an authorized Charder service partner, using original Charder accessories and spare parts. Charder is not liable for any damages arising from improper maintenance or usage.

#### **Disposal**

This product is not to be treated as regular household waste, but should be taken to a designated collection points for electronics. Further information should be provided by local waste disposal authorities.

# **△**Warning

- Only the original adapter should be used with the device. Using an adapter other than the one provided by Charder may cause malfunction.
- Do not touch the power supply with wet hands.
- Do not crimp the power cable, and avoid sharp edges.
- Do not overload extension cables connected to the device.
- Route cables carefully, to avoid tripping.
- Keep device away from liquids.
- Do not remove the plug by yanking on the cable.
- Use only a correctly wired (100-240VAC) outlet, and do not use a multiple outlet extension cable.

- Do not under any circumstances dismantle or alter the device, as this could result in electric shock or injury as well as adversely affect the precision of measurements.
- Do not place the device in direct sunlight, or in close proximity to an intense heat source. Excessively high temperatures may damage the internal electronics.

#### **Incident Reporting**

Any serious incident that has occurred in relation to the device should be reported to the manufacturer, EU representative (if device is used in EU member state), and competent authority of user/subject's member state.

#### B. EMC Guidance and Manufacturer's Declaration

#### Guidance and manufacturer's declaration-electromagnetic emissions

The MS5750/MS5751 Stand-on Floor Scale is intended for use in the electromagnetic environment specified below. The customer or the user of the device should assure that it is used in such an environment.

Emission test	Compliance	Electromagnetic environment-guidance
RF emissions CISPR 11	Group 1	The device uses RF energy only for its internal function. Therefore, its RF emissions are very low and are not likely to cause any interference in nearby electronic equipment.
RF emissions CISPR 11	Class B	The device is suitable for use in all establishments, including domestic establishments and those directly
Harmonic emissions IEC 61000-3-2	Class A	connected to the public low-voltage power supply network that supplies buildings used for domestic
Voltage fluctuations /flicker emissions IEC 61000-3-3	Compliance	purposes.

#### Guidance and manufacturer's declaration-electromagnetic immunity

The MS5750/MS5751 Stand-on Floor Scale is intended for use in the electromagnetic environment specified below. The customer or the user of the device should assure that it is used in such an environment.

Immunity test	IEC 60601 test level	Compliance level	Electromagnetic environment-guidance
Electrostatic discharge(ESD) IEC 61000-4-2	± 8 kV contact ± 2 kV, ± 4 kV, ± 8 kV, ± 15 kV air	± 8 kV contact ± 2 kV, ± 4 kV, ± 8 kV, ± 15 kV air	Floors should be wood, concrete or ceramic tile. If floors are covered with synthetic material, the relative humidity should be at least 30%
Electrical fast transient/burst IEC 61000-4-4	± 2kV for power supply lines + 1kV for input/output lines	+ 2kV for power supply lines + 1kV for input/output lines	Mains power quality should be that of a typical commercial or hospital environment.
Surge IEC 61000-4-5	± 1kV line(s) to line(s) ± 2kV line(s) to earth	+ 1kV line(s) to line(s) + 2kV line(s) to earth	Mains power quality should be that of a typical commercial or hospital environment.
Voltage Dips, short interruptions and voltage variations on power supply input lines IEC 61000-4-11	0% UT for 0,5 cycle 0% UT for 1 cycle 70% UT(30% dip in UT) for 25 cycles 0% UT for 5 s	0% UT for 0,5 cycle 0% UT for 1 cycle 70% UT(30% dip in UT) for 25 cycles 0% UT for 5 s	Mains power quality should be that of a typical commercial or hospital environment. If the user of the device requires continued operation during power mains interruptions, it is recommended that the device be powered from an uninterruptible power supply or a battery.
Power frequency(50/60 Hz) magnetic field IEC 61000-4-8	30 A/m	30 A/m	The device power frequency magnetic fields should be at levels characteristic of a typical location in a typical commercial or hospital environment.

#### Guidance and manufacturer's declaration-electromagnetic immunity

The MS5750/MS5751 Stand-on Floor Scale is intended for use in the electromagnetic environment specified below. The customer or the user of the device should assure that is used in such an environment.

Immunity test	IEC 60601 test	Compliance	Electromagnetic
	level	level	environment-guidance
Conducted RF IEC 61000-4-6	3 Vrms 150 KHz to 80 MHz	3 Vrms 150 KHz to 80	Portable and mobile RF
120 01000 1 0	100 1112 10 00 11112	MHz	communications equipment
Radiated RF IEC	6 V in ISM bands		should be used no closer to any
61000-4-3	between 0,15 MHz and 80 MHz	6 V in ISM bands between	part of the device including
	80 % AM at 1 kHz	0,15 MHz and	cables, than the recommended separation distance calculated
		80 MHz	from the equation applicable to
	3 V/m	80 % AM at 1	the frequency of the
	80MHz to 2,7 GHz	<u>kHz</u>	transmitter.
		3 V/m	
		80MHz to 2,7	Recommended separation
		<u>GHz</u>	distance:
			$d = 1,2 \sqrt{P}$
			$d = 1,2 \sqrt{P}$ 80MHz to 800 MHz
			$d = 2.3 \sqrt{P}$ 800MHz to 2.5 GHz
			Where <b>P</b> is the maximum output
			power rating of the transmitter
			in watts (W) according to the
			transmitter manufacturer and d
			is the recommended separation
			distance in metres (m).
			Field strengths from fixed RF
			transmitters, as determined by
			an electromagnetic site survey <sup>a</sup> ,
			should be less than the
			compliance level in each
			frequency range <sup>b</sup> .
			Interference may occur in the
			vicinity of equipment marked
			with the following symbol:
			4 1
			(((:)))

NOTE1 At 80 MHz and 800 MHz, the higher frequency range applies.

NOTE2 These guidelines may not apply in all situations. Electromagnetic propagation is affected by absorption and reflection from structures, objects and people.

- a Field strengths from fixed transmitters, such as base stations for radio (cellular/cordless) telephones and land mobile radios, amateur radio, AM and FM radio broadcast and TV broadcast cannot be predicted theoretically with accuracy. To assess the electromagnetic environment due to fixed RF transmitters, an electromagnetic site survey should be considered. If the measured field strength in the location in which the device is used exceeds the applicable RF compliance level above, the device should be observed to verify normal operation. If abnormal performance is observed, additional measures may be necessary, such as re-orienting or relocating the device.
- b Over the frequency range 150 kHz to 80 MHz, field strengths should be less than 3 V/m.

# Recommended separation distance between portable and mobile RF communications equipment and the MS5750/MS5751 Stand-on Floor Scale

The device is intended for use in an electromagnetic environment in which radiated RF disturbances are controlled. The customer or the user of the device can help prevent electromagnetic interference by maintaining a minimum distance between portable and mobile RF communications equipment (transmitters) and the device as recommended below, according to the maximum output power of the communications equipment.

Rated maximum output power of	Separation distance according to frequency of transmitter m			
transmitter W	<b>150</b> kHz to 80 MHz d =1,2√P	80 MHz to 800 MHz d =1,2 $\sqrt{P}$	<b>800</b> MHz to <b>2,5</b> GHz d =2,3√ <i>P</i>	
0,01	0,12	0,12	0,23	
0,1	0,38	0,38	0,73	
1	1,2	1,2	2,3	
10	3,8	3,8	7,3	
100	12	12	23	

For transmitters rated at a maximum output power not listed above, the recommended separation distance d in metres (m) can be estimated using the equation applicable to the frequency of the transmitter, where p is the maximum output power rating of the transmitter in watts (W) according to the transmitter manufacturer.

NOTE1 At 80 MHz and 800 MHz, the separation distance for the higher frequency range applies.

NOTE2 These guidelines may not apply in all situations. Electromagnetic propagation is affected by absorption and reflection from structures, objects and people.

#### II. Installation

#### A. Assembly

Device does not require installation, and can be used once power is supplied.

#### **B. Inserting Batteries**

1. Open battery housing cover



2. Take out battery housing



3. Place batteries in compartment (ensure polarity is correct)



4. Insert battery housing




5. Close battery housing cover.



6. Turn on power to confirm that battery is correctly installed.


#### C. Using Adapter

- 1. Connect adapter to indicator before connecting to mains power supply
- 2. Disconnect adapter from mains power supply before unplugging adapter pin from indicator.



#### III. Indicator

#### A. Indicator and Key Functions




#### **Key Function**

HOLD

TARE


- 1. (UNIT): Switch between units. (MS5750 only, inactive on OIML-approved MS5751)
- 2. (Send Data): When printer is connected to the indicator, press this key to send results.
- 3. On/Off/Zero): Power button. Press and hold to turn off. Press once to zero weight.
- 4. (HOLD/BMI): Press once to Hold (determine stable weighing value used when weight is unstable). Press and hold for 3 seconds to enter Body Mass Index (BMI) calculation mode.
- 5. (TARE/BSA): Press once to Tare (deduct weight from reading after measurement). After using BMI function, press once to display Body Surface Area (BSA).

#### **B. Display layout**



**BSA**: Body Surface Area is being displayed **BMI**: Body Mass Index is being displayed **NET**: Net weight appears after tare is activated

HOLD: Weight lock function is in use



**Stable symbol**: Indicates that weight is stable.

**Negative symbol**: Weight under zero.

Zero symbol: Weight is at zero

**Low battery**: Battery needs to be charged or replaced.

### **IV.** Using Device











Unit

On/Off/Zero

Hold/BMI

Tare/BSA

#### A. Basic Operation

Switch on the device using key. (To turn off device, press and hold

key for 3 seconds) The device will automatically perform self-calibration, displaying software version.

Once "0.00 kg" appears on indicator, device is ready for measurement.

**Note**: If "0.00 kg" does not display on indicator, press key to zero the device.



Guide subject to stand on device. After the weight has stabilized, the "stable" symbol will appear on indicator.

**Note**: If subject's weight exceeds scale capacity (including tare), indicator will display "Err" prompt due to overload.

#### B. Hold

The hold function determines average weight, designed to be used if subject's weight will not stabilize (ex: an active child).

Note: if fluctuation is too severe, average weight determination will be difficult and hold may not function correctly

- 1. Switch on the device normally.
- 2. Press the key. "HOLD" will be displayed on the indicator.
- 3. Guide subject to stand on device.
- 4. After a few seconds, the average weight will be displayed on the indicator. This weight will be locked - at this point, subject can leave measurement platform.
- 5. To release the locked weight, press the key again to return to the device to normal mode

**Note**: Hold function can be activated before or after subject stands on

device. However, if subject finds it difficult to stand still, we recommend activating Hold after subject stands on device. Hold function will not function under 2 kg.

#### C. Tare


The tare function allows the user to deduct the weight of objects from the device's measurement result.

- 1. Place object that needs to be tared onto measurement platform.
- 2. Press key after stable symbol appears on indicator. Display will indicate "0.00 kg".
- 3. Guide subject (plus tared object) to stand on device. Conduct measurement.
- 4. To clear tare value, remove all objects from measurement platform, and press key.

#### D. Body Mass Index (BMI)

- 1. In normal mode, press and hold the key to enter BMI mode.
- 2. Display will show last input height. Left-most digit will flash.
- 3. Adjust height value using (increase ↑) and (decrease ↓) keys. Proceed to next digit using key. Press key to confirm.

5. Proceed to weigh subject as usual. Indicator will display weight, height, and BMI after measurement.




Category	BMI (kg/m²)	Risk of obesity-related disease
Under	< 18.5	Low
Normal	18.5-24.9	Average
Over	24.9-29.9	Slightly Increased
Obese I	30.0-34.9	Increased
Obese II	35.0-39.9	High
Obese III	> 40	Very High

(World Health Organization adult BMI standards)

#### E. Body Surface Area (BSA)

1. After calculating BMI, press key. BSA will be displayed on

indicator. Press key to return to BMI mode. Press key to return to normal weighing mode.



#### F. Print

If thermal printer is connected to indicator, results can be printed by pressing key.

### V. Device Setup

When the device is switched on, press and hold the key for 6 seconds, until the display shows the "SETUP", followed by "AOFF" (first option in setting menu).

In device setup menu:

to togale next menu option

to toggle previous menu option

to confirm selection

**Auto Power-Off**: Instruct device to shut off automatically after a certain period of time.

Auto off options: 120 sec / 180 sec / 240 sec / 300 sec / off

to toggle between time options, and Press selection.

Buzzer/Beep:

When function is turned on, beeping noise will be made when: indicator is turned on, keys are pressed, and weight is stable.

to toggle between on/off, and both key to confirm Press selection.

Hold Stop: When Hold Stop is "on", Hold will deactivate after subject leaves measurement platform.

to toggle between on/off, and selection

6 WEF

**Bluetooth (optional)**: If device has Bluetooth module installed, Bluetooth function can be turned on or off.

Press **[HOLD]** to toggle between on/off, and **[TARE]** to confirm selection.



**Wi-Fi (optional)**: If device has Wi-Fi module installed, Wi-Fi function can be turned on or off.

Press **[HOLD]** to toggle between on/off, and **[TARE]** to confirm selection.



**Wi-Fi Setting (optional)**: If device has Wi-Fi module installed, this option will appear.

Press key when appears on indicator to save all settings and return to weighing mode.

## VI. Setup USB Connection to PC

For successful connection, PC hardware connected to device must be compatible with USB 2.0 or above. Operators should select a USB cable length that is most suitable to the operating environment.


1. Charder Smart Data Manager can be used to connect the device to a PC. The software program can be downloaded from the Charder website:

[LINK URL] https://www.chardermedical.com/download.htm

2. Connect USB cable to device indicator and PC. Follow installation instructions.

#### **Program Setup**


1. After installation of Charder Smart Data Manager is complete, software will automatically search for COM port. Press [Connect]. Once connected, [Connect] button will change to [Disconnect].



#### **Conducting Measurement**

1. Input subject's first name, last name, patient ID, date of birth (DD/MM/YYYY), gender, and height (for BMI calculation) into software if needed. Press **[Clear]** to clear all input.

**NOTE**: information can also be input after weight measurement.




2. Conduct measurement. If **[Auto]** is selected, results will be transmitted from device to software automatically and displayed on the left of screen. If **[Manual]** is selected, user must press "Collect".




#### Saving & Printing Results


1. Press **[Save as]** to save measurement results as .csv file on PC. Default file name is same as user ID. (ex: 20190201.csv) To track changes and multiple measurements for the same subject, we recommend not changing the default file name.



2. Result example:




If previous results were saved in "20190201.csv", new results also need to be saved as "20190201.csv" (overwriting old file) in order to save multiple results for the same subject.



Results will be saved in chronological order of measurement.

3. Press the printer icon to print out result using a printer connected to the PC.



**NOTE**: Body Surface Area (BSA) data cannot be transferred to PC. BSA results should be read from device indicator.

#### **VII. Wireless Connection**

If the device has the wireless module installed, the indicator can transmit measurement results wirelessly. Please see Charder wireless software instructions for details.

## **VIII. Troubleshooting**

Before contacting your local Charder distributor for repair service, we recommend considering the following troubleshooting procedures:

#### **Self-inspection**

#### 1. Device will not power on

- If battery power is depleted, replace with new batteries
- If batteries are not used, check if the power adapter is plugged into the device properly. Check if power adapter is plugged into mains properly

#### 2. Indicator showing "0000" ZERO SPAN out of range

- Interference due to factors such as RF disturbance or ground vibration. Relocate device to location without interference and try again
- Unstable platform feet adjust platform feet according to bubble level indication (clockwise to retract, counter-clockwise to extend) and try again
- External objects interfering with measurement platform. Clear platform of objects and try again
- Device may not function properly on soft surfaces such as carpets or lawns. Relocate device to location with solid, stable floor
- If the steps above cannot resolve the problem, re-calibration may be required to correct weighing accuracy

#### 3. Connection failure for data transmission to PC or printer

- Ensure wires are connected correctly between indicator and PC or printer
- Ensure printer is supplied with power. Ensure PC software is set up properly as indicated in this manual

#### **Distributor support required**

If the following errors occur, we recommend contacting your local Charder distributor for repair or replacement services:

#### 1. Device will not power on

- Faulty on/off key
- Broken or damaged wires causing short circuit or faulty connection
- Safety fuse burnout
- Faulty adapter

#### 2. Indicator damage

- Possible hardware defects include: uneven brightness in LCD screen, blurred text, smeared rainbow screen, incorrect decimal display
- Unable to save or read data
- Indicator shows "ERRL" after device is switched on
- Keys not responding
- Buzzer malfunction

**Error Messages** 

Error Message	Reason	Action
LobRE	Low battery warning Voltage of battery is too low to operate device	Replace batteries, or plug in AC adapter
Err	Overload  Total load exceeds device's maximum capacity	Reduce weight on measurement platform and try again
Err5	Counting Error Signal from loadcells too high or low	Error normally caused by faulty loadcell or wiring. Please contact distributor
00000	Zero count over calibration zero range +10% while power on	Remove weight from device and try again. If error persists, please contact distributor
00000	Zero count under calibration zero range -10% while power on	Remove weight from device and try again. If error persists, please contact distributor
ErrAd	Program Error Fault with device software	Please contact distributor
	Negative weight Weight reading below -2 kg.	Press key to return to 0.0.

**IX. Product Specifications** 

1X. Product Specifications			
Model		MS5750/MS5751	
Display		DP4400	
	Capacity	300 kg x 0.1 kg	
Weight	Accuracy	± 1.5e	
Measurement	OIML	Class III (MS5751 only)	
	LCD Screen	1.4-inch LCD screen (5 digits)	
Dimensions	Overall	360(W) x 430(D) x 110(H) mm	
(Standard)	Platform	360(W) x 310(D) x 75(H) mm	
Device	Device Weight 5.5 kg		
Key Functions		Unit (non-functional on MS5751), On/Off/Zero, Send Data, Hold/BMI, Tare/BSA	
Data Transmission		USB NOTE: Device should be connected to network by qualified distributors only	
Power	Supply	6 AA batteries / Power adapter	
Operation Temperature & Humidity		5°C~35°C 15% / 85% RH	
Standard Accessories		User manual x1, Power Adapter x1, USB cable x1	
Optional Accessories Thermal Printer		Thermal Printer	

# **△**Warning

The device is only compatible with the power adapters specified below.

AMP VOLTAG E	DRAWING NO.	CE APPROVED TYPE NO. / MODEL NO.	TYP E	Adapter plug
12V 1A	AD-8095	UE24WCP1-120100SPA	US EU UK AU	90 - degree

## X. Declaration of Conformity

This product has been manufactured in accordance with the harmonized European standards, following the provisions of the below stated directives:

<b>C €</b> 2460	93/42/EEC as amended by 2007/47/EC Medical Device Directive
N 4	2014/31/EU Non-automatic Weighing Instruments Directive (MS5751 only)

Please see separate document showing on sticker of device for above CE marking.

Authorized EU Representative:





CD-IN-00017 REV 005 04/2021