DIMED®

MA LTD821 (MOD. MD300K2) 02_A_06-2015

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DIMED®

MA LTD821 (MOD. MD300K2) 02 A_06-2015

HAND TYPE PULSEOXIMETER WITH ALLARM
INSTRUCTION MAUNAL

Copyright

Our company owns all rights of this unpublished work and intends to maintain it as a confidential work. We may also seek to maintain this work as an unpublished copyright. This publication is to be used solely for the purpose of reference or operation of our software system. No part of this work can be disseminated for other purposes.

In the event of inadvertent or deliberate publication, we intend to enforce its right to this work under copyright laws as a published work. Those having access to this work may not copy, use, or disclose the information in this work unless expressly authorized by our company.

All information contained in this publication is believed to be correct. We shall not be liable for errors contained herein nor for incidental or consequential damages in connection with the furnishing, performance, or use of this material. The information this publication refers to is protected by copyrights or patents and does not convey any license under the patent rights of our company, nor the rights of others. We do not assume any liability arising out of any infringements of patents or other rights of third parties.

Content of the manual is subject to change without prior notice.

All rights reserved.

Edition: Ver3.0K2

Issue Date: September 16, 2013

Contents

1 INTROUDCTUION	1
1.1 Brief Introduction	1
1.2 Intended Use	1
1.3 Measurement Principle.....	1
1.4 Safety Information	1
1.5 Electromagnetism Interference	4
1.7 Product Features.....	4
2 GENERAL DESCRIPTION	5
2.1 Appearance	5
2.2 Power Supply	6
3 TIME AND DATE SETTINGS.....	8
4 TAKE A MEASUREMENT	9
4.1 Install the Probe	9
4.2 Measurement Screens	9
5 SETTINGS	11
5.1 System Setting	11
5.2 Alarm Setting	11
6 DATA MANAGEMENT.....	12
6.1 Data Review	12
6.2 SpO ₂ Trend.....	12
6.3 PR Trend	13
6.4 Delete all data	13
7 ALARM.....	14
8 MAINTENANCE AND REPAIR.....	15
8.1 Maintenance	15
8.2 Safety Checks	15
8.3 Calibration and Verification.....	16
8.4 Troubleshooting.....	17
8.5 Warranty and Repair	18
APPENDIX A SPECIFICATIONS.....	20
APPENDIX B DECALARATION.....	22

1 INTRODUCTION

1.1 Brief Introduction

Thank you for purchasing the MD300K2 pulse oximeter. The main functions of the device include Haemoglobin Saturation (SpO₂) and pulse rate (PR) measurements, visual and audible indication, probe off and finger out indication, data storage, transmitting and review etc. Please read this manual carefully before using the device.

Notes:

The illustrations applied in the manual may differ slightly from the actual device.

The device is designed of handheld structure, please be sure not to turn upside down when using it.

1.2 Intended Use

The MD300K2 handheld pulse oximeter is intended for continuous monitoring, spot-checking of functional pulse oxygen saturation (SpO₂) and pulse rate (PR) of single adult and pediatric patients in hospitals and home care.

1.3 Measurement Principle

Principle of the oximeter is as follows: A mathematical formula is established making use of Lambert Beer Law according to Spectrum Absorption Characteristics of Reductive hemoglobin (RHb) and Oxyhemoglobin (HbO₂) in glow and near-infrared zones. Operation principle of the instrument: Photoelectric Oxyhemoglobin Inspection Technology is adopted in accordance with Capacity Pulse Scanning and Recording Technology, so that two beams of different wavelength of lights (660nm glow and 940nm near infrared light) can be focused onto a human nail tip through a clamping finger-type sensor. A measured signal obtained by a photosensitive element, will be shown on the oximeter's display through process in electronic circuits and microprocessor.

Diagram of Operation Principle

1. Red and Infrared-ray Emission Tube
2. Red and Infrared-ray Receipt Tube


Fig.1-1

1.4 Safety Information

Conception of Warning, Caution and Note

The Warning, caution and Note at this document are special information in favor of users operation.

Warning - Information concerning something that could possibly hurt the patient or operator.

Caution - Reminds the user to pay close attention to device operation, failure of which may cause abnormal function of the instrument.

Note - Informs the user of other important information by suggestion, requirement and supplement.

Warnings!

1. Before use, carefully read the manual. Our company will assume no warranty for using this

equipment improperly.

2. Operation of the handheld pulse oximeter may be affected by the use of an electrosurgical unit (ESU).

3. Sensor malfunction may cause inaccurate data possibly resulting in patient injury or death, so pay close attention to the sensor and inspect it often.

4. Do not use the handheld pulse oximeter in an MRI or CT environment.

5. Although the pulse oximeter has alarms, it is not suggested for long time continuous monitoring.

6. Do not use the handheld pulse oximeter in an explosive atmosphere.

7. The handheld pulse oximeter is intended only as an adjunct in patient assessment. It must be used in conjunction with other methods of assessing clinical signs and symptoms.

8. Check the pulse oximeter sensor application site every 4 hours to determine the positioning of the sensor and circulation and skin sensitivity of the patient.

9. When begin to gauge associated body features please follow the doctor's counseling.

10. When link this equipment to other peripherals, make sure you are sophisticated operator to handle this device. Any peripherals should be in the light of protocol of IEC 950 and IEC 601-1-1. Any input/output device should be following the protocol of IEC 601-1-1.

11. Do not sterilize the device using autoclaving, ethylene oxide sterilizing, or immersing the device in liquid. The device is not intended for sterilization.

12. Follow local ordinances and recycling instructions regarding disposal or recycling of the device and device components, including batteries.

13. This equipment complies with IEC 60601-1-2:2007 for electromagnetic compatibility for medical electrical equipment and/or systems. However, because of the proliferation of radio-frequency transmitting equipment and other sources of electrical noise in healthcare and other environments, it is possible that high levels of such interference due to close proximity or strength of a source might disrupt the performance of this device.

14. You should operate the equipment according to the EMC information provided in the accompanying documents.

15. Portable and mobile RF communications equipment can affect medical electrical equipment.

16. This equipment should not be used adjacent to or stacked with other equipment.

17. This equipment is not intended for use during patient transport outside the healthcare facility

Rx only: "Caution: Federal law restricts this device to sale by or on the order of a physician."

Cautions:

1. The handheld pulse oximeter must be able to measure the pulse properly to obtain an accurate SpO₂ measurement. Verify that nothing is hindering the pulse measurement before

relying on the SpO₂ measurement.

2. Worn-out data cables may also cause inaccurate data, so if the data is used as a reference to treat a patient, pay special attention to data cable and check it more frequently.
3. Do not tangle the SpO₂ cable with the wires of ES (Electrosurgery) equipment.
4. Single use accessories should never be reused.
5. Only use SpO₂ sensors specified by the manufacturer. Other SpO₂ sensors may cause improper performance.
6. Unplug the sensor from the monitor before cleaning or disinfecting to prevent sensor or monitor from being damaged, and to prevent user under safety situation.
7. Alarm must be set up according to different situation of individual patient. Make sure that chime sound can be activated when alarm function begin to work.

Notes:

1. Optical cross talk can occur when two or more sensors are located in adjoining areas. It can be eliminated by covering each site with opaque material. Optical cross talk may adversely affect the accuracy of the SpO₂ readings.
2. Obstructions or dirt on the sensor's red light or detector may cause a sensor failure. Make sure there are no obstructions and the sensor is clean.
3. For routine equipment maintenance, please refer to the service procedures at the associated section as indicated in the manual.
4. As to the other concerns for attention, please carefully look through the specific chapter in this instruction.

Inaccurate measurements may be caused by:

1. Significant levels of dysfunctional hemoglobin (such as carbonyl - hemoglobin or methemoglobin);
2. Intravascular dyes such as indocyanine green or methylene blue;
3. High ambient light. Shield the sensor area if necessary;
4. Excessive patient movement;
5. High-frequency electrosurgical interference and defibrillators;
6. Venous pulsations;
7. Placement of a sensor on an extremity with a blood pressure cuff, arterial catheter, or intravascular line;
8. The patient has hypotension, severe vasoconstriction, severe anemia, or hypothermia;
9. The patient is in cardiac arrest or is in shock;
10. Fingernail polish or false fingernails;
11. Weak pulse quality (low perfusion);
12. Low hemoglobin;

1.5 Electromagnetism Interference

This oximeter is designed and tested in compliance with the EMC standard, complying with the international standard for the EMC of the electronic medical device – IEC 60601-1-2. However, because of the proliferation of radio frequency transmitting equipment and other sources of electrical noise in the health-care and home environments (e.g. cellular phones, mobile two-way radios, electrical appliances) it is possible that high levels of such interference due to close proximity or strength of a source, may result in disruption of performance of this device.

This apparatus complies with the IEC 60601-1-2 international standard. The requirements of this international standard are: CISPR11, GROP1, and CLASS B.

1.6 Explanation of Symbols

Symbol	Explanation	Symbol	Explanation
	Type BF applied part	IPX1	Protected against dripping water
	Caution		Prevent from rain
	Date of Manufacture		Manufacturer's information
	Authorised representative in the European community		European union approval
SN	Serial number	ID	User ID
	Audio alarm inhibition		Battery power indication
	Storage temperature and relative humidity		Follow instructions for use

1.7 Product Features

- Compact, light in weight and simple to use.
- TFT displaying screen with adjustable backlight.
- Up to 127 user ID and 72-hour records storage.
- Visual and three-level audible alarms, low battery alarm.
- Data transfer to PC by USB cable.
- Powered by three AA alkaline / Ni-MH batteries.
- Battery charge function with adapter or charge stand (optional).

2 GENERAL DESCRIPTION

The handheld pulse oximeter adopts 2.4 inch TFT displaying screen. It can display the SpO₂, PR value, pulse bar, SpO₂ waveform and other kinds of parameters, such as the time, ID number, battery status and so on.

2.1 Appearance


Fig.2.1

Description of Fig.2-1:

1. Connector for probe and USB cable: Used to connect with SpO₂ probe and connect USB cable with PC for data transmission.
2. Alarm indicator: it will flash with yellow when alarms occur.
3. Displaying screen.
4. Navigation buttons (UP, DOWN, LEFT, RIGHT): the function of the left button is confirming the items that you choose; the function of the right button is back to the previous or menu interface. In the measuring interface, press the right button to shift between digital mode and wave mode.
5. Menu button: short press to enter the main menu; long press it to silence the audio alarm, the device will be appear the audio alarm indicator icon covered with "X" and countdown, long press it again, the device will recover the audio alarm.
6. Power button: Press and hold the button for 3 seconds to power the device on, and for about 4 seconds to power it off.

7. Connector for adapter.
8. USB interface: Used for update the software, only serves engineers;
9. Side UP/DOWN button: Short press it to increase/decrease the volume of alarm; Press it for extended time to increase/decrease the brightness of the backlight.

NOTE: In the measuring interface, press the menu button for extended time, the alarm pause will be turn on; press it again, the alarm pause will be turn off.

2.2 Power Supply

The device can be powered by 3 AA alkaline batteries, Ni-MH batteries or power adapter.

Rear panel introduction:

- 1: Battery compartment
- 2: Fixing hole
- 3: Speaker for alarm


Fig.2.2 Rear Panel

2.2.1 Powered by batteries

Batteries Installation:

- 1) **Open the battery cover:** Rotate the fixing screw slightly in the rear panel to the position which is marked with “” and then open the cover.
- 2) **Batteries installation:** Install 3 AA alkaline batteries or rechargeable batteries lightly as indicated by the polarity signs as shown in Fig.2-2.
- 3) **Close the battery cover:** Close the battery cover and rotate the screw to the  position. It means that the battery cover is locked. You can open it with attached tool.

Battery life and replacement

When the battery voltage is lower than $3.7 \pm 0.1V$, the battery indicator will become empty and its frame shows in red, which means few of battery capacity remains. You should replace the batteries with new ones in time. When the battery voltage is lower than $3.5 \pm 0.1V$ the unit will shut down.

 **Warnings!**

- If battery fluid should get into your eyes, immediately rinse with plenty of clean water. Consult a doctor immediately.
- Be sure to insert the batteries in the correct polarities, as indicated by polarity markings inside the battery compartment.

Cautions!

- Do not use batteries not specified for this unit.
- Do not dispose of batteries in fire.
- Remove the batteries from this unit when you are not going to use it for a long period.
- Do not use batteries of different types together.
- Do not use new and used batteries together.
- Dispose of batteries in accordance with the local ordinances and regulations.

2.2.2 Charging Batteries through adapter (optional)

To charge, please ensure the installed batteries are Ni-MH ones. Firstly connect the adapter with the oximeter and secondly with the outlet, and then power the device on. A prompt window appears enquiring you "Are the batteries Ni-MH ones?". Select "Yes" to charge them by outlet or "No" to abandon charge.

2.2.3 Charging Batteries by charger stand (optional)

The device can also be charged by charger stand optional with the pulse oximeter. For more information, refer to "Charger Stand Instruction Manual".

 **Warning!**

DO NOT USE Ni-MH rechargeable and alkaline batteries together, otherwise, that may damage the device or injury users.

Notes:

1. The "Brightness" cannot be set during the charging process.
2. It is not recommended to take a measurement during the charging process, for the damages even injuries may be caused to the device or users.
3. Do not charge when the non-rechargeable batteries are installed, otherwise the damages even injuries may be caused to the devices or users.

3 TIME AND DATE SETTINGS

Always set the date and time before using the unit for the first time. The date and time are important indicators when a measurement is taken.

Set the correct time according to the following steps:

1) Power on the pulse oximeter and then press the menu button to enter the main menu, refer to Fig.3-1.

Fig.3.1

Fig.3.2

2) Select the “Data and Time” item in the main menu, and then you will enter its corresponding sub-menu as shown in Fig.3.2.

3) Press the Up or Down button to select it and then press the OK button to confirm it. At last, press the Up or Down button adjust the value, and then press the OK button to confirm the value.

The date is displayed as the order of Year-Month-Day and Time of Hour-Minute

4 TAKE A MEASUREMENT

4.1 Install the Probe

1. Before use, check the pulse oximeter for mechanical damage.


Fig.4.1

2. Insert the SpO₂ probe into the probe socket (refer to Fig.4.1). If the SpO₂ probe is disconnected from the unit, a prompt "Probe Off" will be shown in the status bar in the top left of the display screen.

Note: The connector is also applied to uploading data to MedView software for reviewing. For detailed information, refer to the instruction manual of MedView software.

After finishing the time setting, plug your finger into the sensor as shown in the following illustration.


Fig.4.2 Placement of the sensor

Select the suitable sensor in terms of type and dimension.

4.2 Measurement Screens

There are two display modes as shown in the following pictures.

Fig.4.3 (Digital Mode)


Fig.4.4 (Wave Mode)

Description of measurement screens:

1. Measuring: The pulse oximeter is in the status of measuring.

It shows "Finger off" when there is no finger inserted or no signal is detected.

It shows "Probe off" when the sensor/probe is not connected to the pulse oximeter.

2. SpO₂%: SpO₂ area of display

It shows the oxygen saturation level of functional hemoglobin during normal measurement.

The color of the SpO₂ value will become red when the SpO₂ is beyond the alarm limits.

It shows two dashes throughout probe off and finger out conditions.

3. 100: SpO₂ high alarm limit; 90: SpO₂ low alarm limit.

4. : Pulse bar

5. 100: PR high alarm limit; 060: PR low alarm limit.

6. PR: PR area of display

It shows the pulse rate in beats per minute during normal measurement.

The color of the PR value will become red when the PR is beyond the alarm limits.

It shows three dashes throughout probe off and finger out conditions.

7. ID: 001, the ID number of the current patient is 001.

8. 12:00: The current time.

 **Warnings!**

● Use only SpO₂ probes provided by manufacturer. Other SpO₂ probes may cause improper performance.

● Do not use an SpO₂ probe with exposed optical components.

● Excessive patient movement may cause inaccurate measurements.

● Tissue damage can be caused by incorrect operation or misusing probe; for example, by wrapping the probe too tight. Inspect the probe site to ensure the skin's integrity and the adhesion position of the probe is correct. More frequent inspection should be taken if necessary.

● Loss of pulse signal can occur in any of the following situations:

a) The probe is too tight;

b) There is excessive illumination from light sources such as a surgical lamp, a bilirubin lamp, or sunlight; Pulse sensor should obviate the light source, e.g. radial lamp or infrared lamp.

c) A blood pressure cuff is inflated on the same extremity as the one to which an SpO₂ sensor is attached.

5 SETTINGS

5.1 System Setting

Pick and enter the [System Setting] interface from the main menu. And then press the Navigation buttons to select different item to set.


Fig.5.1

Alarm Volume: you can adjust the value of alarm volume, there are 7 levels, and the default level is 3

Alarm Pause: .there are three modes,60s, 120s, 180s, and the default mode is 120s.

Beep Tone: the level is from 0 to 7, and the default level is 3.

Backlight Setting: the level of brightness is from 1 to 7, and the default level is 3.

Language: English, French, Italian, Deutsch, Japanese and Espanol

5.2 Alarm Setting

Pick and enter the [Alarm Setting] interface from the main menu.


Fig.5.2

You can select the Alarm on or off.

High limit SpO₂ range is71~100

Low Limit SpO₂ range is 70~99

High Limit PR range is 31~235

Low Limit PR₂ range is30~234

6 DATA MANAGEMENT

Pick and enter the [Data Management] interface from the main menu. You can choose continuous measurement or not.

Fig.6.1

6.1 Data Review

Pick and enter the “Data review” interface as shown in Fig.6.1. By pressing the UP / DOWN button, the user can turn to the earlier/later item page.

The pulse oximeter can record the alarming parameter marked with red color. Press the MENU button, it will return to the main menu interface.

Time	SpO2	PR	ID
23/04 06:00:20	98	70	1
23/04 06:00:16	98	70	1
23/04 06:00:12	98	70	1
23/04 06:00:08	98	70	1
23/04 06:00:04	90	60	1
23/04 06:00:00	90	60	1
23/04 05:59:56	90	60	1
23/04 05:59:52	90	60	1
23/04 05:59:48	90	60	1
23/04 05:59:44	90	60	1

Fig.6.2

6.2 SpO₂ Trend

Pick and enter the “SpO₂ Trend” interface as shown in Fig.6.3. Press the Left / Right button; you can turn to the earlier or later item page. Press the Up button, the pulse oximeter returns to the previous interface. Press the MENU button, the pulse oximeter returns to the main menu interface.


Fig.6.3

On the above of the trend, the date and time of the first item are displayed, with year-month-day; hour: minute: second.

6.3 PR Trend

Pick and enter the "PR Trend" interface as shown in Fig.6.4. Press the Left / Right button; you can turn to the earlier or later item page. Press the Up button, the pulse oximeter returns to the previous interface. Press the MENU button, the pulse oximeter returns to the main menu interface.

Press the MENU button, the pulse oximeter returns to the main menu interface.


Fig.6.4

6.4 Delete all data

Pick and enter the "Delete all data" interface as shown in Fig.6.5. You can select "Yes" or "No" by pressing the up or down button, and by pressing the OK button to confirm your selection.

Note: Please take caution to the deletion of data; you will never get the data back once deleted.


Fig.6.5

7 ALARM

ALARM PRIORITY:

There are three-level priorities for selection.

High priority: the highest level alarm, indicates the patient is in the very dangerous situation.

Medium priority: indicates the warning-attention should be paid.

Low priority: indicate the technical alarm caused by the device itself.

Alarms of the oximeter include technical and physiological alarms. All the three priorities divided by built-in module and cannot be changed by user.

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Notes:

1. The alarm sound will go on until alarm disappears or is turned off.
2. After silencing the alarm, the corresponding indicator will indicate this.
3. The power low alarm: the corresponding indication lamp will be flashing with a red frame.

Warnings!

- When the alarm occurs, check patients' conditions immediately.
- Check which parameter is alarming or which alarm is going on.
- Check patient's condition.
- Search for the source of alarm.
- Make the alarm mute if necessary.
- Check the alarm when no warning.

8 MAINTENANCE AND REPAIR

Warnings!

- The advanced circuit inside the oximeter does not require periodic calibration and maintenance, except replacing the batteries.
- Don't open the cover of oximeter or repair electronic circuits. Its open will cause the damage of the device and the annulment of the guarantee.

8.1 Maintenance

Use only the substances approved by us and methods listed in this chapter to clean or disinfect your equipment. Warranty does not cover damage caused by unapproved substances or methods.

We make no claims regarding the efficacy of the listed chemicals or methods as a means for controlling infection. For the method to control infection, consult your hospital's Infection Control Officer or Epidemiologist. Keep your equipment and accessories free of dust and dirt. To avoid damage to the equipment, follow these rules:

- Always dilute according the manufacturer's instructions or use lowest possible concentration.
- Do not immerse part of the equipment into liquid.
- Do not pour liquid onto the equipment or accessories.
- Do not allow liquid to enter the case.
- Never use abrasive materials (such as steel wool or silver polish), or erosive cleaners (such as acetone or acetone-based cleaners).

Caution: If you spill liquid on the equipment or accessories, contact us or your service personnel.

Note: To clean or disinfect reusable accessories, refer to the instructions delivered with the accessories.

8.2 Safety Checks

Before every use, or after your pulse oximeter has been used for 6 to 12 months, or whenever your pulse oximeter is repaired or upgraded, a thorough inspection should be performed by qualified service personnel to ensure the reliability. Follow these guidelines when inspecting the equipment:

- Make sure that the environment and power supply meet the requirements.
- Inspect the equipment and its accessories for mechanical damage.
- Make sure that only specified accessories are applied.
- Inspect if the alarm system functions correctly.
- Make sure that the batteries meet the performance requirements.

- Make sure that the pulse oximeter is in good working condition.

In case of any damage or abnormality, do not use the pulse oximeter. Contact your hospital's biomedical engineers or your service personnel immediately.

Cleaning

Your equipment should be cleaned on a regular basis. If there is heavy pollution or lots of dust and sand in your place, the equipment should be cleaned more frequently. Before cleaning the equipment, consult your hospital's regulations for cleaning the equipment.

- Recommended cleaning agents are:
 - Mild soap (diluted)
 - Ammonia (diluted)
 - Sodium hypochlorite bleach (diluted)
 - Hydrogen peroxide (3%)
 - Ethanol (70%)
 - Isopropanol (70%)

To clean your equipment, follow these rules:

1. Shut down the pulse oximeter and take the batteries out of the battery wharf.
2. Clean the display screen using a soft, clean cloth dampened with a glass cleaner.
3. Clean the exterior surface of the equipment using a soft cloth dampened with the cleaner.
4. Wipe off all the cleaning solution with a dry cloth after cleaning if necessary.
5. Dry your equipment in a ventilated, cool place.

Disinfecting

The applied parts touching the patients' body are required to be disinfected once after each use. The recommended disinfectants include: ethanol 70%, isopropanol 70%, glutaraldehyde-type 2% liquid disinfectants.

Disinfection may cause damage to the equipment and is therefore not recommended for this pulse oximeter unless otherwise indicated in your hospital's servicing schedule. Clean the pulse oximeter before disinfecting it.

Caution: Never use EtO or formaldehyde for disinfection.

8.3 Calibration and Verification

The performance should be checked every one year and after maintenance and repair.

Required Test Equipment: SpO₂ signal Simulator

Note: The simulator cannot be used to assess the accuracy of a pulse oximeter probe or a pulse oximeter.

8.3.1 Control Key Verification

Press Menu key, display the history data.

8.3.2 Sound Verification

- Set the oximeter sound ON.
- Can hear the simulated heart beep sound.

8.3.3 SpO₂ & Pulse Rate Measurement Value Verification

- Connect SpO₂ Probe to the SpO₂ connector on the oximeter.
- Insert the operator's finger into the finger sensor, the SpO₂ measured value of healthy person should be from 95% to 99%, and the pulse rate is the same as heart rate.
- If SpO₂ Simulator is available, verify the accuracy of Oxygen Saturation Value with probes as follows:

Oxygen Saturation	Tolerance
96%	±2%
86%	±2%
70%	±3%

8.3.4 SpO₂ & Pulse Rate Alarm Verification

- Connect SpO₂ Probe to the SpO₂ connector on the oximeter.
- Insert the operator's finger into the finger sensor, the SpO₂ measured value of healthy person should be more than 96%.
- Set the SpO₂ high limit as 90, low limit as 80.
- Verify the SpO₂ visual and auditory alarms, the background color of the SpO₂ data should be red and "dudu" voice should be heard.

8.4 Troubleshooting

Malfunction	Reason	Resloution
SpO ₂ or PR is shown unstably	<ol style="list-style-type: none"> Finger might not be plugged deep enough. Finger is trembling or patient's body is in movement status. 	<ol style="list-style-type: none"> Retry by plugging the finger. Try not to move.
The Oximeter cannot be powered on	<ol style="list-style-type: none"> Power of batteries might be inadequate or not be there at all. Batteries might be installed incorrectly. The Oximeter might be damaged. 	<ol style="list-style-type: none"> Please replace batteries Please reinstall the batteries Please contact with local customer service centre.
Indication lamps are suddenly off	<ol style="list-style-type: none"> The product is automatically powered off when no signal is detected longer than 8 seconds Lower power 	<ol style="list-style-type: none"> Normal Replace the batteries
Others		Please contact with local customer service centre.

8.5 Warranty and Repair

8.5.1 Maintenance Method

(1) Repair response time: AM9:00 to PM17:30 on Monday to Friday except legal holiday.

Repair time: AM9:00 to PM17:30 on Monday to Friday except legal holiday.

(2) Repair service: Including telephone support, field inspecting, fittings replacement.

- Telephone support: we can give guidance to customer's engineer to inspecting the instrument when you dial our service line. Professional repair engineer online provides technical support.

- Field inspecting: we will send engineers to repair the instrument if necessary. Certified engineers of our company or local repair team trained by our company provide this service.

- Fittings replacement: if necessary, we will replace the damaged fittings according to contract. The damaged fittings should be returned to us expect for special reason.

(3) Spare machine for repair: it is used to replace the damaged machine for customer using, customer should send the damaged machine to us to repair.

(4) Repair for sponsoring and contributing machine: customer should send the machine to us to repair.

(5) Updating software is free.

8.5.2 Exempt and limitation

a) Our company isn't responsible for such damage caused by force nature. For example: fire, thunder flash, flood, cyclone, hail, earthquake, house collapse, commotion, plane failing and traffic accident, deliberate damage, lack of fuel or water, labor and capital bother, strike and stop-working etc.

b) No-service offer

- The cost and insurance charge of disassembling, refurbishing, repackaging and conveying of the oximeter or the part of it.

- Damage or loss sustained due to inspected or repaired by other institute that is not certified.

- The damage and failure caused by user or its representative who doesn't use the device according to the operator's manual

c) The damage or lose sustained due to connection to peripheral equipment (such as printer, computer etc.), that are not provided by our company are not covered by the warranty.

d) Responsibility limitation

In the duration of warranty, if user changes the parts manufactured by other manufacturers without our company permission, our company is entitled to cancel contract.

8.5.3 User Guarantees

a) Please read user manual carefully before operation

b) Please operate and make daily maintenance as request of manual and guarantee

c) Power supply and environment must be maintained under manual specifications.

8.5.4 No-guarantee principle

- The device does not remain in original condition.
- The shell of the device is breached or cracked.
- Evidence of water damage.
- Accessories adulterated or appearance of physical abuse.
- Evidence of crushing damage to the probe.
- Original Packaging during transportation is not used.
- Non authorized service is performed on oximeter.
- Damage to a product as a result of not conforming to manual specifications.
- The work environment is not eligible.
- There is smear or marks that are not belong to the instrument and cannot be removed from the outside surface of the instrument.
- The circuit is short and damaged due to liquor or other stuff flow in the instrument or its fittings.
- All probe and its accessories are not free replacement.
- If any code label of parts is damaged or missing, this warranty shall become null and void.
- For example of code label.

Such damage of probe caused by mechanical force doesn't belong to free change range.

- During measurement of SpO₂, principle leads to measuring value difficultly or inaccurate measurement.
- Maintenance seal of oximeter are not opened.

8.5.5 User's Special Request for Guarantee Time

As we stipulate the warranty period according to the relevant electronic regulation of country, which we stipulate is one year, accessory is three months. When customer requires to extend the warranty period, you should consider whether it is reasonable. Because electronic product is quickly replaced, as to the warranty period over three years, purchased accessories may be out of stock. In this case, we will adopt to entirely upgrade or replace the old, you should pay the minimum acceptable cost of renewed device.

8.5.6 Repackage

- Take all the accessories and put them into plastic cover
- Try to use original package and packing material. User will be responsible for such damage caused by bad package during transportation.
- Please offer guarantee list and copy of invoice to standby with the period of guarantee.
- Please describe failure phenomenon in detail and altogether offer oximeter.

Storage and Transportation

Storage: Storage Temperature -20°C~55°C, Relative Humidity ≤93%

Transportation: Transport by airline, train or vessel after packing according to request.

Package: We pack the product with the hard bag. We put the foam between the inner box and the cartoon to alleviate the shake.

APPENDIX A SPECIFICATIONS

Notes:

- Specifications may be changed without prior notice.
- The circuit diagrams, the list of components, the illustration of diagrams, and the detailed rules of calibration, are provided exclusively to professional personnel authorized by our company.

Display

Data: SpO₂, PR, pulse column, SpO₂ plethysmogram

Others: Information such as connection status of probe and power low indication.

Alarm

Alarm: SpO₂ and pulse rate value, probe off, battery exhausted

Alarm mode: audio alarm, visual alarm and information

Alarm limits range: SpO₂ 70%-100%, PR 30-234bpm

Default limits: SpO₂ High 100%, low 95%; PR High 100bpm; low 60bpm

SpO₂

Display range: 0%~100%

Measurement range: 70%~100%

Resolution: 1%

Accuracy: $\pm 2\%$ (80-100%); $\pm 3\%$ (70-79%); Unspecified($\leq 70\%$)

Pulse Rate

Display range: 0~254 bpm

Measurement range: 30~235 bpm

Resolution: 1 bpm

Accuracy: ± 2 bpm or 2%(The larger are applied)

Operation Environment

Operating temperature: 5℃~40℃

Relative humidity: $\leq 80\%$, no condensation

Atmosphere pressure: 86 kPa~106 kPa

Transport and Storage Environment

Storage Temperature: -20~55℃

Relative humidity: $\leq 93\%$, no condensation

Atmosphere pressure: 50~106 kPa

Power Requirement

Power supply: Three AA alkaline or Ni-MH batteries, or power adapter

Working time: work for 24 hours continuously at most

Adapter

Input Voltage: AC 100~240V

Input Frequency: 47~63Hz

Input Current: 0.7A MAX, at 100Vac, full load;

0.4A MAX, at 240Vac, full load.

Output Voltage: DC 5V±5%

Output Current: 2A MAX

Data store and replay

Store and replay 72 hours SpO₂ and PR value, the time interval is 4 seconds.

Outline of product

Dimension: 125mmX60mmX30mm

Weight: 141g(excluding the batteries)

Accessories**Standard accessories:**

1. Three AA alkaline batteries.
2. One instruction manual
3. One adult finger probe:
Model: M-50E.
4. One key for screwing the screw.
5. One piece of software CD
6. One USB cable

Optional accessories:

1. One binding sensor for pediatric: (Pediatric 15-45 Kg)
Model: M-50C
2. One fingertip sensor for pediatric:
Model: M-50B
3. One soft sensor for pediatric (Pediatric 15-45 Kg).
Model: M-50H
4. One soft sensor for adult
Model: M-50G
5. Three Ni-MH rechargeable batteries
6. One power adapter model: FRM015-S05-I
7. One charger stand (with adapter and USB cable)

Equipment Classification

Classification accordind to IEC-60601-1:

Accordind to the type of protection against Electrical shock:

Internal electrical power source equipment and class II equipment;

According to the degree of protection against Electrical shock:

Type BF equipment;

According to the degree of protection against harmful ingress of water:

IPX1;

According to the methods of sterilization or disinfection:

Non-sterilization, use of liquid surface disinfectants only;

According to the mode of operation:

Continuous operation;

APPENDIX B DECLARATION

Guidance and manufacturer's declaration - Electromagnetic emission---- for all EQUIPMENT AND SYSTEM

1	Guidance and manufacturer's declaration - electromagnetic emission		
2	The model MD300K2 Pulse Oximeter is intended for use in the electromagnetic specified below. The customer or the user of the model MD300K2 Pulse Oximeter should assure that it is such an environment.		
3	Emissions test	Compliance	Electromagnetic environment-guidance
4	RF emissions CISPR11	Group 1	The model MD300K2 Pulse Oximeter uses RF energy only for its internal function. Therefore, its RF emissions are very low and are not likely to cause any interference in nearby electronic equipment.
5	RF emissions CISPR11	Group B	
6	Harmonic emissions IEC 61000-3-2	Class A	
7	Voltage fluctuations/IEC 61000-3-3	Complies	

**Guidance and manufacturer's declaration - Electromagnetic Immunity-
For all Equipment and Systems**

Guidance and manufacturer's declaration- electromagnetic immunity			
The model MD300K2 Pulse Oximeter is intended for use in the electromagnetic environment specified below. The customer or the user of the model MD300K2 Pulse Oximeter should assure that it is used in such an environment.			
Immunity test	IEC 60601 test level	Compliance level	Electromagnetic environment-guidance
Electrostatic discharge (ESD) IEC 61000-4-4	±6kV contact ±8kV air	±6kV contact ±8kV air	Floors should be wood, concrete or ceramic tile. If floors are covered with synthetic material, the relative humidity should be at least 30%.
Electrostatic transient/burst IEC 61000-4-4	±2kV for power supply lines ±1kV for input/output lines	±2kV for power supply lines ±1kV for input/output lines	Mains power quality should be that of a typical commercial or hospital environment.
Surge IEC 61000-4-5	±1kV differential mode ±2kV common mode	±1kV differential mode ±2kV common mode	Mains power quality should be that of a typical commercial or hospital environment.
Voltage dips, short interruptions and voltage variations on power supply input lines IEC 61000-4-1	<5% UT (>95%dip in UT) for 0.5 cycles 40% UT (60%dip in UT) for 5 cycles 70% UT (30%dip in UT) For 2.5 cycles <5% UT (>95%dip in UT) for 5 sec	<5% UT (>95%dip in UT) for 0.5 cycles 40% UT (60%dip in UT) for 5 cycles 70% UT (30%dip in UT) For 2.5 cycles <5% UT (>95%dip in UT) for 5 sec	Mains power quality should be that of a typical commercial or hospital environment. If the user of the model MD300K2 Pulse Oximeter requires continued operation during power main interruptions, it is recommended that the model MD300K2 Pulse Oximeter be powered from an uninterruptible power supply or a battery.
Power frequency (50/60Hz) magnetic field IEC 61000-4-8	3A/m	3A/m	Power frequency magnetic fields should be at levels characteristic of a typical location in a typical commercial or hospital environment.
NOTE: UT is the a.c. mains voltage prior to application of the test level.			

**Guidance and manufacturer's declaration- electromagnetic immunity-
For EQUIPMENT and SYSTEM that are not LIFE-SUPPORTING**

Guidance and manufacturer's declaration – electromagnetic immunity

The model MD300K2 Pulse Oximeter is intended for use in the electromagnetic specified below. The customer of the user of the MD300K2 Pulse Oximeter should assure that it is used in such an environment.

Immunity test	IEC 60601 test level	Compliance level	Electromagnetic environment-guidance
Conducted RF IEC 61000-4-6	3Vrms 150kHz to 80kHz	3V	<p>Portable and mobile RF communications equipment should be used no closer to any part of the model MD300K2 Pulse Oximeter, including cables, than the recommended separation distance calculated from the equation applicable to the frequency of the transmitter.</p> <p>Recommended separation distance</p> $d = \frac{3.5}{r_1} \sqrt{P} \quad 150\text{kHz to } 80 \text{ MHz}$ $d = \frac{3.5}{r_1} \sqrt{P} \quad 80\text{MHz to } 800\text{MHz}$ $d = \frac{7}{r_1} \sqrt{P} \quad 800\text{MHz to } 2.5\text{GHz}$ <p>Where P is the maximum output power rating of the transmitter in Watts (W) according to the transmitter manufacture and d is the recommended separation distance in meters (m).</p> <p>Field strengths from fixed RF transmitters, as determined by an electromagnetic site survey, a should be less than the compliance level in each frequency range. b</p> <p>Interference may occur in the vicinity of equipment marked with the following symbol:</p> 
Radiated RF IEC 61000-4-3	3V/m 80MHz to 2.5GHz	3V/m	

NOTE1 At 80MHz and 800MHz, the higher frequency range applies.

NOTE2 These guidelines may not apply in all situations. Electromagnetic propagation is affected by absorption and reflection from structures, objects and people.

a Field strengths from fixed transmitters, such as base situation for radio (cellular/cordless) telephones and land/mobile radios, amateur radio, AM and FM radio broadcast and TV broadcast cannot be predicted theoretically with accuracy. To assess the electromagnetic environment due to fixed RF transmitters, an electromagnetic site survey should be considered. If the measured field strength in the location in which the model MD300K2 Pulse Oximeter is used exceeds the applicable RF compliance level above, the Pulse Oximeter should be observed to verify normal operation. If abnormal performance is observed, the additional measures may be necessary, such as reorienting or relocating the model MD300K2 Pulse Oximeter.

b Over the frequency range 150kHz to 80MHz, field strengths should be less than 3V/m.

Distribuito da / Distributed by:
Moretti S.p.A. Cavriglia (AR)
www.morettispa.com

 **Beijing Choice Electronic Technology Co.,Ltd.**

Room 4104, No. A12 Yuquan Road Haidian District,
100143 Beijing, P.R.China

EC REP Shanghai International Holding
Corp.GmbH(Europe)

Eiffestraße 80, 20537
Hamburg GERMANY

Via Bruxelles, 3 - Meleto
52022 Cavriglia (Arezzo)
Tel. +39 055 96 21 11

Fax. +39 055 96 21 200
www.morettispa.com
info@morettispa.com