MIO-VEIN

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MNPG87 Rev. 4 del 01/07/19

Elettroterapia modello

MIO-VEIN

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Electro therapy model

MIO-VEIN

INDEX	XXVII
TECHNICAL INFORMATION	29
MANUFACTURER	29
DECLARATION OF CONFORMITY	29
CLASSIFICATION	30
PURPOSE AND SCOPE	30
TECHNICAL FEATURES	31
DEVICE AND COMMANDS DESCRIPTION	32
LABELLING	33
<i>Package content</i>	34
HOW TO USE	35
INTRODUCTION TO THE TECHNOLOGY	35
INDICATIONS	35
CONTRAINDICATIONS	36
<i>Side effects</i>	36
WARNINGS	36
DEVICE USE	39
<i>Preliminary settings</i>	39
<i>Therapy settings</i>	40
<i>Key block</i>	42
<i>Battery installation</i>	42
DEVICE'S CARE	43
MAINTENANCE	43
TROUBLESHOOTING	44
<i>Battery charging</i>	45
DISPOSAL	46
WARRANTY	46
<i>Support</i>	48
<i>Spare part</i>	48
ELECTROMAGNETIC INTERFERENCES AND ELECTROMAGNETIC COMPATIBILITY TABLES	48

WARNING. MIO-VEIN IS A MEDICAL DEVICE

Consult your doctor before using MIO VEIN if you suffer from health problems.

Read carefully the user manual and the electrode position manual before using MIO-VEIN. Read the contraindications and warnings carefully.

Manufacturer

I.A.C.E.R. S.r.l.

Via S. Pertini, 24/a • 30030 Martellago (VE)

Tel. 041.5401356 • Fax 041.5402684

IACER Srl. is an Italian manufacturer of medical devices (certified CE no. MED24021 from the Notified Body no. 0476, Kiwa Cermet Italia Spa).

Declaration of conformity

I.A.C.E.R. S.r.l

Via S.Pertini 24/A – 30030 Martellago (Ve), Italia

herewith declares under its own responsibility, that the product

MIO-VEIN

UMDNS Code: **13762**

Batch no.:

Series no.:

has been designed and manufactured according to the European Medical Device Directive 93/4/EEC (transposed in Italy by the D.Lgs. 46/97), as modified by the Directive 2007/47/EC (D.Lgs.37/2010) and further modifications/integrations.

The product has been assigned to class IIa, according to Annex IX, rule 9 of the Directive 93/42/EEC (and further modifications/integrations) and bears the mark

Compliance of the concerned product with the Directive 93/42/EEC has been assessed and certified by the Notified Body

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Certificate no.: MED24021

following the certification procedure according to Annex II (excluding point 4) of the Directive 93/42/EEC.

A handwritten signature in black ink, appearing to read 'Massimo Marcon', is written over a horizontal line.

MASSIMO MARCON

Martellago, 14/09/2018

Place, date

Legal Representative

Classification

The MIO-VEIN has the following classification:

- class IIa (Directive 93/42/EEC, Annex IX, rule 9, 10 and further amendments);
- class II with BF type applied part (classification IEC 60601-1);
- equipment protection level IP22 against liquid and dust penetration;
- equipment and accessories not subject to sterilization;
- equipment unsuitable for use in presence of a flammable anesthetic mixture containing air, oxygen and nitrous oxide;
- equipment suitable for continuous operation;
- equipment unsuitable for outdoors use.

Purpose and scope

Clinical intended use:

Therapeutic

Environmental intended use:

Ambulatory and home

The MIO-VEIN electronic stimulator is a medical device specifically intended for domestic use and is targeted to adults, who acquired the specific knowledge to use the device by reading the present manual.

MIO-VEIN is used to apply electrical micro impulses which stimulates calf muscles and so it helps the correct venous circulation, leading also to an improvement in health conditions in patient who suffer from nervous system pathologies.

Technical features

Characteristic	Specification	
Power supply	Rechargeable batteries AAA Ni-MH 4.8V 800mAh Battery recharger: power supply line AC 100-240V, 50/60Hz, output DC 6.8V/300mAh max*	
Isolation (IEC 60601-1)	II	
Applied part (IEC 60601-1)	BF	
Protection level	IP22	
Applied part to the patient	Electrodes	
Dimensions (length x height x depth)	140x70x28mm	
Weight main body	210g including batteries	
Number output channel	2 independent and galvanically isolated	
Functioning	Continuous	
Output current	60mA at 1kΩ of load	
Impulse	Rectangular, biphasic and compensated	
Frequency	Burst impulses at 250Hz	
Impulse's width	From 30 to 240μs	
Therapy	Time depending on the program, max 20min	
Display	Reflective and illuminated LCD display	
Command	ABS keyboard with 7 keys	
Conditions of use	Environmental temperature	From +5° to +40°C
	Relative humidity	From 30% to 75%
	Atmospheric pressure	From 700 to 1060hPa
Storage and transportation conditions	Environmental temperature	From -10° to +55°C
	Relative humidity	From 10% to 90%
	Atmospheric pressure	From 700 to 1060hPa


WARNIGN. The equipment delivers current over 10mA.

* Use only the battery recharger given by the manufacturer. The use of other recharger could seriously compromise the security and safety both of the patient and of the device.

Expected useful life of the device is set in 3 years, meanwhile the expected useful life of the electrodes is set in 10/15 uses.

Device and commands description


- (1) power supply socket
- (2) channel 1 (CH1) output
- (3) channel 2 (CH2) output
- (4) led indicator CH1
- (5) LCD display
- (6) increase intensity CH1
- (7) decrease intensity CH1
- (8) pause
- (9) led indicator CH2
- (10) increase intensity CH2
- (11) decrease intensity CH2
- (12) lock/unlock key
- (13) ON/OFF
- (14) battery cover

- (1) time indicator
- (2) timer (symbol)
- (3) treatment number or treatment's phase indicator
- (4) lock/unlock (symbol)
- (5) CH1 output intensity (numeric)
- (6) CH1 output intensity (symbol)
- (7) pause (symbol)
- (8) output (symbol)
- (9) battery (symbol)
- (10) CH2 output intensity (numeric)
- (11) CH2 output intensity (symbol)

Labelling

Symbol	Description
	Manufacturer's logo.

Symbol	Description
	Product CE certification released by Notified Body n°0476, Kiwa Cermet.
	Manufacturer
	Manufacturing date (YYYY-MM)
	Read instructions for use
	The product must be disposed as “electronic waste”, in accordance to WEEE Directive on waste electrical and electronic equipment.
	Applied part type BF
	Power supply
IP22	Medical device protected against the penetration of solids (with a diameter $d \geq 12,5mm$) and against the vertical drops when the device is kept at 15° from its normal functioning position.
	Limits of relative humidity (relative humidity of the storage environment, on the package).
	Temperature humidity (temperature of the storage environment, on the package).

Package content

The MIO-VEIN box contains:

- n° 1 mainframe;
- n° 1 batteries pack (in the device);
- n° 2 cable for the electrostimulation (yellow e blue);
- n° 2 set of 2 oval electrodes;
- n° 1 bag for the transportation;
- n° 1 user manual;
- n° 1 user manual of the electrodes' positions.


Introduction to the technology

The electrostimulation consists on the transmission of electric micro impulses to the human body and it can be applied in different fields: pain therapy, muscles stimulation and beauty treatment. Specific electric impulses are used for every one of these applications.

For this reason, MIO-VEIN is an electrotherapy device for deep stimulation of calf muscles. Its stimulation reactivates calf muscular pump that is mainly responsible of venous return.

Calf muscles contraction promotes the compression of superficial and deep veins. Therefore, blood reaches easily the heart.

Gait helps venous return, thanks to muscle pump that moves bloodstream back to the heart. This happens in normal conditions when veins have healthy valves and walls.

As veins walls and valves become weak, blood flow slows down and stagnate (venous stasis). This causes the inversion of venous circulation.

Indications

MIO-VEIN is suggested for all people suffering from venous diseases and also in all situations that represent potential risk factors for venous problems.

The main clinical applications are:

- venous stasis;
- lower limbs edemas;
- varicose veins;
- post-phlebitis syndrome;
- post-thrombotic syndrome;
- pregnancy related venous disorders;
- venous ulcers;
- leg itch after compressive stocking use.

MIO-VEIN helps the correct venous circulation and consequently it promotes healthy and beautiful legs.

The treatment effectiveness depends mostly on the physician choice of therapy, in fact the stimulation intensity and the electrodes positioning should be based on the doctor's prescription.

Contraindications

The device must not be used in presence of cancerous injuries in the area to be treated. The stimulation should not be applied to infected, swollen or inflamed areas and in case of rashes (phlebitis, thrombophlebitis, etc.). The use of the device is not suited for pregnant or in labor women.

It is absolutely forbidden to use MIO-VEIN if the patient suffers from severe arrhythmia or has a pacemaker, suffers from epilepsy, heart-related pathologies, phlebitis in place, thrombophlebitis, in feverish state, anxiety or serious illnesses, abdominal or inguinal hernias or in case of severe injuries in the area to be treated (except in case of medical prescription).

Do not use the device if the source of the pain is unknown or not diagnosed. Use the device **ONLY** after having a diagnosis. In the event of injury, muscle stress or any other health problem consult your doctor before using the device and only use it under medical supervision.

Side effects

No significant side effects are known. In some cases of particularly sensitive people, skin redness occurs at the electrodes after treatment: the redness normally disappears few minutes after treatment. If the redness persists, consult a doctor.

In some rare cases evening stimulation causes some difficulties in falling asleep. In this case, suspend the treatment and consult a doctor.

The inadequate dimensions (too small) of the electrodes in use could provoke skin reactions or bruises.

Warnings

Treatment efficacy depends on the patients' selection by qualified persons. ***The long-term effects of a chronic stimulation are unknown.***

MIO-VEIN has been designed and manufactured to be operated exclusively with the internal, rechargeable batteries.

It is recommended:

- to control position and meaning of all the labels on the equipment;
- not to damage the connection cables to the electrodes and to avoid winding the cables around the device;

- to avoid the use of the device by persons who did not read carefully this manual;
- not to wear metal objects during treatment;
- to consult the physician before using MIO-VEIN, if recently been under surgery: the electrostimulation could slow the healing process;
- to use the electrodes on clean and dry skin. When using the electrodes, follow the instructions given in the manual and on the package of the electrodes. Use only single-patient electrodes, supplied exclusively by the manufacturer, and take care to avoid the exchange of electrodes between different users. MIO-VEIN has been tested and guaranteed for the use only with the electrodes supplied by the manufacturer;
- to use ONLY accessories supplied by device manufacturer.

It is forbidden:

- to use the device in the presence of patient monitoring equipment, of electrosurgical or shortwave or microwave therapy equipment or other equipment that sends electrical impulses into the body and in general in combination to other medical devices;
- to use the device by persons known to be unsound-minded, or suffering from sensibility disorders, permanently or temporarily disabled unless assisted by qualified personnel (e.g. a doctor or therapist); by persons younger than 15 years old or not adequately educated about the device use by an adult person;
- to use the device in presence of signs of deterioration of the device itself, cables and/or electrodes: please contact the dealer or the manufacturer following the instructions given in the paragraph *Support*. Control carefully the integrity of the device before each use;
- to use the device close to flammable substances/gas/explosives, in environments with high concentrations of oxygen, with aerosol-therapy devices or in wet environments (use of the device is prohibited in bathroom or shower areas or while showering/bathing);
- to use the device while driving a vehicle or controlling machineries/equipment;
- to position the electrodes in such a way that the current crosses the heart area (e.g. a black electrode on the chest and a red electrode on the shoulder blade) or the head (i.e. do not position the electrodes on the opposite sides of the head);
- to use the device with electrodes positioned on or close to injuries or cuts;
- to position the electrodes on the carotid sinuses (carotid) or genitals;

- to position the electrodes close to the eyes; make sure that the current delivered does not cross the eyeball (one electrode diametrically opposite to the other in relation to the eye); keep a distance of at least 3cm from the eyeball;
- to use electrodes when damaged, even if they stick to the skin well. A repeated utilization of the same electrodes could compromise the safety of the stimulation, in fact it can cause persisting skin redness also after few hours after stimulation;
- to use pointed or sharp objects on the device keyboard.

Warnings:

- do not use the device while sleeping;
- insufficiently sized electrode sections can cause skin reactions or burns;
- do not use damaged electrodes even if they well adhere to the skin;
- pay attention to use connection cables with children/young people: strangulation danger. Take away from the reach of children;
- do not mix connection cables up with earphones or other devices and do not connect the cable to other equipment;
- keep right distance between electrodes: the contact between electrodes could cause wrong stimulations or irritations/burns;
- ***stimulation intensity and electrodes position should be suggested by the prescriber doctor.***

The manufacturer considers himself responsible for the performances, reliability, safety and security of the device only if:

- any addition, modification and/or repair are carried out by authorized personnel;
- the environmental electrical installation to which MIO-VEIN is connected is compliant to the national laws;
- the instructions for use contained in this manual are strictly followed.

Should any foreign materials penetrate the device contact the retailer or manufacturer immediately. If dropped down, check that the housing is not cracked or damaged in any way; if so, contact the retailer or manufacturer.

Should you notice any changes in the device's performance during treatment, interrupt the treatment immediately and consult the retailer or manufacturer (patients being treated in a center must also be informed).

Muscle heaviness and cramps may be experienced in the event of prolonged muscle stimulation. Should this occur we recommend suspending treatment with MIO-VEIN for a few days.

If the stimulation is uncomfortable decrease intensity. If the problem persists consult a doctor.

Some patient may experience skin irritation or hypertensives due to the electrical stimulation or gel. If rash develops or pain persists, stop the treatment and consult a doctor.

Consult a doctor before using MIO-VEIN with metallic osteosynthesis devices. IF YOU HAVE ANY DOUBTS REGARDING THE DEVICE USE CONSULT YOUR DOCTOR.

Device use

Preliminary settings

ELECTRODES AND CABLE CONNECTIONS

The electrodes supplied are top quality, pre-gelled and ready for use, moreover they are particularly suitable for electrostimulation treatments.

Before positioning the electrodes, make sure that the area to be treated is clean and dry. Then remove the protective plastic from the electrodes and, in sitting or lying position, place them in the treatment area (see images below or follow the medical prescription).

Connect the electrodes to the cable and then the cables to the proper connections' door on the top of the MIO-VEIN (in correspondence to the two channels CH1 and CH2).

After the use restore the electrodes in the protective plastic and therefore in the dedicated space in the bag with the device.

Make sure that MIO-VEIN is switched off **before disconnecting the electrodes** at the end of the treatment.

Electrode duration is determined by the skin's pH value; we recommend no more than 10/15 applications using the same electrode. Repeated use of the same electrodes can compromise stimulation safety, which is why electrodes must not be used when they no longer stick to the skin; in fact, it can cause skin redness that can last for many hours after stimulation (consult a doctor if this occurs).

For a safe use please refer to indications and warnings on electrodes package and in previous *Warning* paragraph.

Therapy settings

1. Switch on MIO-VEIN by pressing ON/OFF button.
2. The display shows in the middle therapy time (20 minutes – ①), on the left side the number of treatments that have been executed (②) and on the right the battery status icon (③).

3. Check battery status (see following paragraph *Battery charging*).
4. Start the program by pressing the ▲ button of the channels CH1 and CH2 and increase the intensity of the output current; the indicator [⊗] will start blinking on the display:

5. The display will show the remaining therapy time, while it will be possible to increase the intensity of the stimulation with the ▲ button of the interested channels until the personal limit tolerance (i.e. stimulation comfort). If the intensity is too high, with the ▼ button it is possible to decrease it. The intensity should be adjusted to obtain a soft muscles contraction. The stimulation intensity is represented in the display for each channel through an increasing

scale from 0 to 60 and through a bar icon; as the intensity increases, the number of the bars increases too next to the corresponding numeric value.

- When the intensity reaches 7mA, the increasing intensity bar icon will be shown for each channel. The maximum intensity that can be set is 60mA, by step of 1mA;

- The program is divided into two phases; an acoustic signal informs the phase changing. The number of the ongoing phase is shown in the central region of the display (**01** and **02** respectively)
- An acoustic signal informs of the therapy end and the stimulator will return to the main menu, decreasing the intensity to zero.
- It is possible to pause the therapy in every moment by pressing the button : the display will show on the upper right side the blinking icon . In order to restart therapy, press again .
- To switch off MIO-VEIN, keep pressing **ON/OFF**.

WARNING! MIO-VEIN will enter into *sleeping mode*, if no key is press for longer than 2 minutes to preserve battery.

WARNING: MIO-VEIN has a safety output control to check the correct connection of cables and electrodes and also the correct transmission of the current to the patient. In case of incorrect **connection is not possible to exceed 9mA intensity: MIO-VEIN will stop providing current and return to main screening.** Check cables and electrodes status (eventually replace them) and the correct connection to the patient. If the problem persists, contact the manufacturer.

The device is provided with a function, which memorizes the treatments that are performed. The treatments' number is shown on the display before starting the session when the device is in standby. If the treatment is interrupted before the end, it will not be counted and therefore memorized.

The maximum number of treatments that can be memorized is 99: over this value the memory is reset and the counting restarts from 1.


N.B. Disconnect all the cables, before placing the stimulator in the bag. If not, the cables will bend excessively near the connectors, which can lead to damage to the cables.

Key block

The device is provided with key lock function to avoid accidental modifications to parameters during the treatment.

- Key lock is automatically set after around 20 seconds from the beginning of the therapy: the display will show the flashing icon  on the left side and the keyboard will be locked.
- To unlock the key lock function, press the button  : all the buttons will be unlocked and the lock-icon on the display will disappear.

Battery installation

Open the battery compartment on the back side of the device, connect the battery wire to the plug inside the battery compartment. Insert the battery and close the battery compartment.

WARNING: remove the battery in case of prolonged inactivity (over 2 month).

Batteries have to be handled by adult persons: keep them out of children's reach.

Maintenance

If used following the instructions given in this user guide, the equipment does not require any kind of maintenance.

It is recommended that the manufacturer carries out a functional test every 24 months.

The manufacturer does not consider the MIO-VEIN device repairable by any personnel outside the company. Each operation of the kind perpetuated by personnel not authorized by the manufacturer will be considered as tampering the device, freeing the manufacturer from granting warranty and from any danger that the user or the operator may be exposed to.

CLEANLINESS

It is recommended to switch off MIO-VEIN at the end of each treatment and to remove the cables from the connectors.

Before cleaning the device, take off ALWAYS the battery. Clean the device using only a dry soft cloth. Resistant stains can be removed using a sponge soaked in solution of water and alcohol.

Device not subject to sterilization.

Note:

- Never use solvents for cleaning. Cleaning agents cause damage to the device.
- Attention to the need for periodic maintenance, especially:
 - inspection of main body for cracks, which may allow the ingress of conductive fluid;
 - inspection of the main cable, hoses and associated connectors.

TRANSPORTATION AND STORAGE

Precaution for the transportation

There is no particular precaution to be taken during transportation of the device, since MIO-VEIN is a portable device. Anyway, it is recommended to keep MIO-VEIN in the supplied carrying bag, together with the rest of the equipment and carefully stored on a secure surface after each treatment. Protect the device from high temperature, direct daylight and liquids. Store the device in a cool, well-ventilated place.

Do not store heavy objects on the device.

Precaution for the storage

The performances of the equipment are granted if it is stored according to the following conditions:

Outside the carrying bag:

temperature:	from +5°C to +40°C
maximum relative humidity:	from 30 to 75%
atmospheric pressure:	from 700hPa to 1060hPa

Inside the supplied carrying bag:

temperature:	from -10°C to +55°C
maximum relative humidity:	from 10 to 90%
atmospheric pressure:	from 700hPa to 1060hPa

Troubleshooting

Any type of work on MIO-VEIN must be carried out exclusively by the manufacturer or by an authorized dealer. In any event, any presumed malfunction of MIO-VEIN must be verified before sending the device to the manufacturer.

Here below are proposed some typical situations:

- **MIO-VEIN cannot be turned on.** Verify the battery is correctly installed. Recharge the batteries and try to switch it on again; if it still does not switch on, check that the power supply socket is working. Check that the ON/OFF button has been pressed correctly (held down for at least one second).
- **MIO-VEIN does not transmit electric impulses.** Check that the cable jacks have been inserted in the electrodes and that the plastic protection has been removed from the electrode. Check that the cables have been connected correctly (connector well inserted in the device). Check the electrodes status: damage or poor-quality electrodes cannot properly transmit electrical impulses. Substitute the electrodes and/or the cables if necessary.
- **MIO-VEIN switches off during operation.** Recharge the batteries and start the treatment again. If the problem persists, contact the manufacturer.

WARNING: MIO-VEIN has a safety output control to check the correct connection of cables and electrodes and also the correct transmission of the current to the patient. In case of incorrect **connection is not possible to exceed 9mA intensity: MIO-VEIN will stop providing current and return to main screening.** Check cables and electrodes

status (eventually replace them) and the correct connection to the patient. If the problem persists, contact the manufacturer.

Battery charging

MIO-VEIN is supplied by internal rechargeable Ni-Mh 800 mAh battery with new long-lasting technology. In case the battery needs to be recharged or has no enough power, the display will show the following blinking icon . In this case it could be not possible to start the therapy or to complete it.

Battery fully
discharged

Battery partially
discharged

Battery
fully
charged

To proceed with the charging follow the steps below:

- make sure that the device is **switched off** before charging;
- make sure that the device is NOT being connected to the patient (disconnect cables and electrodes);
- connect the battery charger to the plug on the upper side of MIO-VEIN and connect the battery charger into the power socket;
- the display will show the flashing icon of battery status on the right side.

After 6/8 hours the recharging will end automatically, so therefore it will be possible to disconnect the charger from the power socket and from the device. Do not immerge the battery in water or other liquids and do not expose it to heat sources.

Do not dispose of dead or defective batteries with domestic waste; dispose of in an authorized waste collection bin or in any case according to the underlying norm (WEEE).

Only adults should be managing the battery. Keep out the reach of children. We suggest to fully discharge batteries before the new charging to guarantee a longer batteries duration.

WARNING: in case of prolonged inactivity it may be necessary to perform a few charge cycles to restore normal battery capacity.

It is recommended to remove the battery in case of prolonged inactivity.

The life cycle of the battery depends indeed on the charge/discharge cycles and from their number; a full recharged battery lasts approximately two weeks, when using 20minutes a day at full strength. It is suggested the following precaution for a longer battery duration:

- recharge the battery once a month, even if the device is not used;

- discharge the battery as much as possible before recharging;
- use only and exclusively the original battery charger or in any case those supplied by the manufacturer/distributor. Do not open or modify the battery charger.

WARNING: after the first recharge, the battery could present a limited capacity. This is a normal characteristic of Ni-MH battery. We suggest recharging the battery after each treatment for 3 or 4 times.

Do not use the battery charger if:

- the plug is damaged or parts of it are broken;
- it has been exposed to rain or any other type of liquid;
- the components have been damaged by a fall.

Use a dry cloth to clean the battery charger.

Do not open the battery charger: it does not contain repairable parts.

Disposal

The electro stimulator MIO-VEIN was designed and engineered to have minimal negative environmental impact, in consideration of its performance and safety requirements, following the disposition given by the European Directive 2012/19/EU, regarding the waste of electrical and electronic equipment.

Rigorous standards were followed in order to minimize the amount of waste, use of toxic materials, noise, non-required radiation and energy consumption. A deep research on the optimization of machine performances guarantees a significant consumption's reduction, in accordance to the saving energy principles.

This symbol means that the product shall not be disposed as domestic waste.

The correct disposal of obsolete equipment, accessories and most of all of batteries contributes in preventing possible negative consequences on human and environmental health.

The user must dispose of scrap equipment by taking it to a recognized center for recycling of electrical and electronic equipment.

For further information on the obsolete equipment disposal please contact the dedicated disposal service or the shop in which the device was bought.

Warranty

IACER Srl guarantees a warranty period from the purchasing date for MIO-VEIN device, unless information contained in this manual regarding installation, use

and maintenance is strictly adhered. The wearing parts (batteries and electrodes) are not included in the warranty, unless of visible manufacturing defects. The warranty is void in case of tampering of the device and in case of intervention on the same by personnel not authorized by the manufacturer or by the authorized dealer.

The warranty conditions are those described in the following paragraph Warranty conditions. As established by the European Medical Devices Directive 93/42/EEC, the manufacturer must maintain the traceability of the equipment, in order to be able to intervene promptly in case needed as a consequence of manufacturing defects.

WARNING: if not sent, the manufacturer declines any responsibility should any corrective interventions on the device be needed.

Should you need to return the goods then please pack the device and all the accessories so that it won't be damaged during transportation. In order to be entitled to the warranty assistance, the purchaser must enclose to the device a copy of the purchasing receipt, proving origin and purchasing date.

For more information on the warranty please contact the distributor or vendor, in order to check the norm and standard in force in your Country or, ultimately, the manufacturer IACER Srl.

Warranty conditions

- 1) Should assistance be needed, enclose the purchasing receipt when sending the device to the manufacturer.
- 2) The warranty period is valid only on the electronic parts. The warranty will be granted by the shop or directly by the manufacturer.
- 3) The warranty covers only the product damages, which causes its malfunctioning.
- 4) Warranty means that only the manufacturing defect components or material are covered by reparation or free substitution, hand work included.
- 5) Warranty is not applied to damages caused by negligence or use not compliant to the given instructions, by intervention on the device from personnel not authorized, accidental causes or negligence from the purchaser.
- 6) Warranty is not applied in case of damages caused by unsuitable power supplies.
- 7) Warranty does not apply to wearing parts.
- 8) Warranty does not include transportation costs which have to be covered by the purchaser.

- 9) After the warranty period the warranty is no more applicable. In this case all the assistance interventions will be performed by debiting the costs of the substitution of the parts, the hand work and the transportations costs.
- 10) The court of Venice has exclusive jurisdiction over any dispute.

Support

The manufacturer is the one and only allowed to operate with technical assistance. For any technical assistance contact:

I.A.C.E.R. S.r.l.
Via S. Pertini, 24/a • 30030 Martellago (VE)
Tel. 041.5401356 • Fax 041.5402684

Technical documentation related to repairable parts could be attached, but only with previous authorization from the manufacturer and only after giving proper training to the staff employed in technical assistance.

Spare part

The manufacturer makes available at any time the original spare parts for the equipment. Please contact:

I.A.C.E.R. S.r.l.
Via S. Pertini, 24/a • 30030 Martellago (VE)
Tel. 041.5401356 • Fax 041.5402684

In order to preserve the warranty, the functionality and the security and safety of the product, it is highly recommended to use exclusively the spare parts given by the manufacturer (see also paragraph *Warnings*).

Electromagnetic interferences and electromagnetic compatibility tables

The MIO-VEIN equipment has been designed and manufactured according to the ELECTROMAGNETIC COMPATIBILITY DIRECTIVE 2014/30/EU with the aim of providing adequate protection from harmful interference when installed in homes and health establishments.

The equipment does not generate significant radio frequency energy and is adequately immune to radiated electromagnetic fields. Therefore, it does not detrimentally interfere with radio-electric communications, electro-medical equipment for monitoring, diagnosis, therapy and surgery, office electronic devices such as computers, printers, photocopiers, fax machines, etc. or any

electric or electronic equipment used in these environments, as long as the equipment complies with the ELECTROMAGNETIC COMPATIBILITY directive.

In general, the use of accessories other than those specified or provided by the manufacturer could result in increased electromagnetic emissions or decreased electromagnetic immunity of the MIO- VEIN and result in improper functioning.

In any case, in order to avoid any interference problems, it is recommended to use the therapy equipment enough far away from critical equipment for monitoring vital patient functions, and to be careful when applying therapy to patients with pacemakers. In any case it is recommended to use the equipment at least at 3 meters away from televisions, monitors, cellphones or any other electronic equipment, in particular portable RF communications equipment (including peripherals such as antenna cables and external antennas) should not be used closer than 30cm (12 inches) to any part of the device, including the cables specified by the manufacturer; otherwise, it could lead to degradation of the performance of the MIO-VEIN .

In conclusion, the use of MIO-VEIN adjacent to or stacked with other equipment should be avoided, since it could cause improper functioning. If such use is necessary, the MIO-VEIN and the other equipment should be constantly observed to verify that they are operating normally.

When the device is used in an environment relatively dry, strong electromagnetic interferences usually occur. At this time, the device may be affected as follows:

- the device stops supplying;
- the device turns off;
- the device restarts.

The above phenomena do not affect the basic safety and essential performance of the device, which can be normally used according to the instructions given in this manual. If you want to avoid the above phenomena, please use MIO- VEIN according to the environment's conditions specified in the manual.

For more details, please see the following EMC TABLES.

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Note: U_T is the A.C. mains voltage prior to application of the test level.

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Portable and mobile RF communications equipment should not be used near any part of the MIO-VEIN (including cables) except when the recommended separation distance is respected. This distance is calculated from the equation applicable to the frequency of the transmitter.

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RF condotta <i>Conducted RF</i> IEC 61000-4-6	$3V_{\text{eff_Vrms}}$ da 150kHz a 80 MHz <i>from 150kHz to 80MHz</i> $6V_{\text{eff_Vrms}}$ da 150kHz a 80MHz bana frequenze amatoriale e per banda ISM <i>amateur radio band from 150kHz to 80MHz and ISM band</i>	$3V_{\text{eff_Vrms}}$ da 150kHz a 80 MHz <i>from 150kHz to 80MHz</i> $6V_{\text{eff_Vrms}}$ da 150kHz a 80MHz bana frequenze amatoriale e per banda ISM <i>amateur radio band from 150kHz to 80MHz and ISM band</i>	$d = 1,2 \sqrt{P}$ da 150kHz a 80MHz e per banda ISM <i>from 150kHz to 80MHz and for ISM band</i>
RF irradiate <i>Radiated RF</i> IEC 61000-4-3	10V/m da 80MHz a 2.7GHz <i>from 80MHz to 2.7GHz</i>	10V/m	$d = 1,2 \sqrt{P}$ da 80MHz a 800MHz <i>from 80MHz to 800 MHz</i> $d = 2,3 \sqrt{P}$ de 800MHz a 2,7GHz <i>from 800MHz to 2,7GHz</i>

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Where P is the maximum output power rating of the transmitter in watts (W) according to the transmitter manufacturer and d is the recommended separation distance in metres (m). Field strengths from fixed RF transmitters, as determined by an electromagnetic site survey^a, should be less than the compliance level in each frequency range^b. Interference may occur in the vicinity of equipment marked with the following symbol:

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The MIO-VEIN is intended for use in an electromagnetic environment in which radiated RF disturbances are controlled. The customer or the user of the MIO-VEIN can help prevent electromagnetic interference by maintaining a minimum distance between portable and mobile RF communications equipment (transmitters) and the MIO-VEIN as recommended below, according to the maximum output power of the communications equipment.

Potenza di uscita massima del trasmettitore specificata (W) <i>Rated maximum output power of transmitter (W)</i>	Distanza di separazione alla frequenza del trasmettitore (m) <i>Separation distance according to the frequency of the transmitter (m)</i>		
	<i>d = 1,2 √P da 150kHz a 80MHz from 150kHz to 800MHz</i>	<i>d = 1,2 √P da 80MHz a 800MHz from 80MHz to 800MHz</i>	<i>d = 2,3 √P da 800MHz a 2,7GHz from 800MHz to 2,7GHz</i>
0.01	0.12	0.12	0.23
0.1	0.38	0.38	0.73
1	1.2	1.2	2.3
10	3.8	3.8	7.3
100	12	12	23

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For transmitters rated at a maximum output power not listed above, the recommended separation distance *d* in metres (m) can be determined using the equation applicable to the frequency of the transmitter, where *P* is the maximum output power rating of the transmitter in watts (W) according to the transmitter manufacturer.

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I.A.C.E.R. Srl