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PLEASE READ THE INSTRUCTIONS CAREFULLY AND THOROUGHLY AND KEEP THEM FOR FUTURE REFERENCE

Orthosis for stabilising and supporting the lumbar spine, with malleable reinforced stays

DECLARATION OF CONFORMITY

As the manufacturer, ORTHOSERVICE AG states under its sole responsibility that this product is a Class I medical device that has been manufactured in accordance with the requirements of Directives 93/42/ EEC and 2007/47/EC. These instructions have been written in conformity with these directives, and their aim is to ensure the safe and proper use of this medical device.

TRADEMARKS OF MATERIALS

Velcro® is a registered trademark of Velcro Industries B.V.

PRECAUTIONS

It is advised that the pressure exerted by the product does not act on parts of the body suffering from wounds, swelling or lumps. The product should not be tightened too firmly so as to avoid the development of pressure points or compression of the nerves or blood vessels underneath. It is advisable to wear the product over a garment, avoiding direct contact with the skin. If in doubt about how to apply, contact an orthopaedic technician. Please read the composition of the product on the inside label carefully. The device should not be worn in the vicinity of open flames or strong electromagnetic fields.

WARNINGS

The product has been designed for the specific indications stated below and must be prescribed by a doctor and fitted by an orthopaedic technician. The orthopaedic technician is responsible for fitting the product and providing information about how to use it safely in accordance with individual requirements. The product must be fitted with the utmost care to ensure its efficacy, tolerability and correct function. Any changes to the structure or adjustment of the device must be prescribed by a physician and performed by an orthopaedic technician. Under no circumstances should the adjustment made by the physician/orthopaedic technician be altered. The manufacturer is no longer liable if the product is used or adapted inappropriately. The orthosis is intended to be used only by a single patient; if it is used for more than one patient, the manufacturer accepts no responsibility or liability in accordance with medical device legislation. Direct contact with the skin may cause reddening or irritation in hypersensitive individuals. In the case of onset of pain, swelling, lumps or any unusual reaction, immediately contact your physician. The orthopaedic efficacy of the product may only be guaranteed when all of its components are used.

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INDICATIONS

- Lombalgia
- Sciatic lumbalgia and lumbar cruralgia
- Spondylarthrosis and lumbar discopathies
- Slight traumas to the lumbosacral column and paravertebral lumbar contractures
- Post-operation period
- Results of vertebral lumbar relaxations
- Treatment of fracture of transverse processes
- Only REF. 1322 - 1322SH - 1354 - 1622 - 1654 - 1321 - 1353:**
- Degenerative lumbar scoliosis
- Spondylosis with or without vertebral listhesis

CONTRAINDICATIONS

Currently not known

CHARACTERISTICS AND MATERIALS

REF.1322-1322SH-1323-1354-1355-1622-1623-1654-1655

- Traction system with elastic straps and loops, for targeted pressure and perfect stabilisation
- 4 posterior guides for insertion of rigid reinforced stays (optional) to replace the ones provided
- Shoulder pads with axillary padding
- Velcro® front fastening
- REF. 1622 - 1623 - 1654 - 1655:**
- Posterior insert of reticular fabric for maximum permeability
- REF. 1622 - 1623 - 1654 - 1655: Main fabric made of cotton with a small amount of elastane
- REF. 1322 - 1322SH - 1323 - 1354 - 1355: Fabric made of cotton with a small amount of elastane
- Individually adaptable to the patient's pathology

REF. 1321 - 1353

- Lateral Velcro® closure with accompanying pocket
- Special shaping to favour comfort and adherence
- Ladder-proof hems
- 4 posterior guides for insertion of rigid reinforced stays (optional) to replace the ones provided
- Traction system of wide elastic stretch straps for targeted pressure and perfect stabilisation
- Traction straps with loops, for support and easy wearing
- Velcro® front fastening
- Fabric made of cotton with a small amount of elastane
- Individually adaptable to the patient's pathology

PUTTING ON THE APPLIANCE