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Index

1. GENERAL INFORMATION.....	4
1.1 MANUFACTURER	4
2. PRODUCT DESCRIPTION AND INFORMATION	4
2.1 DESCRIPTION	4
2.2 INTENDED USE.....	4
2.3 USERS.....	4
2.4 GROUPS OF PATIENTS AND CLINICAL CONDITIONS.....	4
2.5 USAGE ENVIRONMENT	4
2.6 STORAGE AND TRANSPORT	4
2.7 MANUFACTURER'S STATEMENTS.....	4
3. GENERAL WARNINGS.....	5
3.1 SERIOUS ACCIDENTS.....	5
3.2 SYMBOLS ON THE LABELS AND ON THE DEVICE.....	5
3.3 SYMBOLS IN THE MANUAL.....	6
3.4 SYMBOLS PRESENT ON THE PACKAGING.....	6
3.5 SPECIAL WARNINGS	6
3.6 CONTRAINDICATIONS AND SIDE EFFECTS	6
3.7 RESTRICTIONS ON USE	6
3.8 LIST OF MOST USED FUNCTIONS	7
4. AVAILABLE CONFIGURATIONS.....	7
5. PACKAGE CONTENT AND PRODUCT FEATURES	8
5.1 PACKAGE CONTENTS	8
5.2 PRODUCT.....	8
5.2.1 STRUZZO 500 AND 500 RC	8
5.2.2 EASY-UP 100 AND 200.....	8
5.2.3 STAND-UP 100-500.....	9
5.3 THE JOYSTICK	9
5.4 THE HANDSET	10
5.5 THE REMOTE CONTROL.....	10
6. ASSEMBLY AND FIRST USE.....	11
6.1 ASSEMBLY	11
6.2 FIRST USE.....	11
6.3 INSTRUCTION FOR CORRECT USE	11
6.4 NON-COMPLIANT USE	11
6.5 UNSAFE USE	11
6.6 DISCLAIMER.....	11
7. USE.....	11
7.1 WARNINGS BEFORE USE.....	11
7.2 INITIAL CONFIGURATIONS	12
7.2.1 POSITIONING THE JOYSTICK.....	12
7.3 USE OF THE DEVICE	12
7.3.1 SWITCHING ON AND OFF	13
7.3.2 HEIGHT ADJUSTMENT	14
7.3.3 PARKING BRAKE.....	15
7.3.4 TRANSFER TO THE DEVICE.....	15
7.3.5 LIFTING	16
7.3.6 MOVEMENT WITH LIFTED PATIENT.....	17
7.3.7 DESCENT AND CONCLUSION OF THE SESSION.....	18
7.3.8 BATTERY CHARGING	18
7.3.9 REMOTE CONTROLLED MOVEMENT	19
7.4 EMERGENCY DESCENT.....	20
7.5 ADJUSTMENTS.....	20

7.5.1 FOOTPAD ADJUSTMENT	20
7.5.2 ADJUSTMENT OF THE KNEE SUPPORTS	20
7.5.3 ADJUSTING THE PATIENT SUPPORT	21
7.5.4 ADJUSTMENT OF THE LATERAL CONTAINMENT	21
7.5.5 ADJUSTMENT OF THE SUPPORT PLANE.....	21
7.5.6 ADJUSTMENT OF THE ANTI-CRUSHING SYSTEM	22
8. ACCESSORIES.....	23
9. MAINTENANCE AND USEFUL LIFE	24
9.1 ROUTINE MAINTENANCE.....	24
9.2 NON-ROUTINE MAINTENANCE AND REPAIRS.....	25
9.3 BATTERY REPLACEMENT.....	25
9.4 SERVICE LIFE.....	25
10. CLEANING	25
10.1 CLEANING.....	25
10.2 DISINFECTION	25
11. DISPOSAL	26
12. TECHNICAL SPECIFICATIONS	26
12.1 LABELLING.....	27
13. ELECTROMAGNETIC COMPATIBILITY.....	27
13.1 ESSENTIAL PERFORMANCE.....	27
13.2 ELECTROMAGNETIC EMISSIONS.....	27
13.3 RF SIGNAL TRANSMISSION AND RECEPTION	27
13.4 ELECTROMAGNETIC IMMUNITY.....	28
13.5 RECOMMENDED SEPARATION DISTANCES	28
14. TROUBLESHOOTING.....	29
15. WARRANTY.....	29
15.1 GENERAL TERMS AND CONDITIONS	29
15.2 WARRANTY REPAIRS.....	30
15.3 OUT OF WARRANTY REPAIRS	30
15.4 NON-DEFECTIVE PRODUCTS	30
15.5 REPAIRS AT HOME	30
15.6 SPARE PARTS.....	30

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
1. GENERAL INFORMATION

This manual contains information for a practical, correct and safe use of the device and is intended for reading by specialized personnel and the user of the product. It is recommended that you carefully read the entire manual before using the product. If there are any doubts or clarifications, please contact the retailer who will be able to advise you correctly.

The importance of reading and understanding the user manual is highlighted on the product by the following symbol:

1.1 MANUFACTURER

CHINESPORT SPA - Via Croatia, 2 - 33100 Udine, Italy

 +39 0432621621

 <https://www.chinesport.com>

 chinesport@chinesport.it

The company produces in compliance with the UNI EN ISO 13485: 2016 Quality System

2. PRODUCT DESCRIPTION AND INFORMATION

2.1 DESCRIPTION

Standing frames are adjustable in height and have four wheels, at least two with brakes. The decision on which model to use is based on the patient's residual motor and coordination skills and can only be made by physicians with adequate specialization.

Depending on the model, the user can get up manually, clinging to the large side handles and the edges of the tray, or by means of an electric motor with rechargeable battery, operated by hand control, or servo-assisted by a gas spring. It features oversized handles to provide standing support, as well as a large tray to hold items.

The aid is easy to access and can be adapted to the user's morphological characteristics through a wide range of adjustments. In particular, the main body can be adjusted according to the user's height by means of a servo-assisted system by the gas spring or by an electric motor.

2.2 INTENDED USE

Aid intended for people with reduced mobility, muscle strength, movements and stability to assume and maintain an upright posture. Some models also allow autonomous mobility within the home.

It is forbidden to use the device for purposes other than those defined in this manual. The manufacturer declines any responsibility for damage to persons or things deriving from improper use of the device or in any case different from that provided for in this manual. The manufacturer reserves the right to make changes to the product and to the manual without prior notice in order to improve its characteristics and performance.

2.3 USERS

People with reduced mobility, the elderly and the disabled, who can hardly bear their own weight and who cannot independently assume an upright posture.

It is up to the specialist to judge the patient's physical fitness for using the product. It is always recommended to use it under the supervision of the operator.

2.4 GROUPS OF PATIENTS AND CLINICAL CONDITIONS

There are no particular categories of patients that can be excluded except patients who show evident conditions of impossibility of access to the aid.

2.5 USAGE ENVIRONMENT

The aid is suitable for use in a clinical / hospital environment and in the home. Use on flat, dry, stable surfaces with no obstacles. Not usable outside.

Respect the climatic conditions of use

Temperature	Moisture	Pressure
0°C - 40°C	0% - 80%	76kPa - 106kPa

2.6 STORAGE AND TRANSPORT

Keep dry and at normal room temperature. Avoid excessive pressure and contact with discoloured materials. Avoid excessive exposure to direct sunlight.

Never charge the batteries if they are frozen or if they have a temperature higher than 45 °C.

Temperature	Moisture	Pressure
-25°C - 70°C	0% - 85%	76kPa - 106kPa

2.7 MANUFACTURER'S STATEMENTS

The manufacturer declares that

- the device is not a measuring instrument
- The device is not intended for clinical investigations
- the device is not sterile and is not for single use
- for correct operation and for the safety of the user it is necessary to carry out routine maintenance operations as described in the relative paragraph

- the device cannot be used outside the intended use stated in this manual

3.1 SERIOUS ACCIDENTS

In the event that serious accidents involving the device occur, the user is required to promptly notify the manufacturer and the competent authority of the member state in which the device is installed.

3.2 SYMBOLS ON THE LABELS AND ON THE DEVICE

	Follow user instructions.
	CE Certificate
	Dispose of properly
	Manufacturer
	Production date
	Indoor use only.
	Class II appliance
	Type B applied parts
	Type BF, CF applied parts
	Serial number of the appliance
	Unique identifier of the device
	Medical Device
	Identification of the country of production in ISO 3166-1 alpha-3 code (Italy)
IP21	Degree of protection against the ingress of solids and liquids: Protected against dripping water Protected against the ingress of solids > 12mm

3. GENERAL WARNINGS


For correct use of the device, always refer to this manual.

- Store the device in an environment suitable to what is reported on the labels on the packaging and on the specifications given in this manual
- There **useful life of the product is 10 years** in compliance with the correctness of the execution of the ordinary maintenance operations provided for in this manual. Battery life is highly dependent on the number of charge and discharge cycles and on the temperature of use. More details can be found in the paragraph on batteries. Use of the device beyond the declared useful life is strictly prohibited. At the end of the useful life it is possible to proceed according to the procedures illustrated in the relevant paragraph
- The manufacturer is not responsible, within the maximum limits permitted by the law in force, for direct or indirect, particular, incidental or consequential damages caused by:
 - o Incorrect use of the device
 - o Improper use of the device and outside the intended use
 - o Use of the device beyond the useful life stated in this manual
 - o Use of the device in environments not provided for in this manual
 - o Use with unsuitable patients
 - o Distraction of operators or incorrect application of commands and adjustments
 - o Use without first checking the status of the device as reported in the relevant paragraph
 - o incorrect or lack of maintenance
 - o Use with parts or accessories that are incompatible or not approved by the manufacturer
 - o Incorrect disposal or disposal other than that provided for in this manual

The device can be equipped with labels to draw attention to particular dangers such as:

	Danger of hand crushing		Foot crushing hazard
---	-------------------------	---	----------------------

Therefore, pay particular attention when carrying out operations in the areas adjacent to these symbols.

	Lifting of the patient if lifting with an electric motor (handset) is possible
	Lowering of the patient if lifting with an electric motor (handset) is possible
	Increased height of the structure if electric height adjustment is possible (push button panel)
	Reduction of the height of the structure if it is possible to electrically adjust the height (push button panel)
	Forward movement via drive wheels (remote control)
	Rearward movement via drive wheels (remote control)
	Left rotation, counter-clockwise, via drive wheels (remote control)
	Rotation to the right, clockwise, via drive wheels (remote control)
IPx4	Degree of protection against the ingress of solids and liquids (Joystick)
	Charging connector to connect the charger (Joystick)
	Maximum speed control, direction of speed increase (Joystick)
	Horn (Joystick)
	Turning on the guidance system (Joystick)
	Release position for height adjustment by gas spring
	Release position for assisted lifting by gas spring

3.3 SYMBOLS IN THE MANUAL

	Warning! Draws attention to safety-related information
	Joystick manoeuvring directions. Forward, back, left and right.

3.4 SYMBOLS PRESENT ON THE PACKAGING

	Fragile
	Recyclable
	Keep dry
	Keep this side up
	Do not use sharp blades to open

3.5 SPECIAL WARNINGS

Do not use when the workload exceeds the declared capacity.

Do not use with the patient unattended.

Do not use when the patient and / or equipment are not stable enough.

Do not use if the wheels are not firmly secured to the frame.

Do not use if the original product has been modified or tampered with.

Do not use with components even partially defective or following checks with a negative result

Do not use if configured incorrectly or if believed to cause harm to the patient or user

3.6 CONTRAINDICATIONS AND SIDE EFFECTS

There are no known side effects or contraindications. Always read the warnings before use.

It can generate involuntary muscle contractions.

Contraindicated to users with trunk and upper limb control deficits

Contraindicated to users with joint limitation of the knees and / or coxo - femoral and / or ankle.

Contraindicated to users with structural deficiency in the lower limbs.

3.7 RESTRICTIONS ON USE

Always use in accordance with the intended use and with patients who comply with the weight requirements indicated on the product label.

Maximum patient weight: 140 kg.

The physical and mental fitness of the user to use the equipment must be approved by a specialist.

Do not pour liquids directly on the equipment.

Do not use the lifting system continuously, respect the duty-cycle indicated in the TECHNICAL SPECIFICATIONS

Always unplug the charger when using and / or moving the appliance.

Disconnect the batteries before any intervention on the equipment circuits.

Do not connect to an electrical system that does not comply with the Standards in force for the same.

Do not disconnect the battery charger, the handset and the motor by pulling the cable, but always grasping the body of the plug.

Turn off the equipment when not in use to conserve battery power

3.8 LIST OF MOST USED FUNCTIONS

The most frequently used functions of standing stabilizers are listed below. Refer to the respective paragraphs where the functions are explained in detail:

- Lifting (7.3.5)
- Height adjustment (7.3.2)
- Movement (7.3.6)
- Adjustment of the knee supports (7.5.2)
- Wheel brake activation / deactivation (7.3.3)
- Battery charging (7.3.8)

4. AVAILABLE CONFIGURATIONS

The product is available in some configurations linked to specific functions of support, containment, stabilization and movement of the person which differ from each other in:

1. Base.
2. Footboard.
3. Knee supports.
4. Height adjustment.
5. Lifting.
6. Rear containment.
7. Lateral containment.
8. Accompaniment.
9. Table top.

For each element there are different options, which within the limit of compatibility between them produce more than forty configurations with significant differences that meet the most complex needs. For some elements there are accessories that increase the comfort and specificity of the function.

It is possible to fill in the code below in order to have a quick reference regarding the configuration of the purchased product.

A	V	3	*	*	*	*	*	*	*	*	*
01	02	03	04	05	06	07	08	09	10	11	12

The code can consist of the following alternatives:

Location	Val.	Description
BASE Position. 04	1	Wide with wheels, for push movements, easy to access.
	2	Narrow with traction wheels, for autonomous movement, easy driving.
	4	Narrow with traction wheels, for autonomous movements, with advanced control unit
	5	Narrow with traction wheels, for autonomous movement, with advanced control unit and remote control
FOOTPLATES Stand 05	1	In wood adjustable in depth with tools
	2	Metal adjustable in depth without tools.
KNEE SUPPORT Stand 06	1	Adjustable in width and depth.
	2	Depth-adjustable.
HEIGHT ADJUSTMENT Stand 07	A	Height adjustment assisted by a gas spring.
	B	Height adjustment assisted by an electric actuator.
LIFTING Stand 08	0	In autonomy with the strength of the arms.
	2	With the help of a gas spring operated by a lever.
	3	Servo assisted by electric actuator controlled by handset.
REAR CONTAINMENT Position 09	1	Mechanical operated manually.
	3	Standard removable saddle
	4	Long removable saddle
	5	Harness
LATERAL CONTAINMENT Position 10	1	Fixed, not adjustable in width.
	4	Short adjustable in width by sliding
	5	Long adjustable in width by sliding
ACCOMPANIMENT Position 11	2	Dynamic Parallel
	3	Static Fixed
	4	Static Adjustable
TABLE TOP Posture 12.	B	In small polyethylene
	C	In small Plexiglass
	F	In large polyethylene
	G	In large Plexiglass

Not all the obtainable combinations exist, for more information refer to the manufacturer's website or contact the manufacturer himself.

There are standard starting configurations, with a pre-configured code on which it is possible to vary some options such as the patient support, the foot rest platforms or the support surface.

5. PACKAGE CONTENT AND PRODUCT FEATURES

The product is delivered in suitable cardboard packaging so that it can be received intact and in working order. To open the packaging and extract its contents, pay attention to the warnings and symbols shown on the packaging itself. Dispose of packaging and waste material in an appropriate manner and following the information on the packaging and in this manual. The version for independent mobility has the new Joystick holder that is fixed to the handles only for shipping. Before use, find the best position for the joystick, see the paragraph on positioning the joystick.

5.1 PACKAGE CONTENTS

The pack contains

- Instruction manual
- Product already assembled.
- Packaging material to be disposed of

5.2 PRODUCT

Depending on the configuration chosen, the product may have some different characteristics.

5.2.1 STRUZZO 500 AND 500 RC

The variants indicated with the name of Struzzo 500 are those that have a base with drive wheels for autonomous mobility. The variants with RC indicate the availability of the remote control or its predisposition. The letters in the highlighted positions are those that can be configured as reported in the paragraph on the available configurations (4)

Model	Code
STRUZZO 500	AV3221B3442C
STRUZZO 500 RC	AV3521B3442C
STRUZZO 500 RCR	AV3421B3442C

1	Lifting motor
2	Patient Saddle
3	Foot rests
4	Joystick Handling
5	Handset for lifting and height adjustment
6	Lifting motor and height adjustment
7	Drive wheels
8	Knee supports
9	Table
10	Lateral containment
11	Castor Wheels with Brake

5.2.2 EASY-UP 100 AND 200

The variants indicated with the name of Easy-Up 100 and 200 are those that allow assisted lifting with an electric motor and differ from each other for the height adjustment, which can be with an electric motor or with a gas spring.

Model	Code
EASY-UP 100	AV3121A3452G
EASY-UP 200	AV3121B3452G

1	Lifting motor
2	Patient Saddle
3	Foot rests
4	Handset for lifting and height adjustment
5	Height adjustment (Motor or Gas spring)

- 6 Knee supports
- 7 Table
- 8 Lateral containment
- 9 Castor Wheels with Brake


5.2.3 STAND-UP 100-500

The variants indicated with the name of Stand-Up 100, 200, 300, 400 and 500 are those that allow a wider choice of lifting and height adjustment and differ from the easy-up models for the possibility of adjusting independently the height of the table with respect to the rest of the structure and does not depend on the lifting mechanism.

The models are divided into two macro categories:

- 1) Models with saddle
- 2) Models with sling or self-lifting

Model	Code
STAND-UP	AV3121A3454G
STAND-UP	AV3121B3454G
STAND-UP	AV3121A2454G
STAND-UP	AV3121A3554G
STAND-UP	AV3121A0113G


- 1 Lifting motor or spring
- 2 Patient saddle or rear harness or restraint
- 3 Foot rests
- 4 Handset for lifting and height adjustment
- 5 Height adjustment (Motor or Gas spring)
- 6 Knee supports
- 7 Table with independent adjustment
- 8 Lateral containment
- 9 Castor Wheels with Brake

For details regarding all the available combinations, refer to the commercial material or to the manufacturer's website.

5.3 THE JOYSTICK

The models that provide autonomous mobility have a driving system consisting of a joystick for moving and the drive wheels.

There are two Joystick models, based on the purchased aid model. They differ in shape, size and type of user experience, but both have common characteristics.

The compact model comes with the power button and top speed dial set further back than the control handle. Both Joysticks have the charging connector into which the battery charger plug, supplied as standard, will be inserted, and the button for the horn. The control knob is a lever that can be tilted 360° in space and provides an analog steering system that provides an even speed distribution from idle to maximum based on how much tilt the knob is given. To summarise:

	Speed dial of the maximum speed in the standard model. Turning clockwise increases the maximum speed
	Speed dial of the maximum speed in the compact model. Turning counter-clockwise increases the maximum speed
	Audible warning (horn)
	Switch on the steering system
	Charging connector for the battery charger

Since the movement takes place from the front drive wheels, it may be advisable to familiarize yourself with the use of the joystick control lever to understand how movement occurs in space. When switched on, the joystick shows the status of the machine and the battery charge through luminous indications.


The amount of LEDs lit indicates the percentage of battery charge.

	100% charged batteries
	Batteries charged to 80%
	Batteries charged to 60%
	Batteries charged to 40%, recharging is recommended
	Batteries charged to 20%, use is not recommended, charge immediately

Generally it is good practice to always keep the batteries charged, so daily recharging or recharging with low-discharged batteries is recommended compared to less frequent recharging or with very discharged batteries. Short charging cycles extend the life of the batteries. In certain situations, the charge LEDs take on different meanings, which are described in detail in the paragraphs below.

5.4 THE HANDSET

Some models of the aid have motorized height adjustment or motorized lifting support or both. These functions are activated by using the handset connected to the device via a long spiral cable, as can be seen from the descriptive photos. The handset has 4 buttons of different colours. One pair is dedicated to height adjustment (1) and one pair is for lifting (2).

1		Increased height of the structure
		Decrease in the height of the structure
2		Lifting the patient
		Lowering of the patient

5.5 THE REMOTE CONTROL

Some models have the remote steering system, designed exclusively to approach or move away the aid equipped with drive wheels. This function allows the patient both to approach the aid when he needs to stand upright and move independently, and to move the aid away when he has finished using it, but there are no people nearby who can assist him in this regard. In both situations, the situation of departure or arrival is the patient positioned on a seat or bed and the aid a few meters away. The control was not designed for use over long distances, so it has a limited range.

The remote control is coded and has a specific variable code encryption algorithm linked to a unique identification for each remote control. This ensures that only the remote control associated with the aid can work with it and does not interfere with any other similar aid that may be present in the vicinity.

Each remote control has a unique identification stamped on the back and the aid to which it has been associated has a label with the same identification, to quickly identify which remote control it has been associated with.

To make the best use of the remote control and understand the commands more quickly, it is recommended to use it by turning it 45 ° clockwise. In this way, the correspondence with the joystick movement commands is noticeable more quickly.

In this position it is easier to understand in which direction the movement of the aid will be commanded.

	Moving forward
	Moving backwards
	Left rotation, counter-clockwise
	Rotation to the right, clockwise

The remote control only works when the aid is in remote driving mode. This mode can only be activated by entering a specific sequence, to avoid abuse by third parties. This sequence is indicated in the chapters dedicated to using the device (See 7.3.1.4).

6. ASSEMBLY AND FIRST USE

6.1 ASSEMBLY

The device is delivered already assembled and ready for use. The version for independent mobility has the Joystick holder which is fixed to the handles only for shipping. Before use, find the best position for the joystick, see the paragraph on positioning the joystick. If the model is equipped with batteries, it is advisable to carry out a complete recharge before using the device.

6.2 FIRST USE

Before proceeding with the first use it is necessary to read this instruction manual.

Make sure that the product corresponds to what is described in this manual and that the contents of the package coincide with what is described in the relative paragraph.

Before using the product for the first time:

- Check that there are no device components inside the packaging (screws, nuts, knobs ...)
- Check for obvious signs of damage or tampering
- Check that the knee supports show no signs of damage
- Verify that all knobs or pins for adjustments are present
- Check the rotation and the sliding of the wheels
- Check the correct activation and effectiveness of the parking brakes
- Check the operation of the height adjustment control (electric or gas spring)
- Using the height adjustment control, bring the device to the minimum height and then bring it back to the maximum height
- Check the operation of the lifting control (electric or with gas spring if present)
- In the versions with autonomous mobility, turn on the joystick and check that the drive wheels move in the given direction.

6.3 INSTRUCTION FOR CORRECT USE

For a correct use of the device, the user moves on it from a stable and safe wheelchair or from a simple chair or from a sofa or bed with sufficient stability on the floor. The user must assume an upright sitting position and hold firmly to the handles of the AV3 and the chair or wheelchair must not be able to tip over when the user moves into the AV3 seat. It is recommended to be assisted during this phase or in any case to carry out tests in the presence of help to gain confidence with the operations to be done.

AV3 must be handled on a flat and solid surface, INDOOR only.

6.4 NON-COMPLIANT USE

AV3 is not suitable for:

- Diagnoses
- Tracking

- Measure
- competitive training
- Outdoor use

The use of autonomous handling requires responsibility and must be avoided in the presence of uneven surfaces, obstacles on the path, including carpets and the presence of stairs or hatches that are not properly protected.

6.5 UNSAFE USE

Do not use AV3 (electrical versions)

- In combination with other devices that emit ionizing radiation (e.g. radiotherapy, nuclear medicine, etc.)
- In environments where:
 - There are explosive materials
 - There is enriched oxygen
- In combination with the
 - Flammable anaesthetics
 - volatile solvents
- In environments with uneven surfaces, obstacles or hazards that are not protected or signalled.
- Do not use the remote control when the patient is on the aid.

6.6 DISCLAIMER

The manufacturer and its distributors assume no liability as a result of:

- Misuse, incorrect or unintentional use
- Failure to comply with this instruction manual
- Intention or gross negligence
- Too intense training or use, e.g. competitive activity
- Outdoor use or movements
- Use in combination with inadequate unstable wheelchairs, chairs, tables or beds
- Use of the device, although prohibited by the doctor or therapist
- Installation of unauthorized accessories
- Performing repairs or other interventions on the AV3 by persons not authorized by the manufacturer

7. USE

7.1 WARNINGS BEFORE USE

Before each use make sure that:

The product does not show obvious signs of tampering or damage

The product has been sanitized in the parts in contact with the user

The physical and clinical conditions of the user have been evaluated and that they are considered suitable for the use of the device

The environment of use is in accordance with the provisions of this instruction booklet

There are no particular dangers in the areas adjacent to the device (steps, obstacles, flammable materials, etc. ...)

There are no oil stains on the floor caused by the gas spring
If there are batteries, make sure they are always charged

7.2 INITIAL CONFIGURATIONS

When using the device for the first time, it is advisable to adjust the settings in relation to the pathological and physical characteristics of the patient. The setting can be made using all the adjustments provided for each option. Please contact your therapist in order to provide all support for the safest and most comfortable adjustment of the device and supervise the first use of the device to check if the device is used safely. See the paragraphs on adjustments (7.5) for all the adjustments that best suit the model owned.

An incorrect configuration could limit or prevent the correct use of the device

Adjust the knee support for maximum comfort and minimum stress on the joints during the lifting process

Adjust the height of the structure for maximum comfort and the best position of the body and limbs

Seat depth adjustment must be done before starting the transfer from the wheelchair


Before use, check that all adjustment elements are properly tightened

The distributor and the manufacturer are available to provide any support for the initial configuration and to provide all information regarding the safe use of the device. If you have any doubts about the safety information, do not hesitate to contact the manufacturer.

The aids are suitable for customizations that can also be made by the distributor or by the trusted healthcare provider but no changes must be made to the integrity of the product.

7.2.1 POSITIONING THE JOYSTICK

The AV3 models for independent mobility are equipped with a Joystick to control the movements of the device inside homes. The Joystick holder consists mainly of two adjustments that allow two degrees of freedom for the positioning and orientation of the Joystick itself. The first adjustment is used to reach the required position on the table or on the handles. The second adjustment is used for the ergonomic positioning of the Joystick.


The adjustment marked with "A" allows you to open and close the vice to position the Joystick in another point of the table or tubular. Once the point has been chosen, it is necessary to tighten the vice by acting on the adjustment "A" until the vice itself is firmly fixed in the decided point.

The adjustment marked with "B", on the other hand, allows you to orient the Joystick in space by rotating around the two support spheres:

Once the position has been decided it is possible to orient the joystick along the horizontal plane and along the vertical plane.

Once the most suitable position has been found, adjustment B must be tightened as much as possible in order to maintain this position. If you want a certain adaptability, tighten less, while if you want absolute stiffness, tighten as much as possible. Here are some examples of Joystick positioning seen from above and in profile:

The joystick support is intended to support the weight of the joystick only, it is absolutely not intended to support the patient. Therefore it is advisable to always keep a part of the joystick body in contact with the table and never lean with your weight on the joystick or on parts of it or of the support.

7.3 USE OF THE DEVICE

Before each use, check that all adjustment elements are properly tightened and that the batteries are fully charged. If you want to use the remote control, make sure you have the remote control handy and that it is charged

Read the instructions carefully to find the correct use based on the configuration of the device actually owned.

7.3.1 SWITCHING ON AND OFF

Some configurations of the AV3 device contain electric motors for:

- Movement, for independent mobility
- Lifting
- Structure height adjustment


These devices have a main switch that turns on the entire system and allows the use of electric motors. The main switch also activates battery charging when the suitable charger is connected to the device.

7.3.1.1 Main power switch

All models that have batteries have a general switch that serves to electrically disconnect the batteries from the circuit and preserve their charge. In addition, the switch acts as a³⁾ thermal protection in the event of a fault.

To switch on, switch the main switch located under the service table on the right to "I". To switch off, switch to "O".

In the case of AV35x and AV34x models, to activate the lifting⁴⁾ and lifting functions, the guide system must also be switched on, as shown below.⁵⁾

7.3.1.2 Activation of the driving system

The models that have the drive wheels for autonomous mobility have a Joystick to guide the aid inside the house, to be able to go from one room to another. Regardless of the model (see 5.3), on the Joystick there is a button to activate and deactivate the guide marked with the symbol:

To activate, press the power button on the joystick, indicated by JA in the figure. This will light up green and show the battery charge status on LEDs, indicated by LD in the figure.

To deactivate the guide, press the button again and it will switch off, together with the battery status LEDs. The horn button (JC) is used to emit a sound from the Joystick.

7.3.1.3 Joystick safety lock and unlock

Sometimes it is useful to lock the Joystick functions for safety reasons (e.g. children in the house). This way only trained people can unlock the Joystick and start the movement. To lock the joystick follow the instructions below:


- 1) When the device is turned on, press the power button (JA) for more than 4 seconds
- 2) The LEDs (LD) in the odd positions (1, 3, 5) will start flashing for a short time

The Joystick is now locked. If the JA button is pressed again, the LEDs (LD) in the odd positions (1, 3, 5) will flash briefly but the use of the device will not be possible. After 10 seconds the joystick switches itself off.

To unlock the joystick again, follow this procedure:

- 1) Push the power button. The LEDs (LD) in the odd positions (1, 3, 5) flash briefly and the use of the device will still not be possible
- 2) Press the horn button (JC) twice within 10 seconds. The Joystick LEDs will now light up permanently and device use will be enabled again.

7.3.1.4 Activation of the remote guidance system

The AV35x models include a remote steering control system which serves only to approach and remove the aid when assistance personnel are temporarily unavailable. This system allows you to move the aid a few meters using a remote control only when there is no one on it so that you have the possibility to approach or move it away from the position in which you are in order to get on or off the aid.

Activate the remote guide only and exclusively to move the aid a few meters and only when no one is on it

The remote steering system can only be activated if the main switch is switched on and the aid is within range of the remote

control. This range of action is intentionally reduced to a few meters, about 5m, and may be affected by external interference, the presence of walls, the state of charge of the remote control battery. In all ways it is advisable to first test the remote control in the environment of real use in order to be confident in the real range of the remote control.

The remote guide is used exclusively to approach and remove the aid when you are not on it to free your view and is not intended as an emergency system. Any failure of the remote control must not constitute a dangerous situation

To activate the remote guidance system, a specific sequence must be entered on the remote control, so that it cannot be activated by children or by unintentional key presses. Enter the sequence with a constant cadence with a pause of about 1 second between one key and the next.

					
Back	Back	Next	Back	Back	Right

If the sequence is entered correctly, it is possible to see the battery status LEDs on the joystick light up and after a few seconds two beeps from the joystick horn. From now on, the remote guidance system is active, and will automatically deactivate based on one of the following conditions:

- 1) The power button on the joystick is pressed
- 2) 10 seconds pass from the last press of any button on the remote control


The first condition is to bypass the guidance system at any time if someone accidentally enters the activation sequence while the patient is still on the aid.

The second condition is to avoid that the system remains activated for too long and consequently reduces the life of the batteries.

7.3.2 HEIGHT ADJUSTMENT

Before moving to the device, independently or with the help of an assistant, make sure that the height of the structure is already in the correct position and that all adjustments have been made. In the case of assisted standing it is possible to change the height and position of the seat by means of a gas spring or electric motors.

For easy access, if a seat has been chosen as a lifting support (see "Lifting support" - referring to position 9 of the product code), the seat can be aligned with the user's seat surface thanks to the adjustment in height. Such seat alignment as a starting point can simplify transfer and increase safety. In this way the user can move from his wheelchair, bed or sofa with this type of accessibility.

If the model chosen is without a seat, the height can be left unchanged. The height can be adjusted as shown in the following images (please only follow the instructions related to the equipment of your model).


7.3.2.1 Height adjustment by gas spring

In this case the height adjustment is not done very often, so it can be done when the device is first adjusted and when small adjustments are needed.

Identify the release lever of the gas spring for height adjustment, located on the lateral containment handle near the respective symbol:

Pull the height adjustment lever to release the gas spring and change the height of the frame with little force up or down. Release the lever when you have reached the desired setting.

7.3.2.2 Height adjustment by electric motor


When lowering is in progress, the column comes very close to the base. Be careful not to put your feet in the trapping area to avoid the risk of crushing

In this case the device is more versatile and the height adjustment can be done more often without effort with the use of the hand control.

Push the button  on the handset for **to increase** the height of the structure or

Push the button  on the handset for **to decrease** the height of the structure.

When the height of the structure is completely decreased, in case the device has independent mobility motors, the movable wheels are raised and there are small wheels useful for moving the device around (without people on it).

7.3.3 PARKING BRAKE

The parking brake of the rear and front wheels can be used to prevent unwanted movement of the device when it is not in use or when the patient is on it. The devices with drive wheels for autonomous mobility have the parking brake only on the rear wheels as the drive wheels are self-braking. To activate the brake, use your foot to operate the lever on the wheel you want to lock downwards.

To deactivate the brake, use your foot to operate the lever on the wheel to be released upwards.

7.3.4 TRANSFER TO THE DEVICE

Before transferring the patient to the aid, the patient's starting seat must be aligned with the aid seat, or the proximity of the aid in models where a seat is not provided. The device can be approached directly with a wheelchair or can be moved by other people near the patient's seat. In models with remote guidance, the device can be approached using the remote control.

When the patient is in front of the device ready for transfer, follow the instructions below.

7.3.4.1 Preparation for the transfer

Before starting the transfer, the following conditions must be met:

- The rear and front wheels are locked (see parking brake) or if there are drive wheels, check that the latter are well braked
- The wheels of the wheelchair are locked (see the wheelchair user manual) or the seat is firmly on the floor.

7.3.4.2 Transfer in the absence of assisted lifting

In models where there is no assisted lifting (stand-up 500), the transfer must take place in a single solution and at the same time the posterior retention system must be closed.

During lifting, the only support for the patient is his own strength or the strength of the person assisting him. There is no back protection until the back support is closed and locked

Follow the steps below:

- Open the back support until it is fully open.
- Place your feet on the support platforms
- With your own strength or with the help of an assistant, hold onto the side handles and lift yourself up.
- close and lock the rear support in place.

Bad starting position

Good starting position

Once the aid seat is within reach, follow the steps below to transfer to the device:

- Place your feet on the support platforms
- With your own strength or with the help of an assistant, hold on to the side handles and transfer to the device

If you have a model with a sling as a patient support, instead of transferring to the device, you must hold yourself so that the person providing assistance can pass the sling under the patient and hook it to the relative supports.

Only when you are in the seat or harness can you move the wheelchair away.

7.3.5 LIFTING

Once the patient is on board of the aid, he can be placed in an upright posture using the lifting assistance present on the device. Only in the model with no lifting (stand-up 500) the patient must verticalize with his own strength or with the help of another person.

7.3.5.1 Lifting with partial assistance via gas spring

In this case there is a gas spring which provides partial lifting aid and reduces the patient's weight. With your own strength or with the help of an assistant, hold on to the side handles and prepare for movement, looking for the release lever of the lifting gas spring, located on the lateral containment handle near the respective symbol:

7.3.4.3 Transfer with assisted lifting

In models where lifting is assisted with a motor or gas spring, the patient support, whether it be the seat or the harness, must be lowered to the minimum level in order to facilitate the transfer of the patient to the aid. In the case of lifting assisted by a gas spring, it is necessary to release the lifting spring and guide it with the hand towards the lowest position.


In the case of motor-assisted lifting, the handset must be used to lower the patient support to a minimum, using the lifting descent button.

Push the button  on the handset for **lower** patient support (saddle or harness)

Press the buttons  or  on the handset for **decrease or increase** the height of the structure if necessary

The lowering step of the structure should be performed only if necessary, i.e. when the seat level is much higher than that of the chair where the patient is sitting:

The gas spring is not powerful enough to support the entire weight of the patient, so once the weight is released, albeit reduced by the assistance of the spring, it rests on the patient's limbs

When ready:

Pull the lifting lever to release the gas spring and at the same time carry out the lifting movement of the body. The gas spring will help to partially lift the patient while reducing the body weight. The lifting phase can be performed in several steps since the lifting movement can be stopped by releasing the lifting lever. In this way the patient can exert the effort later and become less tired. The same rules apply to the descent.

7.3.5.2 Lifting completely assisted by electric motor

In this case there is an electric motor that completely assists the lifting of the patient.

Do not use the lifting motor for a long time. Observe the duty cycle indications on the label. Incorrect stage setting could generate compression on the user during lifting if the anti-crush attachment is not present. Please monitor the patient's condition and reactions while the lift is in progress. A wrong configuration of the knee supports can put stress on the joints during the standing process.

When Ready:

On the handset, locate the buttons for lifting:

Push the button  on the handset for **raise up** patient support (saddle or harness)

In case of need,

Push the button  on the handset for **lower** patient support (saddle or harness).

It is advisable to take breaks during lifting to check the patient's state of health. Moving from a sitting to a vertical position often brings about some sudden changes in the patient's vascular system and also has repercussions on the urinary and digestive systems.

If the anti-crushing accessory is present, a possible collision of the patient's body immediately stops the ascent function. The function will be inhibited as long as the button is pressed. The same recommendations also apply if the patient support consists of a sling. What differs is the final position which is less upright:

7.3.6 MOVEMENT WITH LIFTED PATIENT

Once the vertical position has been reached, whatever the model owned, it is possible to move the aid. If the model is with drive wheels for autonomous mobility, the patient himself can manoeuvre the device so that he can reach the desired place in the house. In the case of all the other models,

the movement is carried out manually by the person assisting the patient.

7.3.6.1 Assisted mobility

In the case of manual movement, the person assisting the patient takes care of moving the structure with the patient on board. The wheel brakes must be deactivated before carrying out the manoeuvre. Deactivate the parking brakes as indicated in the relevant paragraph (7.3.3) and carefully move the aid using the swivel castors.

Move with caution and without significant acceleration.
Only move the aid indoors, on flat and obstacle-free surfaces.
The patient should use the lateral restraint structures to hold on while moving.

Once the desired position has been reached, it is necessary to reactivate the parking brakes, reversing the previous operation.


7.3.6.2 Mobility in autonomy

In the case of motorized movement, the patient himself who is on board the aid can manoeuvre the movement using the Joystick.

Be extremely careful when using the motorized movement
Use independent mobility only indoors and only on flat, unobstructed surfaces
Please always set the minimum required speed
Do not use height adjustment or lifting during movement

Regardless of the joystick model (see details in paragraph 5.3) once the guidance system has been activated (see paragraph 7.3.1.2) to start the movement simply move the control knob (JM) in the desired direction.

In detail, the steps are as follows:

- If this has not been done before, turn on the joystick by pressing the power button (JA) to activate the guidance system

- The Joystick turns on and the status LEDs (JE) indicate the battery charge level. When all the LEDs are lit, the battery is fully charged. If instead of lighting only the odd LEDs flash, the Joystick is locked, it must be unlocked (see 7.3.1.3)
- To move, use the control knob (JM) by tilting it in the desired direction. The speed and acceleration are adjusted based on the inclination. The maximum speed is obtained when the joystick control lever is tilted to its maximum extent
- To change the maximum reachable speed, and therefore the adjustment range of the knob, use the speed dial (see paragraph 5.3)

7.3.7 DESCENT AND CONCLUSION OF THE SESSION

When the patient wants to end the lifting or autonomous movement session or no longer needs to use the aid, to get off the aid he can carry out the operations he has carried out in reverse. All the operations that require the assistance of third parties are the same and it is necessary to make sure that the destination seat is stable and does not move with the weight of the patient. To transfer from the device to the seat, once the patient support has been lowered, the person can transfer independently or with the help of third parties. All the recommendations set out in the respective paragraphs are valid. At this point, if you do not plan to use the aid further, it is advisable to charge it or turn off the main switch if you do not need to use additional functions, such as moving by remote control.

7.3.8 BATTERY CHARGING

When the device is not in use, it is recommended that you recharge the battery for future use. Daily recharging of the batteries improves their expected life and prevents the device from discharging during use.


Never leave the device with dead batteries or never use the device when the batteries are not fully charged.
Use only the charger supplied by the manufacturer


To recharge the batteries, the main switch must remain on. The charging method depends on the model.

7.3.8.1 Charging models with drive wheels

Models equipped with drive wheels for autonomous mobility have the charging connector located on the Joystick body (see also 5.3).

Standard Model

Compact Model

Once you have located the XLR connector, insert the charging plug and connect the charger to the power socket to start charging.

The charging process takes about 6 hours but sometimes a complete cycle of 10 hours is required. The charger has an orange light on it, which lights up when the charging process is in progress. When the battery is fully charged, a green light replaces the orange light.

If the orange light does not come on, check that the main switch is in the "I" position.

7.3.8.2 Charging models without drive wheels

All models that have at least one electric motor, be it lifting / lifting or height adjustment, must be recharged using the charging connector located on the front of the device.

Once the jack is inserted into the device socket, connect the charger to the power outlet to start charging.

The charging process takes about 10 hours. The charger has a yellow light that lights up when the charging process is in progress. When the battery is fully charged, a green light replaces the yellow light.

7.3.9 REMOTE CONTROLLED MOVEMENT

Some models have the possibility of being moved by remote control only to store the aid when it is no longer needed or to resume use of the aid when needed. An example could be the following:

The patient wants to sit in the chair or chair after using the aid and watch TV. After getting off the aid on an appropriate seat, he needs to move the aid a few meters in order to free the view from the TV.

Similarly, the patient finishes watching TV and wants to approach the aid in order to be able to stand upright and move independently.

Activate the remote guide only and exclusively to move the aid a few meters and only when no one is on it

Although the speed of movement via the remote control is limited, attention must always be paid to every slightest movement because any collision with the aid can cause injury or trauma

To be able to move the device using the remote control, it is necessary to activate the remote guide using the activation sequence (see paragraph 7.3.1.4). Once the sequence has been entered, the device can be moved using the remote control. No more than 10 seconds should pass between one command and the next, otherwise the device will automatically deactivate remote guidance and remain in stand-by.

Backward	Backward	Forward	Backward	Backward	Right

	Moving forward
	Moving backwards
	Left rotation, counter-clockwise
	Rotation to the right, clockwise

It is advisable to proceed with extreme caution when handling the remote control and to familiarize yourself with the controls by practicing in a large room without obstacles.

The remote control is not intended as a safety function and any malfunction of the remote control or the discharge of the battery of the remote control or the aid must not constitute a dangerous situation.

7.4 EMERGENCY DESCENT

These indications must necessarily be carried out by the person assisting the patient on board the aid

If for some reason the device stops working with the patient still in an upright position (models with electric motor assisted lifting only), it is necessary to establish if the problem depends on whether the battery is completely discharged or not. If the batteries are completely discharged, you can put the device in charge using the appropriate charger. After a few minutes the charge will be sufficient to lower the patient support to allow it to descend from the device.

If it is not dependent on the batteries, follow this procedure to free the patient from the device.

The following emergency manoeuvre must be performed by the person caring for the patient. If the patient is not able to support himself autonomously on the lateral containment structures, it is advisable to have an additional person present to support the patient. From the moment the support pin is removed the patient support is free to fall, be careful to stop the fall

- Approach a suitable seat on which the patient will be placed after the descent and position it so as to facilitate the descent of the patient support
- Have the patient hold on the lateral containment structures or on the table or if this is not possible, have it supported by a person
- While holding the patient support (if present), open the spring and remove the top safety pin from the actuator (be careful, the patient support is now free to fall down). If present, remove the patient sling.

- Fully lower the patient support to the floor

- Help the patient to reach the support seat.

- Contact the technical service to identify the cause of the malfunction

7.5 ADJUSTMENTS


The adjustments below do not apply to all models, please take those that apply to your device as valid.

7.5.1 FOOTPAD ADJUSTMENT

In most of the models there are platforms for foot support, in metal or wood. Both are adjustable in depth using a tool or simply by unscrewing a knob in the case of the metal footrests. The latter also, unlike the wooden platform, are individually adjustable,

An incorrect setting of the footplates could generate compression on the user's lower limbs during lifting

- Loosen the 4 screws (S and S') with the wrench supplied for the wooden platform, or
- Loosen the 2 knobs (K or K') of the metal footplate you want to adjust. Do not unscrew completely, it is enough to loosen the tightening.
- Slide the footplate in the desired direction in or out of the starting position
- Once the position has been reached, re-tighten the screws (S and S') or the knobs (K or K') to lock the footboard again.


7.5.2 ADJUSTMENT OF THE KNEE SUPPORTS

The knee supports can be adjusted in height, width, depth and rotation.


An incorrect setting of the knee supports can interfere with the kneecap (too high) or can cause the feet to slide off the footplates (too low).

Furthermore, when lifting, the joints can be overstressed and comfort is compromised

The height adjustment (green arrow) takes place as a whole and affects both knee supports. The other adjustments (red arrows) instead take place individually on each knee support. To adjust the height, loosen the T-shaped knob in the center of the bar that holds the knee supports (T), while to adjust the depth, width and rotation of the single knee support you need to loosen the quick release lever (L).

Once the position that guarantees maximum patient comfort has been found, to be checked especially during the lifting phase, the knob or the quick release levers must be tightened again.

7.5.3 ADJUSTING THE PATIENT SUPPORT

Some models have a short seat as patient support, other models have a short seat and others have a harness. Both seats can be adjusted in depth.

Do not pull the seat beyond the red limit indicated on the seat stem
Before moving to the device, check the correct tightening of the fixing handle

To adjust the seat carry out the following steps:

- Loosen the quick release lever

- Adjust the desired depth
- Tighten the quick release lever

7.5.4 ADJUSTMENT OF THE LATERAL CONTAINMENT

All the models are equipped with lateral containment. Regardless of the shape, you can adjust the width of the side containment handles.

Before moving to the device, check the correct tightening of the adjustment knobs of the lateral containment handles

To adjust the width of the side containment handles carry out the following steps:

- Loosen the respective T-knob under the table (T or T')
- Achieve the desired handle width by sliding the handles
- Tighten the T-knob

7.5.5 ADJUSTMENT OF THE SUPPORT PLANE

The support surface in most models can be adjusted in inclination. In some models it can also be adjusted in height and depth.


Before use, check that the pins are properly locked
When adjusting the depth of the stage, check that it could not cause compression in the raised position

To change the inclination of the table follow the steps below:

- Press the blue button and completely remove the locking pin (P)

- Tilt and hold the stage in position until the desired inclination is reached within the four available positions (visible from the fixing holes)
- Reinsert the pin and release the blue button



In the models where the configuration provides for fixed or adjustable accompaniment, the height of the support surface does not depend on either the height adjustment or the lifting, and can be adjusted independently. In these models the support surface can be adjusted, in addition to the inclination, also in height and depth.

To carry out these adjustments there are two clamping elements, that is a T-knob to tighten the movement (T) and a positioning pin which serves to give preset positions to the adjustments (P).


To adjust the height of the support surface, follow the steps below:

- Loosen the T-knob near the right side handle (T)
- Firmly support the support surface and at the same time unlock the positioning pin near the left side handle (P). Be careful because all the weight of the top will be free to slide
- Once the desired height has been reached, insert and lock the pin by snapping it into one of the holes available for adjustment
- Tighten the T-knob

Similarly, to adjust the depth of the table, follow the steps below:

- Loosen the T-knob on the right side (T)
- Unlock the positioning pin on the left side (P)
- Adjust the depth of the table
- Once the desired depth is reached, insert and lock the pin by snapping it into one of the holes available for adjustment
- Tighten the T-knob

7.5.6 ADJUSTMENT OF THE ANTI-CRUSHING SYSTEM

The aid models that have the motorized lifting support, have an anti-crushing support as a basic accessory that prevents the patient's abdomen from being crushed in the event of a collision during the motorized lifting. The button immediately stops the lifting function and only allows descent as long as the button is pressed.

The button can be adjusted in depth to adapt to the patient's body and better calibrate the collision point.

After the adjustment, check that the button works, starting the lifting and pressing the button to stop it.

Find the position so that any collision with the abdomen does not cause any crushing

To adjust the anti-squeeze button in depth follow the steps below:

- Press the blue button and completely remove the locking pin (P)
- Slide the support forward or backward to find the correct position of the anti-crush system

- Insert the locking pin and look for the exact hole so that the pin is fully inserted again and release the blue button

AC0691 - SHORT SEAT - ADJUSTABLE

Seat length 24 cm, adjustable in depth. Weight 4.9kg
Dimensions: 52.4 x 17.4 x 14.2 h cm
Pad sold separately (AC0045)

AC0045 - SHORT SEAT BASE PADDING

Basic padding for the AC0691 adjustable short seat for added comfort

AC0543 - SHORT SEAT PADDING

Breathable padding, easy to wash / sterilize, guarantees flexibility and comfort while maintaining its shape, for the short seat AC0691.

AC0871 - HYGIENE SEAT

Hygienic seat in fiberglass, easily sanitized. Depth-adjustable. Weight 2.8 kg.
Dimensions: 61.7 x 32.2 x 18 h cm

AC0690 - LONG SEAT - ADJUSTABLE

Long saddle. Seat length 36 cm, adjustable in depth.
Weight: 5.5 KG
Dimensions: 58 x 20 x 14.2 h cm. Padding included (AC0046)

AC0046 - LONG SEAT BASE PADDING

Basic padding for the AC0690 adjustable long seat for added comfort

8. ACCESSORIES

The aids have a wide range of accessories available, many of which are compatible with all models.

For the complete list of accessories and compatibility with your device model, contact your distributor or refer to the website www.chinesport.it

Use only original Chinesport accessories

AC0686 - HEEL REST

Individually adjustable in width for metal footrests. The height of the heel stops is 2.5 cm from the footrests. Stabilization strap included

AC1300 - ADJUSTABLE HEEL RESTS

Individually adjustable in width and rotation for wooden platforms. The height is 4 cm. Applicable in two positions, one central and one more external

AC1271 - PADDED KNEE SUPPORTS KIT

In foam with leatherette cover. Independently adjustable in width and depth, as well as simultaneous height adjustment


AC0544 - LONG SEAT PADDING

Breathable padding, easy to wash / sterilize, guarantees flexibility and comfort while maintaining its shape, for the long seat AC0690.

AC0048 - POCKET FOR LONG SEAT

Useful for the user with difficulty in transferring from the wheelchair to the seat. It is recommended to use with the long adjustable seat (AC0690)

AC0049 - LATERAL SUPPORTS

In self-skinning polyurethane foam. Resistant, shockproof and easy to sanitize. Applicable only to short type handles, in fixed position

AC0702 - ADJUSTABLE LATERAL SUPPORTS

In self-skinning polyurethane foam. Resistant, shockproof and easy to sanitize. Applicable only to long type handles, adjustable along the handles

AC0014 - LEG STRAP

For greater stabilization of the lower limb. It is applied at the level of the knee support

AC0783 - BACK STRAP

For containment during lifting. Applied by the assistant to ensure the user in cases of difficulty in controlling the trunk. One size adjustable in length

AC0693 - FRONT SUPPORT

Front support to be used in the absence of adequate control of the trunk by the user. Dimensions: 30 x 38.5 x 26 h cm

AC1216 - TRUNK SUPPORT

Wide support to be used in the absence of adequate control of the trunk by the user. Only Plexiglass, to be applied to the existing safety anti-crushing system. Dimensions: 30 x 38.5 x 26 h cm

AC1124 - STORAGE BAG

Bag with velcro straps for easy installation on both sides. Useful for carrying magazines, newspapers, tablets, phones, remote controls

0168 - HAND ANCHOR

Ergonomic grip, can be applied to obtain greater stability

9. MAINTENANCE AND USEFUL LIFE

9.1 ROUTINE MAINTENANCE

For the correct functioning of the device, for the correctness of use, for the safety of the patient and the operator and to guarantee the performance of the device over time, it is necessary to carry out scheduled maintenance interventions. Ordinary maintenance interventions can be carried out by specialized operators or by authorized technical personnel.

FREQUENCY	CATEGORY	CHECK
Every 6 months	Visual inspection of the device	Integrity of the knee supports
		Check the condition of the driving wheels
		Integrity of the frame and absence of rust or unsoldered elements

		Check the braking wheels and brake conditions
		Check the tightening of screws and adjustment elements
		Check the integrity of the plastic covers
		Check the condition of connectors and plugs
	Functional test	Check that the height adjustment is working correctly
		If present, check the operation of the gas springs, the thrust, the sliding and the seal
		Correct rotation of the braking wheels
		Correct rolling of the wheels
		Effectiveness of the parking brakes
		Experience the comfort of the padding of the knee supports
2 years	Instrumental and functional checks	Electrical controls test
		Check the operation of the anti-crush button
		Checking the general operation of the device
	ELECTRICAL CHECKS	Control of electronic devices and wiring
Check for wear and, if necessary, replace: motors, wheels, electronic components		
		Electrical safety tests

It is recommended to carry out the electrical safety tests (CEI EN 62353 / EN 60601-1) of the device at least once every two years.

9.2 NON-ROUTINE MAINTENANCE AND REPAIRS

Extraordinary maintenance can only be performed by personnel authorized by the manufacturer, under penalty of immediate termination of any warranty conditions. The manufacturer declines all responsibility if tampering with the original product is ascertained.

Any malfunctions encountered by the user must be promptly reported to the distributor or directly to the manufacturer and inhibit the use of the device.

Repairs can only be carried out by technical personnel authorized by the manufacturer and may include the withdrawal of the device to carry out the necessary repairs.

The list of spare parts is provided directly by Chinesport SpA on request and includes ordering codes and price lists.

9.3 BATTERY REPLACEMENT

For battery replacement, contact your local distributor or Chinesport customer service. If the instructions on the care of the device and the management of the recharge are strictly followed, the life of the batteries can even reach 4-5 years.

9.4 SERVICE LIFE

The useful life of the device is also defined at the beginning of this manual, and is approximately 10 years, provided that all the indications in this manual are respected.

At the end of the useful life it is possible to proceed in the following ways:

- 1) Dispose of the device as described in the "disposal" paragraph
- 2) Request the manufacturer to recondition and re-certify the device so that it can continue to be used

As reported in the "General warnings" paragraph, the manufacturer declines all responsibility for the use of the device beyond the useful life established in this manual.

10. CLEANING

The device must be cleaned and sanitized at the end of each use if the device is intended for different users.

Cleaning must be carried out with the device disconnected from the mains and with the main switch set to 0

10.1 CLEANING

Remove dust from metal parts using non-aggressive detergents and drying immediately after washing.

Clean the Plexiglass and plastic parts with a damp cloth with a neutral product, immediately drying the surfaces.

10.2 DISINFECTION

To disinfect the product, use disinfectants with a low chlorine content, such as AMUCHINA® 10% or equivalent solutions with 0.1% sodium hypochlorite concentration and with the following spectrum of action:

- bactericidal in the presence of interfering substances according to EN 1276: 1997;
- fungicide in the presence of interfering substances according to EN 1650: 1997;
- bactericidal surface test according to EN 13697: 2001;
- surface test fungicide according to EN 13697: 2001;
- active on Salmonella typhimurium according to EN 13697: 2001;
- active on: HIV, HAV - HBV, HCV;
- virucidal according to EN14476: 2005;
- active on H1N1 influenza virus according to EN14476: 2005.

For safe use, refer to the instruction knob included in the product.

If you have the Plexiglass top DO NOT use aggressive products such as alcohol or bleach. They could irreparably damage the plan

To disinfect the Plexiglass top, use the specific product for Plexiglass or methacrylate, such as the SANI-PLEX PRO product. Contact the manufacturer for more information.

11. DISPOSAL

The symbol shown on the equipment label indicates that the waste must be disposed of in a "separate collection".

Therefore, the user must give the waste to the differentiated collection centers set up by the local administrations, or hand it over to the retailer against the purchase of new equipment of an equivalent type. The differentiated collection of waste and the subsequent treatment, recovery and disposal operations favor the production of equipment with recycled materials and limit the negative effects on the environment and health that may be caused by improper waste management.

Illegal disposal of the product by the user involves the application of the administrative sanctions referred to in article 50 and following of Legislative Decree no. 22/1997.

12. TECHNICAL SPECIFICATIONS

Below are the technical specifications of the standard pre-configured models.

TECHNICAL SPECIFICATIONS	STRUZZO	EASY-UP	EASY-UP	STAND-UP
	500 - 500 RC	100	200	100
Base - accessibility	466 mm	630 mm	630 mm	630 mm
Overall dimensions (length x width x height)	90 x 60 x 85-127 cm	103.5 x 75 x 85-127 cm	103.5 x 75 x 85-127 cm	103.5 x 75 x 98-146 cm
Wheels	Ø 150 mm + Ø 75 mm	Ø 75 mm	Ø 75 mm	Ø 75 mm
Height of the frame from the ground	31 mm	31 mm	31 mm	31 mm
Parking brake	Motor brake	individual	individual	individual
Safe workload	140 kg	140 kg	140 kg	140 kg
Weight	95 kg	81 kg	79 kg	95 kg
Batteries	2x 12V 12Ah	2x12V 2Ah	2x12V 2Ah	2x12V 2Ah
Insulation class	II	II	II	II
Applied Parts	BF	BF	BF	BF
Power supply voltage	100-240V~	100-240V~	100-240V~	100-240V~
Battery charger	50/60Hz	50/60Hz	50/60Hz	50/60Hz
Degree of protection	IP21	IP21	IP21	IP21
Device class 93/42	I	I	I	I

AID MOVEMENT

Handling	independent	assisted	assisted	assisted
Movements type	power-driven.	manual	manual	manual
Maximum step surmountable	20 mm	-	-	-
Maximum practicable slope	2°	-	-	-
Command Type	Joystick - remote control	-	-	-

Max. speed	2.5 km/h	-	-	-
Footpad type	Metal	Metal	Metal	Metal
Knee supports	anatomical	anatomical	anatomical	anatomical
Structure height adjustment	Electric	Gas spring	Electric	Gas spring
Lifting	Electric	Electric	Electric	Electric
Patient ID support	long seat	long seat	long seat	long seat
Lateral containment	Short handles	Long handles	Long handles	Long handles
Superior structure	dynamic	dynamic	dynamic	fixed
Support table	small	large	large	large

ELECTRICAL PART

Autonomy of continuous use	90 min	90 min	90 min	90 min
Lifting cycles (maximum load)	20	20	20	20
Average battery recharge time	8-10 hours	6-8 hours	6-8 hours	6-8 hours
Battery charge level	5 LEDs	3 LEDs	3 LEDs	3 LEDs
Auto-off	included	included	included	included
Fault diagnosis and maintenance	On Joystick	On LED	On LED	On LED
Lifting stop button	included	included	included	included

TECHNICAL SPECIFICATIONS	STAND-UP	STAND-UP	STAND-UP	STAND-UP
	200	300	400	500
Base - accessibility	630 mm	630 mm	630 mm	630 mm
Overall dimensions (length x width x height)	103.5 x 75 x 98-146 cm	103.5 x 75 x 98-146 cm	103.5 x 75 x 98-146 cm	103.5 x 75 x 98-146 cm
Wheels	Ø 75 mm	Ø 75 mm	Ø 75 mm	Ø 75 mm
Height of the frame from the ground	31 mm	31 mm	31 mm	31 mm
Parking brake	individual	individual	individual	individual
Safe workload	140 kg	140 kg	140 kg	140 kg
Weight	87 kg	73 kg	83 kg	72 kg
Batteries	2x 12V 2Ah	-	2x12V 2Ah	-
Insulation class	II	-	II	-
Applied Parts	BF	-	BF	-
Power supply voltage	100-240V~	-	100-240V~	-
Battery charger	50/60Hz	-	50/60Hz	-
Degree of protection	IP21	IP21	IP21	IP21
Device class 93/42	I	I	I	I

AID MOVEMENT

Handling	assisted	assisted	assisted	assisted
Movements type	manual	manual	manual	manual
Maximum step surmountable	-	-	-	-
Maximum practicable slope	-	-	-	-
Command Type	-	-	-	-
Max. speed	-	-	-	-
Footpads type	Metal	Metal	Metal	Metal
Knee supports	anatomical	anatomical	anatomical	anatomical
Structure height adjustment	Electric	Gas spring	Gas spring	Gas spring
Lifting	Electric	Gas spring	Electric	Manual
Patient ID support	long seat	long seat	Harness	-
Lateral containment	Long handles	Long handles	Long handles	Fixed handles
Superior structure	fixed	fixed	fixed	fixed
Support table	large	large	large	large

ELECTRICAL PART

Autonomy of continuous use	90 min	-	90 min	-
Lifting cycles (maximum load)	20	-	20	-

Average battery recharge time	6-8 hours	-	6-8 hours	-
Battery charge level	3 LEDs	-	3 LEDs	-
Auto-off	included	-	included	-
Fault diagnosis and maintenance	On LED	-	On LED	-
Lifting stop button	included	-	included	-

12.1 LABELLING

The following labels are present on the device and on the packaging. All symbols are explained in the paragraph 3.2 . The serial number can be identified on the label should it be necessary to speak to the manufacturer's service center.

13. ELECTROMAGNETIC COMPATIBILITY

The manufacturer declares that the AV3 device including all original accessories and the power cable supplied with the device complies with the requirements of the EN 60601-1-2 standard for electromagnetic compatibility. As this is medical electrical equipment, it requires special precautions regarding electromagnetic emissions (EMC) and must be installed and put into service according to the EMC information provided in this manual.

The use of accessories, transducers or cables other than those specified or supplied by the device manufacturer could increase the electromagnetic emissions and decrease the electromagnetic immunity of the device. Use of this device near or stacked (on top or under) other devices should be prohibited as it may cause improper operation. If this type of

use is necessary, the device in question and the others involved must be kept under control to verify its correct functioning.

The device is not compatible with high frequency surgical equipment.

13.1 ESSENTIAL PERFORMANCE

When used in all intended environments and in accordance with this manual, the device ensures basic safety and essential performance. In the presence of electromagnetic interference, the device may undergo operational degradation but basic safety remains guaranteed and the following performances remain guaranteed:

- The device remains safe for the patient and the operator
- The device does not make unexpected movements
- The handling parameters of the system do not change
- The charger may shut down in the event of a voltage dip, but it will become operational again when the voltage is restored (LED on)
- During handling with the remote control, movement interruptions may occur but unwanted movements (departures or deviations) cannot occur

13.2 ELECTROMAGNETIC EMISSIONS

AV3 is intended to be used both in the home and in the clinical setting, with the exception of places near active high frequency (HF) surgical devices or in Radio Frequency (RF) shielded environments, in an electromagnetic environment described in the following table. The user of the AV3 aid must ensure that it is used in the appropriate environment.

Emission test	-	Recommended electromagnetic environment
Radiated / Conducted RF Emissions CISPR 11 / EN 55011	Group 1	AV3 devices use RF energy only for its internal functions. Therefore, its RF emissions are very low and are unlikely to cause any interference in nearby electronics. Remote control versions emit RF energy from the remote control
Radiated / Conducted RF Emissions CISPR 11 / EN 55011	Class B	The AV3 devices are suitable for use in all environments, including domestic environments and all those directly connected to the low voltage power supply network that supplies buildings to homes.
Harmonic emissions IEC 61000-3-2	Class A	
Voltage fluctuations and flickering EN 61000-3-3	Compliant	

13.3 RF SIGNAL TRANSMISSION AND RECEPTION

AV3 models that can be remote controlled have a receiver and an FCC certified commercial RF remote control. RF energy is transmitted from the remote control and received by the receiver inside the AV3.

The RF signal has the following characteristics

- Transmission and reception frequency: 433.92 MHz
 - Modulation: FM.
 - Encoding: Rolling code FOBLOQ
- The remote control complies with the essential requirements of the 2014/53 / EU regulation, identified as a narrow range SRD device and complies with the harmonized standards:
- EN 60950-1: 2006 + A2: 2013
 - ETSI EN 301 489-1 V2.1.1 (2017-02)
 - ETSI EN 301 489-3 V2.1.1 (2017-03)
 - ETSI EN 300 220-2 V3.1.1 (2017 (02)
 - EN 50581:2012

13.4 ELECTROMAGNETIC IMMUNITY

AV3 is intended for use both in the home and in the clinical setting, in the electromagnetic environment described in the following table. The user of the AV3 aid must ensure that it is used in the appropriate environment. Portable and Mobile RF communications equipment can affect medical electrical equipment.

Immunity test	IEC / EN 60601-1-2 test levels	Electromagnetic environment - Recommendations
Electrostatic discharge (ESD) IEC/EN 61000-4-2	±8 kV contact ±2 kV, ±4 kV, ±8 kV, ±15 kV air	Floors should be of wood, concrete or ceramic tiles. If the floors are covered with synthetic material, the relative humidity must be at least 30%
Fast Transient / Electric Burst IEC/EN 61000-4-4	± 2 kV for power supply lines ± 1kV for input/output lines PRF AT 100kHz	Mains power quality should conform to typical commercial or hospital applications.
Overvoltage surge IEC/EN 61000-4-5	±1 kV differential mode 1.2 / 50us ±2 kV common mode 0 °, 90 °, 180 °, 270 °	Mains power quality should conform to typical commercial or hospital applications.
Voltage dips, short interruptions and voltage variations on IEC / EN 61000-4-11 power lines	0% UT (100% dip in UT) for 0.5 cycles 0°, 45°, 90°, 135°, 180°, 225°, 270°, 315° 0% UT (100% dip in UT) for 1 period 70% UT (30% hole in UT) for 25 cycles 0% UT (> 100% dip in UT) for 5s	Mains power quality should conform to typical commercial or hospital applications. If the user of the device needs continuous operation of the system even when the power is interrupted, it is recommended to power it through a UPS. (battery charger)
Magnetic fields at mains frequency (50/60 Hz) IEC 61000-4-8	30A/m 50Hz	Network frequency magnetic fields must be tuned to levels typical of a standard network used for commercial or hospital use.

Immunity test	IEC / EN 60601-1-2 levels	Electromagnetic environment - Recommendations
RF Conduits IEC/EN 61000-4-6 RF Radiate IEC/EN 61000-4-3	3 Vrms Mod. AM 80% 150 kHz - 80 MHz 6 Vrms in ISM band 3 V/m Mod. AM 80% 1kHz 80 MHz - 2,7 GHz	Portable and mobile RF communications equipment should not be used near any part of the AV3 device, including cables, within the recommended separation distance calculated from the equation applicable to the transmitter frequency. Recommended separation distance $d = 1.16 \sqrt{P}$. $d = 1.16 \sqrt{P}$. 80 MHz - 800 MHz $d = 2.33 \sqrt{P}$. 800 MHz - 2.7 GHz Where P is the maximum output power rating of the transmitter in watts (W) according to the transmitter manufacturer and d is the recommended separation distance in meters (m). The energy field of fixed RF transmitters, as determined from an electromagnetic field study ¹⁰ , must be less than the compliance level in each frequency range ^b . Interference may occur near devices marked with the following symbol: 

13.5 RECOMMENDED SEPARATION DISTANCES

AV3 is intended for use in an electromagnetic environment in which radiated RF disturbances are controlled. The customer or the user of the AV3 device can help prevent electromagnetic interference by maintaining a minimum distance between portable and mobile RF communications transmitters and the AV3 device as recommended below, in relation to the maximum output power of the device. communication.

Mobile, radio frequency communication devices (including antenna cables and external antennas) should not be used closer than 30 cm (12 inches) to any part of the AV3 device,

including cables specified by the manufacturer. Otherwise, degradation of the performance of the device may occur

Recommended separation distances between portable and mobile RF communications equipment and the monitor			
Estimated Maximum Output Power of the Transmitter (W)	Separation distance according to frequency of transmitter (m)		
	150 kHz~80 MHz $d = 1.16\sqrt{P}$	from 80 MHz to 800 MHz $d = 1.16\sqrt{P}$	from 800 MHz to 2.7 GHz $d = 2.33\sqrt{P}$
0.01	0.116→ 0.3	0.116→ 0.3	0.233→ 0.3
0.1	0.368	0.366	0.737
1	1,166	1,166	2,333
10	3,689	3,689	7,378
100	11,666	11,666	23,333

For transmitters considered at maximum output power not listed above, the recommended separation distance in meters (m) can be calculated using the equation applicable to the transmitter frequency, where P is the estimate of the transmitter's maximum output power in Watts (W) according to the transmitter manufacturer.

NOTE 1 At 80 MHz and 800 MHz, the separation distance is applied for the higher frequency range.

NOTE 2 These guidelines cannot be applied in all situations. Electromagnetic propagation is subject to absorption and reflection by structures, objects and people.

Electric height adjustment does not work	LEDs under the chassis are off)	
	Motor connector is disconnected or not connected properly	Check the connections and connect the motor correctly
	The batteries are flat (red LED on)	Recharging your batteries
	Main fuse blown (3 LEDs off) o Control box or manual override or lift actuator is faulty	Contact your local distributor or Chinesport customer service
Gas spring lift does not work or Height adjustment with gas spring does not work	Emergency button pressed (Electric lift only works when you press the button to lower)	Stop pressing the emergency button or remove any obstacles that prevent the button from loosening
	The gas spring is broken	Contact your local distributor or Chinesport customer service
The device does not move (assisted mobility)	The cable from the lever to the gas spring is broken	Replace the cable
	The wheels are blocked or there are obstacles on the wheels	Unlock the wheels or remove obstacles from the wheels
The device does not move (autonomous mobility)	The main switch is set to "0", (The 3 LEDs under the chassis are off)	Turn on the main switch by positioning it on "I"
	Joystick is turned off (power button not lit)	Turn on the Joystick (See7.3.1)
	The wheels are blocked or there are obstacles on the wheels	Unlock the wheels or remove obstacles from the wheels
	Batteries are flat	Recharging your batteries
The device does not move (autonomous mobility)	Joystick is stuck (power button is on but odd LEDs are blinking)	Unlock the Joystick (See7.3.1.3)
	Control unit or joystick are broken	Contact your local distributor or Chinesport customer service

14. TROUBLESHOOTING

If you encounter a problem or malfunction, try the following troubleshooting table. If none of the proposed solutions work, contact your local distributor or Chinesport customer service.

PROBLEM	CAUSE	SOLUTIONS
The batteries do NOT charge when the charger is connected	The main switch is set to "0", (Green LED on the battery charger)	Disconnect the battery charger, turn on the device with the main switch on 'I' and connect the battery charger again
	The batteries died (Green LED on the battery charger, Main switch on "I" and device not working)	Try charging the batteries for 12 hours. If it has no effect, replace the batteries
	The battery charger is not connected correctly or is not powered (the LED on the battery charger is off)	Check connections Check the mains voltage or change the socket Check the power cable condition regularly.
	The battery charger is broken (connected correctly but the LED is off)	Replace the battery charger
Electric lift does not work or	The main switch is set to "0", (The 3	Turn on the main switch by positioning it on "I"

15. WARRANTY

15.1 GENERAL TERMS AND CONDITIONS

All Chinesport products are guaranteed against defects in materials or manufacturing for a period of 24 months from the date of sale of the product, subject to any exclusions, limitations or conditions defined in the supply phase of the product itself.

The warranty is not to be considered valid in case of improper use, tampering with the device, abuse or modification of the

product or for any use or operation not explicitly indicated in this manual.

The guarantee is not valid if the device has not been subjected to correct and documented maintenance as described in this manual or if the instructions regarding storage, cleaning and sanitation are not respected.

The manufacturer is not responsible for damage or injury or any situation caused by incorrect installation or configuration of the device or by use of the equipment that does not comply with what is reported and provided for in the installation, assembly and use manuals.

The manufacturer does not guarantee its products against defects or damage in the presence of extraordinary conditions such as: natural disasters, unauthorized maintenance and repairs, non-compliant power supply (where provided), use of non-original parts or components or accessories, non-original shipping damage directly managed by the manufacturer, lack of maintenance, manifest negligence on the part of the user or operator.

The warranty does not cover consumables, rechargeable batteries, and in general all materials subject to wear, failures caused by shocks, falls, incorrect or improper use, accidental events, damage caused by transport. If the equipment is tampered with, the warranty automatically lapses.

15.2 WARRANTY REPAIRS

In the event of a report of material or manufacturing defects, the manufacturer assesses whether the defect is covered by the warranty

The repairs under warranty must be requested expressly and are intended at our laboratory, subject to authorization and with the issue of the return number.

For products sent with the original packaging, the return shipment will be made carriage paid.

For warranty repairs, a fiscal document is required where the purchase date falls within the warranty period (sales note, purchase invoice, receipt).

The costs related to labor for the repair under warranty (when the validity of the warranty conditions is ascertained) are borne by the manufacturer.

The repair of a product or the replacement of the same does not renew or extend the terms and expirations of the warranty.

15.3 OUT OF WARRANTY REPAIRS

Products out of warranty can be repaired by the manufacturer by returning them after authorization of the same by the technical assistance service. The repair costs including shipping, material and labor are to be understood as the responsibility of the customer or the retailer. The parts and components to be repaired are to be considered covered by a guarantee for 24 months from the date of receipt of the repaired device.

15.4 NON-DEFECTIVE PRODUCTS

In the event that the manufacturer does not find malfunctions or defects in the returned products, it is

concluded that the product is not to be considered defective. The costs of shipping and handling the device will be charged to the customer or distributor.

15.5 REPAIRS AT HOME

In the case of repairs at the customer's, a written request must be made indicating the complete personal details of the applicant, the type of machine and the fault.

The mileage cost of the technician's transfer is to be agreed in relation to the customer's urgency.

In the event that the machine in question is under warranty, only the costs of the transfer will be charged.

The time is counted from the departure of the technician from our laboratory until his return, the return time will be estimated on the basis of the time taken for the outward journey

15.6 SPARE PARTS

You can ask the manufacturer for a detailed list of all spare parts.

Spare parts are sold following a formal request for an offer to the technical assistance service. The processing times are related to the availability of the parts. Returns are not accepted for replacement parts.

Payment will be cash on delivery, unless otherwise agreed.

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