

## Let's Go Out

Outdoor rollator

Black/Silver Art. No. 0120

Black Night Art. No. 0140

Beige/Silver Art. No. 0150

Red Art. No. 0130



Let's Go Out has a rubber suspension in the rear wheels for smoother walk.

Let's Go Out rollator is light and steady. It weighs only 6.4 kg and has a width of only 58 cm. It is easy to fold and the bag is folded together with the rollator. Let's Go Out is available in four colours: Black, Black Night, Red and Beige. The large front wheels make it easy to turn and climb kerbs. Let's Go Out has fenders on the rear wheels and built-in spring suspension that provides comfortable walking. The rubber handles are ergonomically shaped and can be adjusted in height from 78 to 96 cm to fit different users.



When Let's Go Out is folded it is easy to bring with you in the car or on a bus as the width is only 20 cm.



Bag – black lid. Art. No. 0121  
Bag – brown lid. Art. No. 0151



Transportation bag. Art. No. 0125



Crutch holder. Art. No. 0126



Tray. Art. No. 0479



Black backrest band. Art. No. 0123  
Brown backrest band. Art. No. 0153

### Technical specifications:

Tested according to ISO 11199-2:2005

Max. user weight	130 kg
Width	58 cm
Length	65 cm
Weight	6.4 kg
Adjustable handle height	78 – 96 cm
Seat height	59.5 cm
Front wheel diameter	24.5 cm
Rear wheel diameter	20 cm
Max. bag load	10 kg

