



SL52045325397001TX

Date: January 06,2021

Page 1 of 3

ELDO CARE CORPORATION LIMITED

RM.1902, EASEY COMM. BLDG., 253-261 HENNESSY ROAD, WANCHAI, HONGKONG

The following sample(s) was/were submitted and identified on behalf of the client as:

Sample Description	:	(A)GIMA BRANDED 20701 - 3 PLY SURGICAL MASK-DISPOSABLE - EARLOOP-BLACK
Composition Sample Color	:	(A)67% non-woven, 33% melt-blown (A)BLACK
Style No. Lot No. Manufacturer	: : :	ELFM3600 Not provided ELDO CARE CO., LIMITED
Test Performed	:	Selected test(s) as requested by applicant
Sample Receiving Date Testing Period Test Result(s)	: : :	Dec 22, 2020 Dec 23, 2020 - Jan 06, 2021 Unless otherwise stated the results shown in this test report refer only to the sample(s) tested, for further details, please refer to the following page(s).

Signed for and on behalf of SGS-CSTC Standards Technical Services (Shanghai) Co., Ltd Testing Center

Sara Guo (Account Executive)

Juan Helen

Helen Xuan (Authorized Signatory)



Unless otherwise agreed in writing, this document is issued by the Company subject to its General Conditions of Service printed overleaf, available on request or accessible at http://www.sgs.com/en/Terms-and-Conditions.aspx and, for electronic format documents, subject to Terms and Conditions for Electronic Documents at http://www.sgs.com/en/Terms-and-Conditions/Terms-e-Document.aspx. Attention is drawn to the limitation of liability, indemnification and jurisdiction issues defined therein. Any holder of this document is advised that information contained hereon reflects the Company's findings at the time of its intervention only and within the limits of Cilent's instructions, if any. The Company's sole responsibility is to its Cilent and this document cannot be reproduced except in full, without prior written approval of the Company. Any unauthorized alteration, forgery or faisification of the content or appearance of this document is unlawful and offenders may be prosecuted to the fullest extent of the law. Unless otherwise stated the results shown in this test report reter only to the sample(s) tested and such sample(s) are retained tor 30 days only. Attention: To check the authenticity of testing /inspection report & certificate, please contact us at telephone: (86-755) 8307 1443, or email: CN\_Doccheck@esps.com

3rd Building, No. 889, Yishan Road, Xuhui District Shanghai, China 200233 中国・上海・徐汇区宜山路889号3号楼 邮编: 200233 t (86-21) 61402666 f (86-21) 64958763 t (86-21) 61402666 f (86-21) 64958763



SL52045325397001TX

Test Result

## EN 14683:2019+AC:2019 Medical Face Masks-Requirements and Test Methods

#### Clause 5.2 Performance Requirement

## Clause 5.2.2 Bacterial Filtration Efficiency (BFE)

(EN 14683:2019+AC:2019 Annex B)

Dimensions of test specimen:~174mm x 162mmPositive Control Average:2636 CFUNegative Monitor Count:< 1 CFUMean Particle Size:3.0 ±0.3µmTest bacteria:Staphylococcus aureus ATCC 6538	Positive Control Average Negative Monitor Count Mean Particle Size	:	2636 CFU < 1 CFU 3.0 ±0.3μm
---	--	---	-----------------------------------

Test Item	Specimen No.	Result
	1	99.8%
Destarial Filtration Efficiency	2	99.8%
Bacterial Filtration Efficiency (BFE)	3	99.8%
(BFE)	4	99.7%
	5	99.8%

Remark:

- 1) Performance Requirement: Type I≥95%, Type II≥98%, Type IIR ≥98%
- 2) The number of specimens that shall be tested is minimum 5, but can be greater and shall be increased if necessary to allow for an AQL(Acceptable Quality Level) of 4%.



Unless otherwise agreed in writing, this document is issued by the Company subject to its General Conditions of Service printed overleaf, available on request or accessible at http://www.gg.com/en/Terms-and-Conditions.aspx and, for electronic format documents, subject to Terms and Conditions for Electronic Documents at http://www.gg.com/en/Terms-and-Conditions/Terms-e-Document.aspx. Attention is drawn to the limitation of liability, indemnification and jurisdiction issues defined therein. Any holder of this document is advised that information contained hereon reflects the Company's findings at the time of its intervention only and within the limits of Client's instructions, if any. The Company's sole responsibility is to its Client and this document cannot be reproduced except in full, without prior written approval of the Company. Any unauthorized alteration, forgery or faisification of the sand the essues the start company. Any unauthorized alteration, forgery or faisification of the sand the solution of the law. Unless otherwise stated the results shown in this test report refer only to the sample(s) tested and such sample(s) are catalized to advised.

3<sup>2\*</sup>Building,No.889,Yishan Road,Xuhui District Shanghai,China 200233 t (86–21) 61402666 f (86–21) 64958763 www.sgsgroup.com.cn 中国・上海・徐江区宜山路889号3号楼 邮编: 200233 t (86–21) 61402666 f (86–21) 64958763 e sgs.china@sgs.com



SL52045325397001TX

Date: January 06,2021

**Sample Photo** 



The statement of conformity in this test report is only based on measured values by the laboratory and does not take their uncertainties into consideration.

\*\*\*End of Report\*\*\*



Unless otherwise agreed in writing, this document is issued by the Company subject to its General Conditions of Service printed overleaf, available on request or accessible at http://www.sgs.com/en/Terms-and-Conditions.aspx and, for electronic format documents, subject to Terms and Conditions for Electronic Documents at http://www.sgs.com/en/Terms-and-Conditions/Terms-e-Document.aspx. Attention is drawn to the limitation of liability, indemnification and jurisdiction issues defined therein. Any holder of this document is advised that information contained hereon reflects the Company's findings at the time of its intervention only and within the limits of Client's instructions, if any. The Company's sole responsibility is to its Client and this document does not exonerate parties to a transaction forour price of this document cannot be reproduced except in full, without prior written approval of the Company. Any unauthorized alteration, forgery or falsification of the content or appearance of this document tast report refer only to the sample(s) tested and such sample(s) are retained for 30 days only. Attention: To check the authenticity of testing /inspection report & certificate, please contact us at telephone: (86-755) 8307 1443, or email: CM.Doccheck@sgs.com

3<sup>rd</sup>Building,No.889,Yishan Road,Xuhui District Shanghai,China 200233 t (86-21) 61402666 f (86-21) 64958763 中国・上海・徐汇区宜山路889号3号楼 邮编: 200233 f (86-21) 64958763 t (86-21) 61402666





SL52045325396201TX

ELDO CARE CORPORATION LIMITED RM.1902, EASEY COMM. BLDG., 253-261 HENNESSY ROAD, WANCHAI, HONGKONG

The following sample(s) was/were submitted and identified on behalf of the client as:

Sample Description	:	(A)GIMA BRANDED 20701 - 3 PLY SURGICAL MASK-DISPOSABLE - EARLOOP-BLACK
Composition Sample Color Style No. Lot No. Manufacturer	::	(A)67% non-woven, 33% melt-blown (A)BLACK ELFM3600 Not provided ELDO CARE CO., LIMITED
Test Performed	:	Selected test(s) as requested by applicant
Sample Receiving Date Testing Period Test Result(s)	: : :	Dec 22, 2020 Dec 23, 2020 - Jan 06, 2021 Unless otherwise stated the results shown in this test report refer only to the sample(s) tested, for further details, please refer to the following page(s).

Signed for and on behalf of SGS-CSTC Standards Technical Services (Shanghai) Co., Ltd Testing Center

Sara Guo (Account Executive)



Unless otherwise agreed in writing, this document is issued by the Company subject to its General Conditions of Service printed overleaf, available on request or accessible at http://www.sgs.com/en/Terms-and-Conditions.aspx and, for electronic format documents, subject to Terms and Conditions for Electronic Documents at http://www.sgs.com/en/Terms-and-Conditions/Terms-e-Document.aspx. Attention is drawn to the limitation of liability, indemnification and jurisdiction issues defined therein. Any holder of this document is advised that information contained hereon reflects the Company's findings at the time of its intervention only and within the limits of Client's instructions, if any. The Company's sole responsibility is to its Client and this document cannot be reproduced except in full, without prior written approval of the Company Any unauthorized alteration, forgery or faisification of the stated the sease stated the full, stated on the advise at the test report refer only to the sample(s) test and such sample(s) are retained for 30 avays only.

t (86-21) 61402666

t (86-21) 61402666

3<sup>rd</sup>Building,No.889,Yishan Road,Xuhui District Shanghai,China 200233 中国 · 上海 · 徐汇区宜山路889号3号楼 邮编: 200233 f (86–21) 64958763 f (86–21) 64958763

Member of the SGS Group (SGS SA)



SL52045325396201TX

Test Result

# EN 14683:2019+AC:2019 Medical Face Masks-Requirements and Test Methods

## Clause 5.2 Performance Requirement

## Clause 5.2.3 Breathability

(EN 14683 :2019+AC:2019 Annex C)

Sample: A Test Side

Test Side	:	Randomly test in different location (1 around and 4 away from the centric point) on each of the 5 masks
Pre-Conditioning	:	Minimum of 4 hours at 21±5°C and 85±5% R.H.
Test Area	:	4.9 cm <sup>2</sup>
Flow Rate	:	8 l/min

Specimen No.	Test Area No.	Different Pressure for each tested area (Pa/cm <sup>2</sup> )	The average value for each test specimen (Pa/cm <sup>2</sup> )
	1-1	49.2	specifien (r a/citr )
	1-2	58.0	-
1	1-3	48.1	52
	1-4	48.5	
	1-5	54.6	
	2-1	41.5	
	2-2	51.9	
2	2-3	45.6	48
	2-4	46.4	
	2-5	54.8	
	3-1	54.1	
	3-2	47.1	
3	3-3	54.9	50
	3-4	47.5	
	3-5	46.8	
	4-1	58.1	
4	4-2	58.0	
	4-3	56.6	57
	4-4	52.1	
	4-5	58.3	
	5-1	58.1	
	5-2	47.6	
5	5-3	53.2	51
	5-4	41.3	
	5-5	57.2	

Remark:

1) Performance Requirement: Type I<40 Pa/cm<sup>2</sup>, Type II<40 Pa/cm<sup>2</sup>, Type IIR<60 Pa/cm<sup>2</sup>

中国・上海・徐汇区宜山路889号3号楼 邮编: 200233

2) The number of specimens that shall be tested is minimum 5, but can be greater and shall be increased if necessary to allow for an AQL(Acceptable Quality Level) of 4%.



Unless otherwise agreed in writing, this document is issued by the Company subject to its General Conditions of Service printed overleaf, available on request or accessible at http://www.sgs.com/en/Terms-and-Conditions.aspx and, for electronic format documents, subject to Terms and Conditions for Electronic Documents at http://www.sgs.com/en/Terms-and-Conditions/Terms-en-Document.aspx. Attention is drawn to the limitation of liability, indemnification and jurisdiction issues defined therein. Any holder of this document is advised that information contained hereon reflects the Company's findings at the time of its intervention only and within the limits of client's instructions, if any. The Company's sole responsibility is to its Client and this document cannot be reproduced except in full, without prior writen approval of the Company. Any unauthorized alteration, forgery or falsification of the sample(s) tested and such sample(s) are retained for 30 days only. Attention: To check the authenticity of testing /inspection report & certificate, please contact us at telephone: (86-755) 8307 1443, or email: <u>CN. Doccheck@sgs.com</u>. [3\*BuildingNo.889,/tshan Roa/,Xuhu Dishid Shanghai,China 200233 t (86-21) 61402666 f (86-21) 64958763 www.sgsgroup.com.cn

t (86-21) 61402666

Member of the SGS Group (SGS SA)

e sgs.china@sgs.com

f (86-21) 64958763





The statement of conformity in this test report is only based on measured values by the laboratory and does not take their uncertainties into consideration.

\*\*\*End of Report\*\*\*



Unless otherwise agreed in writing, this document is issued by the Company subject to its General Conditions of Service printed overleaf, available on request or accessible at <a href="http://www.sgs.com/en/Terms-and-Conditions.aspx">http://www.sgs.com/en/Terms-and-Conditions.aspx</a> and, for electronic format documents, subject to Terms and Conditions for Electronic Documents at <a href="http://www.sgs.com/en/Terms-and-Conditions">http://www.sgs.com/en/Terms-and-Conditions.aspx</a> and, for electronic format documents, subject to Terms and Conditions for Electronic Documents at <a href="http://www.sgs.com/en/Terms-and-Conditions/Terms-e-Document.aspx">http://www.sgs.com/en/Terms-and-Conditions/Terms-e-Document.aspx</a>. Attention is drawn to the limitation of liability, indemnification and jurisdiction issues defined therein. Any holder of this document is advised that information contained hereon reflects the Company's findings at the time of its intervention only and within the limits of Client's instructions, if any. The Company's sole responsibility is to its Client and this document cannot be reproduced except in full, without prior written approval of the Company. Any unauthorized alteration, forgery or falsification of the content or appearance of this document cannot be reproduced and such sample(s) tested and such sample(s) are retained tor 30 days only. Attention: To check the authenticity of testing /inspection report & certificate, please contact us at telephone: (86-755) 8307 1443, or\_email: CN.Doccheck@sgs.com</a>

3<sup>rd</sup>Building,No.889,Yishan Road,Xuhui District Shanghai,China 200233 中国・上海・徐汇区宜山路889号3号楼 邮编: 200233 t (86-21) 61402666 f (86-21) 64958763 t (86-21) 61402666 f (86-21) 64958763

**Test Report** 





SL52045325398401TX

## ELDO CARE CORPORATION LIMITED RM.1902, EASEY COMM. BLDG., 253-261 HENNESSY ROAD, WANCHAI, HONGKONG

The following sample(s) was/were submitted and identified on behalf of the client as:

Sample Description	:	(A)GIMA BRANDED 20701 - 3 PLY SURGICAL MASK-DISPOSABLE - EARLOOP-BLACK
Composition	:	(A)67% non-woven, 33% melt-blown
Sample Color	:	(A)BLACK
Style No.	:	ELFM3600
Lot No.	:	Not provided
Manufacturer	:	ELDO CARE CO., LIMITED
Test Performed	:	Selected test(s) as requested by applicant
Sample Receiving Date	:	Dec 22, 2020
Testing Period	:	Dec 23, 2020 - Jan 07, 2021
Test Result(s)	:	Unless otherwise stated the results shown in this test report refer only to the sample(s) tested, for further details, please refer to the following page(s).

Signed for and on behalf of

SGS-CSTC Standards Technical Services (Shanghai) Co., Ltd Testing Center

Sara Guo (Account Executive)



Unless otherwise agreed in writing, this document is issued by the Company subject to its General Conditions of Service printed overleaf, available on request or accessible at http://www.gg.com/en/Terms-and-Conditions.aspx and, for electronic format documents, subject to Terms and Conditions for Electronic Documents at http://www.gg.com/en/Terms-and-Conditions/Terms-e-Document.aspx. Attention is drawn to the limitation of liability, indemnification and jurisdiction issues defined therein. Any holder of this document is advised that information contained hereon reflects the Company's findings at the time of its intervention only and within the limits of Client's instructions, if any. The Company's sole responsibility is to its Client and this document cannot be reproduced except in full, without prior written approval of the Company. Any unauthorized alteration, forgery or faisification of the sance stated the results shown in this test report refer only to the sample(s) test and such sample(s) are contact to a stelephone: (86-755) 8307 1443, or email: CN.Doccheck@ggs.com

3<sup>st</sup>Building,No.889,Yishan Road,Xuhui District Shanghai,China 200233 中国・上海・徐汇区宜山路889号3号楼 邮编: 200233

port & certificate, please contact us at t (86–21) 61402666 f (86–21) 64958763 t (86–21) 61402666 f (86–21) 64958763



SL52045325398401TX

Test Result

# EN 14683:2019+AC:2019 Medical Face Masks-Requirements and Test Methods

:

:

## Clause 5.2 Performance Requirement

## Clause 5.2.4 Splash Resistance

(ISO 22609 :2004)

Sample: A Test Blood Pressure Pre-Conditioning Distance of the mask to the tip of cannula

16.0kPa Minimum of 4 hours at 21±5°C and 85±5% R.H. 300±10mm

Test Specimen#	Penetration on inside surface	Conclusion	Test Specimen#	Penetration on inside surface	Conclusion	
1	None Seen	Pass	17	None Seen	Pass	
2	None Seen	Pass	18	None Seen	Pass	
3	None Seen	Pass	19	None Seen	Pass	
4	None Seen	Pass	20	None Seen	Pass	
5	None Seen	Pass	21	None Seen	Pass	
6	None Seen	Pass	22	None Seen	Pass	
7	None Seen	Pass	23	None Seen	Pass	
8	None Seen	Pass	24	None Seen	Pass	
9	None Seen	Pass	25	None Seen	Pass	
10	None Seen	Pass	26	None Seen	Pass	
11	None Seen	Pass	27	None Seen	Pass	
12	None Seen	Pass	28	None Seen	Pass	
13	None Seen	Pass	29	None Seen	Pass	
14	None Seen	Pass	30	None Seen	Pass	
15	None Seen	Pass	31	None Seen	Pass	
16	None Seen	Pass	32	None Seen	Pass	
Number	Number of Pass:		32			
Overal	Overall result:		Acce	ptable		

Remark:

- 1) Performance Requirement Type I: N/A, Type II: N/A, Type IIR: ≥16.0kPa
- 2) Test was conducted within 60s after removal from conditioning chamber.
- 3) An acceptable quality limit of 4.0% is met for a single sampling plan when 29 or more of the 32 tested specimens show pass results.



Unless otherwise agreed in writing, this document is issued by the Company subject to its General Conditions of Service printed overleaf, available on request or accessible at http://www.sgs.com/en/Terms-and-Conditions.aspx and, for electronic format documents, subject to Terms and Conditions for Electronic Documents at http://www.sgs.com/en/Terms-and-Conditions/Terms-e-Document.aspx. Attention is drawn to the limitation of liability, indemnification and jurisdiction issues defined therein. Any holder of this document is advised that information contained hereon reflects the Company's findings at the time of its intervention only and within the limits of Client's instructions, if any. The Company's sole responsibility is to its Client and this document cannot be reproduced except in full, without prior written approval of the Company any unauthorized alteration, forgery or falsification of the sate stated the results shown in this test report refer only to the sample(s) tested and such sample(s) are contact to adv asy only.

3<sup>rd</sup>Building,No.889,Yishan Road,Xuhui District Shanghai,China 200233 t (86-21) 61402666 中国・上海・徐江区宜山路889号3号楼 邮编: 200233 t (86-21) 61402666

f (86–21) 64958763 www.sgsgroup.com.cn f (86–21) 64958763 e sgs.china@sgs.com



The statement of conformity in this test report is only based on measured values by the laboratory and does not take their uncertainties into consideration.

\*\*\*End of Report\*\*\*



Unless otherwise agreed in writing, this document is issued by the Company subject to its General Conditions of Service printed overleaf, available on request or accessible at http://www.sgs.com/en/Terms-and-Conditions.aspx and, for electronic format documents, subject to Terms and Conditions for Electronic Documents at http://www.sgs.com/en/Terms-and-Conditions/Terms-e-Document.aspx. Attention is drawn to the limitation of liability, indemnification and jurisdiction issues defined therein. Any holder of this document is advised that information contained hereon reflects the Company's findings at the time of its intervention only and within the limits of Client's instructions, if any. The Company's sole responsibility is to its Client and this document does not exonerate parties to a transaction forour price of this document cannot be reproduced except in full, without prior written approval of the Company. Any unauthorized alteration, forgery or falsification of the content or appearance of this document tast report refer only to the sample(s) tested and such sample(s) are retained for 30 days only. Attention: To check the authenticity of testing /inspection report & certificate, please contact us at telephone: (86-755) 8307 1443, or email: CM.Doccheck@sgs.com

3<sup>rd</sup>Building,No.889,Yishan Road,Xuhui District Shanghai,China 200233 t (86-21) 61402666 中国・上海・徐汇区宜山路889号3号楼 邮编: 200233 t (86-21) 61402666 f (86-21) 64958763

f (86-21) 64958763 www.sgsgroup.com.cn e sgs.china@sgs.com





# SL52045325391701TX

Date: January 06,2021

ELDO CARE CORPORATION LIMITED

RM.1902, EASEY COMM. BLDG., 253-261 HENNESSY ROAD, WANCHAI, HONGKONG

The following sample(s) was/were submitted and identified on behalf of the client as:

Sample Description	:	(A)GIMA BRANDED 20701 - 3 PLY SURGICAL MASK-DISPOSABLE - EARLOOP-BLACK
Composition	:	(A)67% non-woven, 33% melt-blown
Sample Color	:	(A)BLACK
Style No.	:	ELFM3600
Lot No.	:	Not provided
Manufacturer	:	ELDO CARE CO., LIMITED
Test Performed	:	Selected test(s) as requested by applicant
Sample Receiving Date	:	Dec 22, 2020
Testing Period	:	Dec 23, 2020 - Jan 06, 2021
Test Result(s)	:	Unless otherwise stated the results shown in this test report refer only to the sample(s) tested, for further details, please refer to the following page(s).

Signed for and on behalf of SGS-CSTC Standards Technical Services (Shanghai) Co., Ltd Testing Center

Sara Guo (Account Executive)

Vorijin

Dongjing Liu (Authorized Signatory)



Unless otherwise agreed in writing, this document is issued by the Company subject to its General Conditions of Service printed overleaf, available on request or accessible at http://www.sgs.com/en/Terms-and-Conditions.aspx and, for electronic format documents, subject to Terms and Conditions for Electronic Documents at http://www.sgs.com/en/Terms-and-Conditions.aspx and, for electronic format documents, Attention is drawn to the limitation of liability, indemnification and jurisdiction issues defined therein. Any holder of this document is advised that information contained hereon reflects the Company's findings at the time of its intervention only and within the limits of Cilent's instructions, if any. The Company's sole responsibility is to its Cilent and this document does not exonerate parties to a transaction from exercising all their rights and obligations under the transaction documents. This document cannot be reproduced except in full, without prior written approval of the Company. Any unauthorized alteration, forgery or faisification of the content or appearance of this document is unlawful and offenders may be prosecuted to the fullest extent of the law. Unless otherwise stated the results shown in this test report reter only to the sample(s) tested and such sample(s) are retained tor 30 days only. Attention: To check the authenticity of testing /inspection report & certificate, please contact us at telephone: (86-755) 8307 1443, or email: CN\_Doccheck@egss.com

3<sup>rd</sup>Building, No.889, Yishan Road, Xuhui District Shanghai, China 200233 中国・上海・徐汇区宜山路889号3号楼 邮编: 200233 t (86-21) 61402666 f (86-21) 64958763 t (86-21) 61402666 f (86-21) 64958763



SL52045325391701TX

Test Result

# EN 14683:2019+AC:2019 Medical Face Masks-Requirements and Test Methods

#### Clause 5.2 Performance Requirement

#### Clause 5.2.5 Microbial Cleanliness

(EN 14683:2019+AC:2019 Annex D and EN ISO 11737-1:2018)

Sample: A

Test Specimen#	Mask Weight(g)	Total Bioburden, (CFU/mask)	Total Bioburden, (CFU/g)
1#	3.53	36	10.20
2#	3.52	15	4.26
3#	3.52	51	14.49
4#	3.51	48	13.68
5#	3.61	42	11.63

Remark: Performance Requirement: Type I≤30 CFU/g, Type II≤30 CFU/g, Type IIR≤30 CFU/g

# 

The statement of conformity in this test report is only based on measured values by the laboratory and does not take their uncertainties into consideration.

\*\*\*End of Report\*\*\*



Unless otherwise agreed in writing, this document is issued by the Company subject to its General Conditions of Service printed overleaf, available on request or accessible at http://www.gg.com/en/Terms-and-Conditions.aspx and, for electronic format documents, subject to Terms and Conditions for Electronic Documents at http://www.gg.com/en/Terms-and-Conditions/Terms-e-Document.aspx. Attention is drawn to the limitation of liability, indemnification and jurisdiction issues defined therein. Any holder of this document is advised that information contained hereon reflects the Company's findings at the time of its intervention only and within the limits of Client's instructions, if any. The Company's sole responsibility is to its Client and this document cannot be reproduced except in full, without prior written approval of the Company. Any unauthorized alteration, forgery or faisification of the sand the essues the start company. Any unauthorized alteration, forgery or faisification of the sand the solution of the law. Unless otherwise stated the results shown in this test report refer only to the sample(s) tested and such sample(s) are catalized to advised.

3<sup>16</sup>Building,No.889,Yishan Road,Xuhui District Shanghai,China 200233 中国·上海·徐汇区宜山路889号3号楼 邮编: 200233 port & certificate, please contact us at 1 t (86-21) 61402666 f (86-21) 64958763 t (86-21) 61402666 f (86-21) 64958763