



AG19

Materasso in fibra cava

Il materasso in fibra cava AG19 deve essere utilizzato su una base di appoggio (per esempio il materasso principale del letto) ed è formato da:

- un rivestimento in tessuto di cotone diviso in 13 cilindri distinti contenenti la fibra cava siliconata
- una cerniera sul lato lungo del materasso che consente l'accesso alla parte interna
- due lembi testa - piedi che assicurano l'ancoraggio al materasso principale del letto

Hollow fibre mattress

The AG19 hollow fibre mattress must always be used on a supporting base (for example, the main mattress of the bed) and is composed of:

- a cotton fabric lining divided into 13 separate cylinders containing silicone hollow fibre
- a zipper along the long side of the mattress to access the inside of it
- two head - feet flaps to ensure the main mattress is anchored to the bed

Dati Tecnici / Technical Specifications

Lunghezza / Length	cm 195
Larghezza / Width	cm 85
Altezza / Height	cm. 13
Peso / Weight	kg 5
Portata max / Max. capacity	kg 120



AG19/E

Materasso in fibra cava

Il materasso in fibra cava AG19 deve essere utilizzato su una base di appoggio (per esempio il materasso principale del letto) ed è formato da:

- un rivestimento in tessuto di cotone diviso in 13 cilindri distinti contenenti la fibra cava siliconata
- L'imbottitura che è inserita in un secondo cilindro di TNT che può essere inter cambiato con le altre sezioni
- Una cerniera sul lato lungo del materasso che consente l'accesso alla parte interna
- Due lembi testa - piedi che assicurano l'ancoraggio al materasso principale del letto

Hollow fibre mattress

The AG19 hollow fibre mattress must always be used on a supporting base (for example, the main mattress of the bed) and is composed of:

- a cotton fabric lining divided into 13 separate cylinders containing silicone hollow fibre
- The padding that is inserted in a second Non-Woven cylinder that can be replaced with the other sections
- A zipper along the long side of the mattress to access the inside of it
- Two head - feet flaps to ensure the main mattress is anchored to the bed

Dati Tecnici / Technical Specifications

Lunghezza / Length	cm 195
Larghezza / Width	cm 85
Altezza / Height	cm. 13
Peso / Weight	kg 7
Portata max / Max. capacity	kg 120