

LIBRA Modern design and compact dimensions

Libra bedside cabinet provides all that's necessary in daily care. It combines wide amount of easily accessible storage space with modern design, compact dimensions and solid structure. Raised cabinet sides over table top prevent items from falling down and integrated slide-out table ensures increase of space when necessary. Storage space divided between upper drawer and lower container drawer – both supported by strudy and smooth-running pull-out system.



Key features:

- Compact dimensions
- Raised cabinet sides preventing items form falling of the table top
- Integrated slide-out table with maximum working load of 15kg
- Sliding drawer and sliding container drawer
- Four castors (two with brakes)
- Wide range of colours
- Guarantee: 2 years











LIBRA Modern design and compact dimensions

Raised cabinet sides - to prevent items form falling of the table top



BASIC DIMENSIONS	
Ext. dim. (width x height x depth)	42,5 x 83 x 47cm
Slide-out table (width x depth)	39x36cm
Castors with brakes	2




Adjust design to individual needs

Wide range of wood colours:





REHA-BED Sp. z o.o. ul. Spacerowa 1, 41-253 Czeladź





) biuro@rehabed.com.pl

