

L'eccellenza per la mobilità di anziani e disabili su lunghi percorsi
The excellence for mobility of elderly and disabled people on long distances




SURACE ON AIR Forniture per aeroporti, ferrovie e turismo Supplies for airports, railways and tourism

Pedana per conducente
 chiudibile
 36 cm
 Folding driver's platform
 36 cm



Entra in tutti gli ascensori per disabili.
Enters into all disabled elevators.



SKYWAY ONE
 Sedia elettronica da transito con conducente a bordo

Standard: Motori da 330 watt - velocità 6 km/h. Batterie 100 AH senza manutenzione. Comando con clacson, regolatore di velocità. Chiave di sicurezza e fungo di emergenza. Luci regolamentari con lampeggiante. Vano portabagaglio più supporto bagaglio supplementare e 2 ganci porta zaino. Tasca sullo schienale. Pedana per conducente chiudibile. Braccioli ribaltabili e regolabili in larghezza. Appoggiatesta girevoli e regolabili. Seduta ergonomica igienizzabile con rivestimenti ignifughi. Caricabatterie a bordo con cavo auto avvolgente metri 8. Cintura di sicurezza autorilavante.

Optional: Motori maggiorati da 480 watt, velocità max 10 km/h. Batterie al litio (+60% durata -18 kg peso). Capote. Appoggiatesta. Divaricatore. Sostegni per il tronco.

Guarda il video sul link <http://youtu.be/tv075MyjYPM>

Codice LEA 12.23.06.009 (ricond. 015)

SKYWAY ONE
 Electronic travel chair with driver on board

Standard: 330 watt motors - 6 km/h speed. 100 AH battery without maintenance. Controller with horn and speed regulator. Security key and emergency mushroom head button. Regulatory lights flashing. Luggage compartment, extra luggage support and 2 hooks for backpack. Pocket at the back. Folding driver's platform. Tipping and adjustable in width armrests. Swinging and adjustable legrests. Ergonomic, sanitized seat with fireproof upholsteries. On-board battery charger with 8 meters auto rewind cable. Auto rewind security belt.

Optional: 480 watt increased motors, maximum speed 10 km/h. Lithium batteries (life 60% more, weight 18 kg less). Capote. Headrest. Thigh divider. Trunk supports.

See video on link <http://youtu.be/ROGNBIIFHIS>

ISO Class 12.21.27.006 (ricond. 009)

cm Ø 33	cm Ø 22	L 50 - 65	L 70	L 119 - 180	Weight 100 kg	150 kg Passeggero + Bagaglio Passenger + Luggage	100 kg Conducente Driver