


**Bella, svelta e pratica
larghezza esterna solo 58 cm
Nice, easy and practical
only 58 cm external width**



715 MAGIC
Sedia elettronica pieghevole

Standard: Motori da 330 watt con batterie da 100 A senza manutenzione. Comando con clacson, regolatore di velocità. Telaio pieghevole. Leggere imbottiture con tessuti ignifughi. Freni elettromagnetici. Braccioli desk regolabili in altezza ed estraibili. Pedane girevoli ed estraibili. Cintura di sicurezza. Ruote antiribaltamento. Caricabatteria. Velocità max 7 km/h. Autonomia 6h - 40 km.

Optional: Comando allo schienale.

715 MAGIC
Electronic folding chair

Standard: 330 watt motors with 100 A batteries without maintenance. Control box with horn, speed governor. Folding frame. Light padding with fire retardant upholsteries. Electromagnetic brakes. Adjustable and removable desk armrests. Swinging and detachable footrests. Security belt. Anti-overturn small wheels. Battery charger. Max speed 7 km/h. Range: 6h - 40 km.

Optional: Control box on backrest.



+ kg 30
BATTERY



Codice LEA 12.23.06.009
ISO Class 12.21.27.006

cm Ø 33	cm Ø 20	L 43	L 58	L 75 - 98	h 52	Standard	20%	Portata 120 kg