

L'indistruttibile *The unbreakable*



◀ 63 - 65,5 - 68 ▶

Larghezza braccioli
variabile

*Adjustable width
between armrests*



🇮🇹 PANZER

Poltrona bariatrica per utenti oltre 150 kg

Standard: Struttura super rinforzata a 6 ruote, in tubolare 22 e 26 x2, porta fino a 220 kg. WC igienizzabile. Super imbottiture anatomiche. Imbottiture e tessuti ignifughi. Schienale regolabile con pistoni a gas. Appoggiamambe contenitivi maggiorati, girevoli e regolabili. Pedane regolabili ed estraibili. Fianchi estraibili con paravesti. Gomme soft. Maniglione e maniglie singole di spinta allo schienale.

Optional: Tavolino sagomato. Appoggiatesta. Freni a tamburo.

🇬🇧 PANZER

Bariatric armchair for users over 150 kg

Standard: Super reinforced structure with 6 wheels, in 22 and 26 x 2 tubular, carries up to 220 kg. Easy cleanable WC. Super anatomical padding. Fire-retardant padding and upholsteries. Adjustable backrest with gas springs. Containing, oversized, swinging and adjustable in height legrests. Adjustable and removable footrests. Removable sides with clothes guards. Soft wheels. Big push handle and single push handles at the back.

Optional: Contoured table. Headrest. Drum brakes.

Codice LEA 12.22.18.015 ricond.

ISO Class 12.21.03.006 + 12.24.21.103 + 12.24.06.124

cm Ø 20	L 63 - 65,5 - 68	L 73	L 100	h 54 - 73	h 50	Standard	Optional	Optional	Optional	Portata 220 kg