


GRAVIS

available colours:



4 Wheeled Mobility Scooter **GRAVIS** is our heavy-duty disability scooter, which can be used in any road conditions thanks to its great suspension system. On top of that it is equipped with modern LED lighting which provides a lot of safety in the night and are highly visible in any conditions. Model **GRAVIS** is a perfect solution to everyone who looks for full stability and a sense of security when driving a mobility scooter.

The scooter comes **assembled and ready to go!** 

SPECIFICATIONS

OVERVIEW

Max user weight: 160 kg / 25.2 st
Scooter weight: 125 kg / 19.7 st
Hill climbing ability: 17°
Maximum range on single charge:
45-60 km / 28-37.3 mi
Speed of response: Low
Maximum speed: 8 mph / 12 km/h
Acceleration system: Right-hand
scooter throttle
Seat adjustment: Yes
Suspension: Full
Compliant with DVLA requirements:
Yes
Scooter class: Class-3
Vehicle autonomy: Level 0
Optimal air humidity for storage
and use: 30%
Working temperature of electronic
parts: -10 °C to +40 °C
Turning radius: 220 cm / 86.6"

MOTOR

Power: 1000 W
Voltage: 60 V

BATTERY

Type: Lead-Acid Gel
Capacity: 20 Ah
Voltage: 60 V
Life cycle: 600 full charges
Charge time: ≈8 h

BRAKES

Front: Hydraulic disc
Rear: Drum
Magnetic brake: No

WHEELS

Amount: 4
Type of tyres: Pneumatic tubeless
Size: F 3.00-8" / R 3.00-10"
Tyre pressure: 2.25 bar / 225 kPa /
32.5 PSI

LIGHTS

Low beam headlight: LED (integrated with lamp)
High beam headlight: LED (integrated with lamp)
Tail lights: LED (integrated with lamp)
Daytime running lights: LED (integrated with lamp)
Signal lights: LED (integrated with lamp)
Brake light: LED (integrated with lamp)
Hazard lights: LED (integrated with lamp)

FEATURES

Horn
Rearview mirrors
Steering wheel lock
Parking brake
Electric LED display
USB Port
Alarm
Remote control
Comfortable seat

GRAVIS

