

# ZT63



available colours:



Our **ZT63** retro Vespa electric scooter is a perfect choice for any customer who fancies the look straight from early and middle 1900s. The Vespa style mobility scooters are a perfect balance between old body and brand-new heart of the vehicle. This model is favoured by everyone who prefers to have their mobility vehicle not look like a mobility aid. The vehicle was designed in a way that makes it possible for two people to get on it.

The scooter comes **assembled and ready to go!** 

# SPECIFICATIONS

## OVERVIEW

Max user weight: 199 kg / 31.3 st  
Scooter weight: 94 kg / 14.8 st  
Hill climbing ability: 17°  
Maximum range on single charge:  
40-60 km / 24.9-37.3 mi  
Speed of response: Medium  
Maximum speed: 8 mph / 12 km/h  
Acceleration system: Right-hand  
scooter throttle  
Seat adjustment: No  
Suspension: Full  
Compliant with DVLA requirements:  
Yes  
Scooter class: Class-3  
Vehicle autonomy: Level 0  
Optimal air humidity for storage  
and use: 30%  
Working temperature of electronic  
parts: -10 °C to +40 °C  
Turning radius: 195 cm / 76.8"

## MOTOR

Power: 650 W  
Voltage: 60 V

## BATTERY

Type: Lead-Acid Gel  
Capacity: 20 Ah  
Voltage: 60 V  
Life cycle: 600 full charges  
Charge time: ≈8 h

## BRAKES

Front: Hydraulic disc  
Rear: Drum  
Magnetic brake: No

## WHEELS

Amount: 3  
Type of tyres: Pneumatic tubeless  
Size: 3.00-10"  
Tyre pressure: 2.25 bar / 225 kPa /  
32.5 PSI

## LIGHTS

Low beam headlight: S2  
High beam headlight: S2  
Position lights: W5W  
Tail lights: P21/5 W  
Daytime running lights: S2  
Signal lights: R10W 12V  
Brake light: P21/5 W  
Hazard lights: R10W 12V

## FEATURES

Horn  
Rearview mirrors  
Steering wheel lock  
Parking brake  
Electric LED display  
Alarm  
Remote control

# ZT63



