

ZT15



available colours:



Produced by renowned manufacturer Veledo company, **ZT15** mobility scooter is known for its high quality. Assembled in Poland and delivered fully ready-to-drive, it is a perfect solution for everyone who has to deal with daily mobility problems of any kind. Being one of the first models manufactured by the company, it is a great vehicle for first time electric mobility scooter users.

The scooter comes **assembled and ready to go!** 

SPECIFICATIONS

OVERVIEW

Max user weight: 160 kg / 25.2 st
Scooter weight: 104 kg / 16.4 st
Hill climbing ability: 17°
Maximum range on single charge:
45-60 km / 28-37.3 mi
Speed of response: Medium
Maximum speed: 8 mph / 12 km/h
Acceleration system: Right-hand
scooter throttle
Seat adjustment: Yes
Suspension: Full
Compliant with DVLA requirements:
Yes
Scooter class: Class-3
Vehicle autonomy: Level 0
Optimal air humidity for storage
and use: 30%
Working temperature of electronic
parts: -10 °C to +40 °C
Turning radius: 190 cm / 74.8"

MOTOR

Power: 900 W
Voltage: 48 V

BATTERY

Type: Lead-Acid Gel
Capacity: 20 Ah
Voltage: 48 V
Life cycle: 600 full charges
Charge time: ≈8 h

BRAKES

Front: Drum
Rear: Drum
Magnetic brake: Yes

WHEELS

Amount: 3
Type of tyres: Pneumatic tubeless
Size: 3.00-10"
Tyre pressure: 2.25 bar / 225 kPa /
32.5 PSI

LIGHTS

Low beam headlight: S2
High beam headlight: S2
Position lights: W5W
Front and rear signal lights: P21/5
W LED / R10W 12V
Brake light: P21/5 W LED
Hazard lights: P21/5 W LED / R10W
12V
Registration plate light: W5W

FEATURES

Horn
Rearview mirrors
Steering wheel lock
Parking brake
Electric LED display
Alarm
Remote control

ZT15



